PAXgene™ Blood RNA Tube

English	1-4	Nederlands	25-28
Español	5-8	Svenska	29-32
Português	9-12	Dansk	33-36
Français	13-16	Ελληνικά	37-40
Deutsch	17-20	Polski	41-44
Italiano	21-24	Česky	45-48

English

I. Intended Use

The PAXgene™ Blood RNA System consists of a blood collection tube (PAXgene™ Blood RNA Tube) and nucleic acid purification kit (PAXgene™ Blood RNA Kit). It is intended for the collection, storage, and transport of blood and stabilization of intracellular RNA in a closed tube and subsequent isolation and purification of intracellular RNA from whole blood for RT-PCR used in molecular diagnostic testing.

Performance characteristics for the PAXgene™ Blood RNA System have only been established with FOS and IL1B gene transcripts. The user is responsible for establishing appropriate PAXgene™ Blood RNA System performance characteristics for other target transcripts.

Product Features

PAXgene™ Blood RNA Tube	2.5 mL • 16 x 100 mm
762165	Plus • BD Hemogard™
	6.9 mL Additive
	100 tubes/case

II. Summary and Explanation

The collection of whole blood is the first step in many molecular assays used to study intracellular RNA. A major challenge in this type of testing is the instability of intracellular RNA, which rapidly degrades within hours after blood collection. Furthermore, certain species of RNA, through the process of gene induction, increase *in vitro* after blood collection. Both *in vitro* RNA degradation and gene induction can lead to an under- or overestimation of *in vivo* relative gene transcript number.

The PAXgene™ Blood RNA Tube contains an additive that stabilizes the *in vivo* gene transcription profile by reducing *in vitro* RNA degradation and minimizing gene induction. When used in conjunction with the PAXgene™ Blood RNA Kit, the PAXgene™ Blood RNA Tube provides for accurate detection and quantitation of gene transcripts.

III. Warnings

- Contents of this tube are irritating to skin.
 - After inhalation, supply fresh air: consult doctor in case of complaints.
 - After skin contact, immediately wash with water and soap, and rinse thoroughly.
 - After eye contact, rinse opened eye for fifteen minutes under running water, then consult a doctor.
 - After swallowing, immediately call a doctor.
- A blood collection set must be used with the PAXgene™ Blood RNA Tube. See Ordering Information.


IV. Specimen Collection and Preparation for Analysis

A. Required Blood Collection Accessories

(Not included with the PAXgene™ Blood RNA Tube)

- Blood collection set such as the BD Vacutainer® Safety-Lok™ Blood Collection Set. See Ordering Information.
- A BD Vacutainer® Needle Holder must be used to ensure proper function. See Ordering Information.
- A "Discard Tube" if the PAXgene™ Blood RNA Tube is the only tube being drawn. See Ordering Information.
- Labels for positive donor identification of samples
- Alcohol swab for cleansing site
- Dry sterile gauze
- Tourniquet
- Needle disposal container for used needle or needle/holder combination

B. Procedure for Specimen Collection

Draw Order:

Tube A - Discard Tube

Tube B - PAXgene™ Blood RNA Tube

- Ensure that the PAXgene™ Blood RNA Tube is at 18°C to 25°C prior to use and properly labeled with patient identification.
- If the PAXgene™ Blood RNA Tube is the only tube to be drawn, blood should be drawn into a "Discard Tube" prior to drawing blood into the PAXgene™ Blood RNA Tube so the interior volume of the blood collection set used during phlebotomy can be primed. Otherwise the PAXgene™ Blood RNA Tube should be the last tube drawn in the phlebotomy procedure.
- Using a blood collection set and a holder, collect blood into the PAXgene™ Blood RNA Tube using your institution's recommended procedure for standard venipuncture technique. The following techniques shall be used to prevent possible backflow:
 - Place donor's arm in a downward position.
 - Hold tube in a vertical position, below the donor's arm during blood collection.
 - Release tourniquet as soon as blood starts to flow into tube.
 - Make sure tube additives do not touch stopper or end of the needle during venipuncture.

4. Allow at least 10 seconds for a complete blood draw to take place. Ensure that the blood has stopped flowing into the tube before removing the tube from the holder. The PAXgene™ Blood RNA Tube with its vacuum is designed to draw 2.5mL of blood into the tube.

5. Immediately after blood collection, gently invert the PAXgene™ Blood RNA Tube 8 – 10 times.

6. Store the PAXgene™ Blood RNA Tube upright at room temperature (18°C to 25°C) for a minimum of 2 hours and a maximum of 72 hours before processing or transferring to refrigerator (2°C to 8°C) or freezer (-20°C). See Procedure for Freezing and Thawing Specimens in the PAXgene™ Blood RNA Tubes for details if storage temperature of -70°C/-80°C is desired.

C. Procedure for Freezing and Thawing Specimens Collected in the PAXgene™ Blood RNA Tubes

1. Stand the PAXgene™ Blood RNA Tube upright in a wire rack. **Do not freeze tubes upright in a styrofoam tray as this may cause the tubes to crack.**
2. The PAXgene™ Blood RNA Tubes can be stored at -20°C and below. If tubes are to be kept at temperatures below -20°C, freeze them first at -20°C for 24 hours, then transfer them to -70°C or -80°C.
3. Thaw the PAXgene™ Blood RNA Tubes in a wire rack at ambient temperature (18°C – 25°C) for approximately two hours. **Do not thaw the PAXgene™ Blood RNA Tubes, at temperatures above 25°C.**
4. Carefully invert the thawed PAXgene™ Blood RNA Tubes 10 times.

Note: If the two-hour ambient temperature incubation was not carried out prior to freezing, incubate tubes at least another two hours after reaching room temperature before processing.

Note: The frozen PAXgene™ Blood RNA Tubes are subject to breakage upon impact. To reduce the risk of breakage during shipment, frozen tubes should be treated in the same manner as glass tubes. Users must validate their own freezing and shipping protocol for the PAXgene™ Blood RNA Tubes.

D. Procedure for Specimen Preparation for Analysis


Blood samples must be processed according to the instructions provided with the PAXgene™ Blood RNA Kit. For catalog numbers, see Ordering Information.

V. Performance Characteristics

When filled correctly, the additive to blood ratio in the PAXgene™ Blood RNA Tube is 2.76 mL of additive per mL of blood. After blood is introduced into the tube, the intracellular RNA profile remains stable for 3 days at 18°C to 25°C (Figures 1A and 1B), 5 days at 2°C to 8°C (Figures 2A and 2B), or 6 months at -20°C, and for 6 months at -70°C / -80°C.

RNA Stability in Blood Samples at 18–25°C

A. FOS

B. IL1B



Figure 1

Blood was drawn from 10 donors, with duplicate samples, and stored at 18–25°C for the indicated number of days, followed by total RNA purification. Blood was collected and stored in PAXgene™ Blood RNA Tubes, and total RNA was purified using the PAXgene™ Blood RNA Kit. Relative transcript levels of FOS and IL1B were determined by real-time, duplex RT-PCR, using 18S rRNA as an internal standard. The values for all samples are plotted (20 data sets for each gene, blue lines), with means (red lines) and standard deviations (black bars) of all samples shown. The dashed lines indicate the $\pm 3x$ total precision of the assays (FOS: 2.34 C_T ; IL1B: 1.93 C_T).

RNA Stability in Blood Samples at 2–8°C

A. FOS

B. IL1B


Figure 2

Blood was drawn and total RNA purified, after storage at 2–8°C, as described in Figure 1. Relative transcript levels of FOS and IL1B were determined by real-time, duplex RT-PCR, using 18S rRNA as an internal standard. The values for all samples are plotted (20 data sets for each gene, blue lines), with means (red lines) and standard deviations (black bars) of all samples shown. The dashed lines indicate the $\pm 3\times$ total precision of the assays (FOS: 2.34 C_T ; IL1B: 1.93 C_T).

The actual duration of RNA stabilization may vary depending upon the species of intracellular RNA.

VI. Limitations

1. Under-filling of the PAXgene™ Blood RNA Tubes will result in an incorrect blood-to-additive ratio and may lead to incorrect analytic results or poor product performance.
2. The PAXgene™ Blood RNA System is not suitable for the collection and purification of viral RNA.
3. Blood samples obtained with the PAXgene™ Blood RNA Tube must be prepared only with the PAXgene™ Blood RNA Kit. For catalog numbers, see Ordering Information.
4. The quantity of blood drawn should be 2.5mL per PAXgene™ Blood RNA Tube, but this volume might vary with altitude, ambient temperature, barometric pressure, tube age, venous pressure, and filling technique.
5. The PAXgene™ Blood RNA Kit is intended for purification of intracellular RNA from human whole blood with leukocytes counts between $4.8 \times 10^6 - 1.1 \times 10^7$ leukocytes/mL.

VII. Precautions

1. Practice universal precautions. Use gloves, gowns, eye protection, other personal protective equipment, and engineering controls to protect from blood splatter, blood leakage, and potential exposure to bloodborne pathogens.
2. Handle all biologic samples and blood collection sets "sharps" according to the policies and procedures of your facility. Obtain appropriate medical attention in the event of any exposure to biologic samples (for example, through a puncture injury), since they may transmit viral hepatitis, HIV (AIDS), or other infectious diseases. Utilize any built-in needle protector, if the blood collection set provides one. PreAnalytiX does not recommend resheilding used needles. However, the policies and procedures of your facility may differ and must always be followed.
3. Discard all blood collection tubes in biohazard containers approved for their disposal.
4. Do not re-use the PAXgene™ Blood RNA Tubes.
5. Do not use the PAXgene™ Blood RNA Tubes after the expiration date printed on the tube label.
6. Since the PAXgene™ Blood RNA Tube contains a chemical additive, a blood collection set must be used for blood collection to prevent possible backflow from the tube. See *Specimen Collection and Preparation for Analysis* section.
7. Excessive centrifugation speed (over 10,000 RCF) may cause the PAXgene™ Blood RNA Tube breakage, exposure to blood and possible injury.
8. Do not transfer a sample from a syringe to a tube.

VIII. Storage

1. Store the unused PAXgene™ Blood RNA Tubes at 18°C to 25°C. Limited excursion temperatures up to 40°C are permitted.

The PAXgene™ Blood RNA Tube is made in the UK by BD for PreAnalytiX GmbH. PAXgene and PreAnalytiX are trademarks of PreAnalytiX GmbH all other trademarks are property of Becton, Dickinson and Company.

U. S. Patent Nos. 4,741,446, 4,991,104, 6,602,718, and 6,617,170

PreAnalytiX GmbH, 8634 Hombrechtikon, CH
©2005 PreAnalytiX GmbH www.PreAnalytiX.com

Ordering Information

PAXgene™ Products

	Contents
PAXgene™ Blood RNA Tube	100 Blood collection tubes. Use with the PAXgene™ Blood RNA Kit (50)
762165	
PAXgene™ Blood RNA Kit (50)	50 PAXgene™ Spin Columns, Processing Tubes, RNase-free Reagents and Buffers
762164 N.A.	Use with the PAXgene™ Blood RNA Tubes (100)
762174 Ex. N.A.	

BD Products*

BD Vacutainer® Safety-Lok™ Blood Collection Set	21G 3/4 inch needle, 12 inch tubing with luer adapter.
367281 US	50/box, 200/case
367286 CE	
BD Vacutainer® One Use Holder	Case only for 13mm and 16mm diameter
364815	1000/case
BD Vacutainer® Plus Serum Tubes	13 X 75mm 4.0 ml draw with Red BD Hemogard™ closure and paper label
367812 US	100/box, 1000/case
368975 CE	

* These catalog numbers represent typical products that can be used with the PAXgene™ Blood RNA Tube.

To order BD Blood Collection Accessories:

Go to: <http://www.bd.com/support/contact/international.asp>

- Select country
- Scroll down for BD Diagnostics - Preanalytical Systems
- Dial number

BD – Customer Service / Orders

Location	Phone	Fax	Email
Australia BD: North Ryde	612.8875.7000	612.8875.7000	bd_anz@bd.com
Canada BD: Oakville	800.268.5430	800.565.0897	custsvccan@bd.com
Belgium BD: Erembodegem	32.53720408	32.53720558	
Poland BD: Warsaw	48.22.651.75.88	48.22.651.75.89	
USA BD: Franklin Lakes	888.237.2762	800.847.2220	

QIAGEN – Customer Service / Orders

Location	Phone	Fax
Australia	03.9489.3666	03.9489.3888
Canada	888.387.2381	800.713.5951
France	01.60.920.920	01.60.920.925
Germany	02103.29.12000	02103.29.22000
Italy	02.33430411	02.33430426
Japan	03.5547.0811	03.5547.0818
United Kingdom	01293.422.911	01293.422.922
United States	888.387.2381	800.718.2056

Technical Services

BD Technical Services

<http://www.bd.com/vacutainer>
 Vacutainer_TechServices@bd.com
 800.631.0174 (USA Only)

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