

Helping all people
live healthy lives

BD Biosciences Stem Cell Source
— from isolation to analysis

Collagen Composite and OPLA® Scaffolds are proprietary biomaterials of Kensey Nash Corporation. Calcium Phosphate Scaffold is a proprietary biomaterial is Millenium Biologix, all manufactured for BD Biosciences. OPLA is a registered trademark of Kensey Nash Corporation. OPLA material is manufactured under one or more of the following US Patents: 4,186,448; 5,736,160; 5,855,608; 5,935,594; 5,981,825; 6,005,161; and 6,264,701. Alexa Fluor is a trademark of Molecular Probes, Inc. OPLA is a trademark of Kensey Nash Corp. PuraMatrix is a registered trademark of 3D Matrix, Inc. Calcium Phosphate Scaffolds is a proprietary biomaterial of Millenium Biologix Corp., Ontario, Canada, manufactured for BD Biosciences. BD flow cytometers are class I (1) laser products. ©2006 Becton, Dickinson and Company. All rights reserved. No part of this publication may be reproduced, transmitted, transcribed, stored in retrieval systems, or translated into any language or computer language, in any form or by any means: electronic, mechanical, magnetic, optical, chemical, manual, or otherwise, without prior written permission from BD Biosciences. For research use only. Not for use in diagnostic or therapeutic procedures. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton Dickinson and Company is strictly prohibited. Prices are subject to change without notice. BD, BD Logo, and all other trademarks are the property of Becton, Dickinson and Company. © 2006 BD

Table of Contents

Introduction	5
Stem Cells and their Functionality	6
Embryonic Stem Cells	6
Adult Stem Cells	6
Summary of Similarities and Differences Between Embryonic and Adult Stem Cells	6
Hematopoietic and Stromal Stem Cell Differentiation	7
Potential Uses of Human Stem Cells	8
Stem Cell Markers	8
Extracellular Matrix Molecules, Growth Factors, and other Factors used to Promote Stem/Progenitor Cell Proliferation and Differentiation <i>In Vitro</i>	13
Highlighted Products for Stem Cell Biology and Tissue Engineering Research	15
BD™ PuraMatrix™ Peptide Hydrogel	15
BD™ Three Dimensional Scaffolds	16
CD34 Human Hematopoietic Progenitor Cell Marker	17
BD Procount™ Progenitor Cell Enumeration Kit	18
BD™ Stem Cell Control Kit	18
BD FACSCalibur™ Flow Cytometry System	18
Human Hematopoietic Stem Cell Research Reagents	19
Alexa Fluor® Dyes	19
Mouse Hematopoietic Stem (HSC) Cell Research Reagents	20
Tools for Stem Cell Research	21
Instrumentation	21
Next Generation Flow Cytometry High-Speed Sorter	21
Revolutionary Optics System of the BD FACSAria™ Flow Cytometer	22
BD FACSVantage™ SE System	22
BD Biosciences Fluorochrome Reference Chart	23
Typical Instrument Configuration	23
Identification of a Primitive Progenitor Cell Subset Using Hoechst 33342	24
BD IMag™ Magnetic Cell Separation System	26
Mouse Hematopoietic Progenitor (Stem) Cell Enrichment Set - DM	26
Depletion Flow Chart for Mouse Hematopoietic Progenitor (Stem) Cells	27
Imaging – Seeing is Believing	28
BD CARV II™ Confocal Imager	28
BD Pathway™ Bioimager	28
Product Listing	30
Other Reagents for Stem Cell Research	41
References	44
Additional Resources	45
bdbiosciences.com/stemcellsource	45

BD Biosciences Stem Cell Source — from isolation to analysis

- Over 20 years of cell separation leadership
- Dependable innovation in cell culture technology
- Advanced functional assays
- Thousands of proven reagents and systems for cell analysis



Introduction

Stem cell research is a rapidly expanding area of investigation, with the ultimate goal to prevent, diagnose, and treat human diseases, including heart disease, diabetes, cancer, stroke, and neurological disorders, such as Parkinson's and Alzheimer's disease.

BD Biosciences has provided integrated, high-value products and services for cell culture and analysis for more than two decades with hallmark brands such as BD Falcon™ and BD FACST™. BD FACST™ brand instruments are considered the gold standard in the field of flow cytometry and, along with BD Pharmingen™ and BD Transduction Laboratories™ brand antibodies, we offer one of the largest collections of reagents for applications including flow cytometric analysis, western blotting, immunofluorescence microscopy, and ELISA.

With the addition of new instrument platforms such as the BD Pathway™ high-content imager and BD CARV II™ confocal imager, we now provide one of the most comprehensive portfolios of cell culture and cell analysis tools to the research community. BD Biosciences - Discovery Labware and its BD Falcon™, BD BioCoat™, and BD™ brands offer tools for researchers investigating the mechanisms that dictate stem cell growth and differentiation, as well as stem cell use for tissue engineering applications.

In stem cell research, BD Biosciences was first to commercialize systems and reagents for characterizing and isolating hematopoietic stem cells for research use. This, fueled by research developments, has led to the now well-accepted therapeutic approaches to leukemia and lymphoma. It is one of our continuous goals to build our successes and maintain our leadership as a provider of research and diagnostic tools for somatic and embryonic stem cell research.

This brochure provides an overview of our current product portfolio for stem cell investigation. Part of the BD Biosciences core focus has been to foster collaborations to develop or license new reagents and technologies. If you would like to contribute to our expanding portfolio of stem cell research products and services, please contact one of our BD Biosciences Marketing or Business Development Departments.

Stem Cells and Their Functionality

Stem cells are unspecialized cells that have the remarkable ability to differentiate into the many cell types required for tissue and organ development.

Stem cells are known to proliferate over very long periods. An important area of research is the study of signals in an organism that cause a stem cell to proliferate and remain unspecialized until differentiated cells are needed for the growth or repair of a specific tissue. Understanding this regulation of cell differentiation caused by internal or external factors is critical to ultimately developing therapies to fight a multitude of diseases.

Embryonic Stem Cells

Embryonic stem cells (ESCs) are pluripotent stem cells derived from the embryo blastocyst. They renew themselves indefinitely and develop into different mature cell types, like the brain, heart, bones, muscles, and skin cells. Embryonic stem cells are the early building blocks of each human being.

Adult Stem Cells

Adult stem cells, somatic stem cells, or tissue-specific progenitor stem cells, are undifferentiated cells found in tissues or organs. Their role is to maintain and repair the tissues in which they exist. Adult stem cells can only proliferate for a limited number of cycles. Their ability to differentiate is inversely correlated to the number of proliferation cycles.

Recently, researchers have determined that adult stem cells are present in tissues previously thought to be devoid of these cells. These tissues include brain, bone marrow, peripheral blood, blood vessels, skeletal muscle, skin, adipose, and liver.

Three kinds of stem cells are in adult bone marrow. Hematopoietic stem cells (HSCs), form all blood cells in the body. Bone marrow stromal stem cells, or mesenchymal stem cells (MSCs), generate the cells of bone, cartilage, fat, and fibrous connective tissue.

Central nervous system (CNS) stem cells are capable of self-renewal and give rise to the three main cell types of the CNS: neurons, astrocytes, and oligodendrocytes.


Embryonic and adult stem cells have advantages and disadvantages regarding their potential use for cell-based regenerative therapies:

- Embryonic stem cells can become any cell type in the body, while adult stem cells are generally limited to the different cell types within the tissue or organ.
- Embryonic stem cells can be grown in culture relatively easily, while adult stem cells are rare and difficult to expand for stem cell therapies. Exceptions are MSCs and hepatocyte progenitor cells that are relatively easy to expand.
- An advantage to adult stem cells is that the patient's own cells can be used, grown in cell culture and then reintroduced, thus avoiding the difficulties caused by immune rejection. A tremendous amount of research is focused on the optimization of methods for growing and differentiating adult stem cells.

Summary of Similarities and Differences Between Embryonic and Adult Stem Cells

- Embryonic stem cells are pluripotent cells that can differentiate into all cell types of the body.
- Embryonic stem cells can reproduce themselves indefinitely.
- Adult stem cells have limited potential for reproduction and differentiation.
- Embryonic stem cells can be relatively easily grown in cell culture, while adult stem cells are more rare, and methods for growing them in culture are more complicated. This limits their use since a large number of cells are needed for stem cell therapies.
- The advantage to adult stem cells is that the patient's own cells can be expanded in culture and reintroduced to the patient. The cells would not be rejected by the immune system.

Hematopoietic and Stromal Stem Cell Differentiation



Stem Cells and Their Functionality *(continued)*

Potential Uses of Human Stem Cells

- **Basic studies of the role of signal transduction and gene expression in the regulation of stem cell growth and differentiation**

Conditions such as cancer and birth defects are due to abnormal cell division and differentiation. An understanding of the mechanisms that dictate these conditions will be required to develop strategies for therapy.

- **For use in drug development**

Stem and progenitor cells may enable the identification of new drug targets for a variety of indications.

- **For clinical use**

Stem cells have great promise for use in cell therapy and regenerative medicine. For example, stem cell transplantation into the heart may assist functional tissue regeneration in patients that suffer from heart disease.

While embryonic and adult stem cells may ultimately be used for cell and tissue engineering-based therapies in humans, it is well recognized that considerable basic research will be required to achieve this goal. This work will include both *in vitro* and *in vivo* analyses of stem cell growth and differentiation.

Stem Cell Markers


Molecular markers are critical to the study of cellular differentiation and cell-fate specification. To study the transition from one cell type to another, a researcher must conclusively identify distinct cell types.

While stem cells are best defined functionally, a number of molecular markers have been used to characterize various stem cell populations. In many cases, a combination of multiple markers is used to identify a particular stem cell type.

For example, CNS stem cells express a range of markers, including nestin, vimentin, Sox2, Ki-67, GD2 ganglioside, and CD15 (Lewis X), as well as the tetraspanins CD9 and CD81, CD95 (Fas), and major histocompatibility complex (MHC) class I antigens.

A convenient approach for using markers as a research tool is flow cytometry and cell sorting. A cell sorter can separate the rare stem cells from the millions of other cells present in a tissue by selecting for more than one type of marker at one time. These methods also help determine the purity of separated cells.

Although a number of tissue engineering-based therapies have been developed, a major obstacle is availability of suitable cells. While undifferentiated stem and progenitor cells are a promising source for key differentiated cell types, there is no method for expanding many undifferentiated cell populations. Ultimately, we will learn how to control the proliferation and differentiation of cells appropriately *in vitro* and *in vivo*, by controlling their microenvironment. This knowledge could enable tissue engineering to emerge as a viable therapeutic option for the treatment of many diseases.



Stem Cell Markers

Marker Name	Cell Type	Localization	Significance	Hu	Ms
Blood Vessel					
ABCG2 (CDw338)	Primitive stem cells	Cell surface	A multi-drug resistance protein that is a member of the ATP binding cassette (ABC) transporters. It is highly expressed on primitive stem cells as identified by the "side-population" (SP) phenotype.	√	
CD144 Vascular endothelial cell cadherin (VE-Cadherin)	Smooth muscle	Cell surface	Identifies smooth muscle cells in the wall of blood vessels.	√	√
Fetal liver kinase-1 (Flk1)	Endothelial	Cell surface	Cell-surface receptor protein that identifies endothelial cell progenitors; marker of cell-cell contacts.		√
Smooth muscle cell-specific myosin heavy chain	Smooth muscle	Intracellular	Identifies smooth muscle cells in the wall of blood vessels.		
Bone Marrow (BM) and Blood					
Bone morphogenetic protein receptor (BMPR)	MSC and progenitor cells	Cell surface	Important for the differentiation of committed mesenchymal cell types from mesenchymal stem and progenitor cells. BMPR identifies early mesenchymal lineages (stem and progenitor cells).	√	√
Bone Morphogenic Protein (BMP-2 and BMP-4)	ESC, endoderm	Secreted	BMP plays an important role in induction of mesoderm, hematopoiesis, and epidermis formation.		
CD3	Leukocyte lineage	Cell surface	Specific for developing and mature T lymphocytes.	√	√
CD4	Leukocyte lineage	Cell surface	Specific for developing and mature T lymphocytes.	√	√
CD8	Leukocyte lineage	Cell surface	Specific for developing and mature T lymphocytes.	√	√
CD19	Absent on HSC, present on leukocyte lineage	Cell surface	Specific for developing and mature B lymphocytes.	√	√
CD34	HSC, satellite, endothelial progenitor	Cell surface	Indicative of a HSC or endothelial progenitor; CD34 also identifies muscle satellite, a muscle stem cell.	√	√
CD38	Absent on HSC, present on leukocyte lineage	Cell surface	Cell-surface molecule that identifies leukocyte lineages. Selection of CD34 ⁺ /CD38 ⁻ cells enables purification of HSC populations.	√	√
CD44	HSC, MSC	Cell surface	Cell-adhesion molecule expressed by HSC, MSC, most developing and mature leukocytes, and many mesenchymal cell types.	√	√
CD45 (Leukocyte common antigen)	Absent on HSC, present on leukocyte lineage	Cell surface	Cell-surface protein on leukocytes and their progenitors.	√	√
CD45R/B220	Absent on HSC, present on leukocyte lineage	Cell surface	Specific for developing and mature B lymphocytes.		√
CD59	HSC	Cell surface	Glycoprotein expressed on hematopoietic and non-hematopoietic cells. CD59 also participates in spontaneous T-cell/erythrocyte adhesion, interacts with CD2, and plays a role in T-cell activation.	√	√
CD90 (Thy-1)	HSC, MSC	Cell surface	Cell-surface protein; negative or low detection is suggestive of HSC.	√	√
CD90.1 (Thy-1.1)	HSC	Cell surface	Cell-surface protein; negative or low detection is suggestive of HSC.		√
CD90.2 (Thy-1.2)	HSC	Cell surface	Cell-surface protein; negative or low detection is suggestive of HSC.		√
CD117 (c-Kit)	HSC, MSC	Cell surface	Cell-surface receptor that identifies HSC and MSC in BM. Binding by fetal calf serum (FCS) enhances proliferation of ES, HSCs, MSCs, and hematopoietic progenitor cells.	√	√
CD135 (Flt-3)	HSC, progenitor cells	Cell surface	Receptor protein tyrosine kinase that regulates hematopoiesis.	√	√
CDCP1	HSC, progenitor cells	Cell surface	Expressed on a subset of CD34 ⁺ stem/progenitor cells in the bone marrow, cord blood and mobilized peripheral blood. Studies describe the detection of CDCP1 on CD34 ⁺ /CD38 ⁻ bone marrow stem/progenitor cells, but not on mature peripheral blood leukocytes. CDCP1 is expressed on lung, colon and breast malignant cells, and is coexpressed with CD34 and CD133 on myeloid leukemic cells.	√	
Hoechst dye	Absent on HSC	Intracellular	Fluorescent dye that binds DNA; HSCs extrude the dye and stain lightly compared with other cell types.	√	√
Lineage surface antigen (Lin)	Absent on HSC, MSC, present on differentiated erythrocytes and leukocyte lineage	Cell surface	Markers of mature blood cell lineages; detection of Lin-negative cells assist in the purification of HSC and hematopoietic progenitor populations.	√	√
Ly-6G and Ly-6C (Gr-1)	Absent on HSC, MSC, present on differentiated leukocyte lineage	Cell surface	Specific for developing and mature granulocytes.		√

Key: **BM** – Bone Marrow; **EB** – Embryoid Bodies; **MSC** – Mesenchymal Stem Cell; **ECC** – Embryonal Carcinoma Cell;

ESC – Embryonic Stem Cell; **HSC** – Hematopoietic Stem Cell; **Hu** – Human; **Ms** – Mouse; √ – Available from BD Biosciences

Stem Cells and Their Functionality (continued)

Stem Cell Markers (continued)

Marker Name	Cell Type	Localization	Significance	Hu	Ms
Ly-6A/E (Sca-1)	HSC, MSC	Cell surface	Indicative of HSC and MSC.		√
Mac-1 (CD11b/CD18)	Leukocyte lineage	Cell surface	Specific for mature granulocyte and macrophage.	√	√
Muc-18 (CD146)	Bone marrow fibroblasts, endothelial cells, MSC	Cell surface	May be important in hematopoiesis.	√	
PU.1	HSC	Cell surface	Transcription factor that is exclusively expressed on cells with hematopoietic lineage.	√	
Stem cell factor (SCF, mast cell growth factor, or c-Kit ligand)	Stromal cells, fibroblasts	Cell surface, Secreted	The two forms of the protein enhance proliferation of ES and EC and bind the receptor CD117.	√	√
Stro-1 antigen	MSC	Cell surface	Selection of Stro-1+ cells assist in isolating MSC.		
TER-119	Erythrocytes	Cell surface	Expressed on developing and mature erythrocytes.		√
Liver					
Albumin	Hepatocyte	Secreted	Principal protein produced by the liver; indicates functioning of maturing and fully differentiated hepatocytes.	√	√
CD29 (β-1 integrin)	Hepatocyte	Cell surface	Cell-adhesion molecule important in cell-cell interactions; marker expressed during development of liver.	√	√
Nervous System					
CD56 (N-CAM)	Ectoderm	Cell surface	Cell-surface molecule that promotes cell-cell interaction; indicates primitive neuroectoderm formation.	√	√
CD133 (AC133)	CNS stem cell, HSC	Cell surface	Marker for pluripotent stem cells. The CD133 positive fractions of human bone marrow, cord blood and peripheral blood have been shown to efficiently engraft in xenotransplantation models, and have been shown to contain the majority of the granulocyte/macrophage precursors, NOD/SCID repopulating cells and CD34 ⁺ dendritic cell precursors.		
CNPase	NSC, Oligodendrocyte	Intracellular	Expressed at high levels by oligodendrocytes in the central nervous system and by Schwann cells in the peripheral nervous system and is virtually absent in other cell types.		
FOXD3 (Genesis)	ESC	Intracellular	The WH gene FOXD3 (previous Hfh2 or Genesis) is expressed in the early preimplantation embryo, in ES cells, and then later in the neural crest.		
GABA	NSC	Secreted	Release of the neurotransmitter GABA by adult neuronal precursor cells that develop into neurons. Limits stem cell proliferation.		
GAD65	Neurons, oligodendrocyte, astrocyte	Intracellular	GAD65 is the predominant form of Glutamic Acid Decarboxylase found in human islets and has been shown to be the major target of antibodies in human type 1 diabetes.	√	√
GAP-43	NSC	Cell surface	GAP43 is expressed by developing and regenerating neurons, and to a lesser extent, reactive glial cells.		
Glial fibrillary acidic protein (GFAP)	Ectoderm, Astrocyte	Cell surface	Protein specifically produced by astrocytes.	√	√
HLXB9 (homeo box HB9)	Ectoderm, Mesoderm	Intracellular	Mutation in the HLXB9 transcription factor causes an autosomal dominant form of sacral agenesis.		
Microtubule-associated protein-2 (MAP-2)	Neuron	Cytoplasm	Dendrite-specific MAP; protein found specifically in dendritic branching of neuron.	√	
MSI1	Ectoderm	Intracellular	MS/1 Musashi homolog 1. Preferentially expressed in neural progenitors capable of generating both neurons and glia during embryonic and postnatal development of the CNS..		
Myelin basic protein (MBP)	Oligodendrocyte	Cell surface	Protein produced by mature oligodendrocytes; located in the myelin sheath surrounding neuronal structures.		√
Nestin	Neural progenitor	Intracellular	Intermediate filament structural protein expressed in primitive neural tissue.	√	√
Neural tubulin	Ectoderm, Neuron	Intracellular	Important structural protein for neuron; identifies differentiated neuron.		
NeuroD1	Ectoderm, Neuron	Intracellular	Important in initiating neuronal differentiation and maintaining the nervous system.		
Neurofilament (NF)	Neuron	Intracellular	Important structural protein for neuron; identifies differentiated neuron.	√	√
Neuron Specific Enolase (NSE)	Neuron	Intracellular	Useful marker for neuronal lineage.	√	
Neuron Specific Nuclear Protein (NeuN)	NSC	Intracellular	NeuN is a neuron-specific, DNA-binding nuclear protein, which makes it an important marker for neurons in vertebrates.		
P75 Neurotrophin R	NSC	Cell surface	Suggested to serve critical functions during the development of the nervous system.		
Noggin	Neuron	Intracellular	A neuron-specific gene expressed during the development of neurons.		

Key: **BM** – Bone Marrow; **EB** – Embryoid Bodies; **MSC** – Mesenchymal Stem Cell; **ECC** – Embryonal Carcinoma Cell;

ESC – Embryonic Stem Cell; **HSC** – Hematopoietic Stem Cell; **Hu** – Human; **Ms** – Mouse; √ – Available from BD Biosciences

Stem Cell Markers *(continued)*

Marker Name	Cell Type	Localization	Significance	Hu	Ms
O4	Oligodendrocyte	Cell surface	O4 is present on immature, developing oligodendrocyte.		
O1	Oligodendrocyte	Cell surface	O1 characterizes mature oligodendrocyte.		
Pax-5	Ectoderm	Intracellular	A B cell-specific activator protein (BSAP) is present in the early stages of B lymphocyte development.	√	√
Pax-6	Ectoderm, Neuron	Intracellular	Human homologues of the <i>Drosophila melanogaster</i> gene <i>prd. Pax-6</i> is expressed in the developing nervous system and in developing eyes.		
PSA-NCAM (Polysialic acid-neural cell adhesion molecule)	NSC	Cell surface	An embryonic form of N-CAM mainly expressed in the developing nervous system.		
S100	NSC, Astrocyte	Intracellular	Involved in microtubule dissociation and inhibition of microtubule assembly		
SOX1	Ectoderm	Intracellular	SOX1 is essential for ventral telencephalic development.		
Synaptophysin	Neuron	Cell surface	Neuronal protein located in synapses; indicates connections between neurons.		√
Tau	Neuron	Intracellular	Type of MAP; helps maintain structure of the axon.		
Tyrosine Hydroxylase (TH)	Ectoderm, Neuronal lineage	Intracellular	Rate-limiting enzyme in the synthesis of catecholamines, tyrosine hydroxylase has a key role in the physiology of adrenergic neurons.		
Wnt 1	NSC, Endoderm	Secreted	Signaling molecule that helps control the differentiation of neurons that release the transmitter dopamine and serotonin		
Pancreas					
Cytokeratin 19 (CK19)	Pancreatic epithelium	Intracellular	CK19 identifies specific pancreatic epithelial cells that are progenitors for islet cells and ductal cells.	√	
Glucagon	Pancreatic islet, Endoderm	Intracellular	Expressed by α -islet cell of pancreas.		
Insulin	Pancreatic islet, Endoderm	Intracellular	Expressed by β -islet cell of pancreas.	√	√
Insulin-promoting factor-1 (PDX-1)	Pancreatic islet, Endoderm	Intracellular	Transcription factor expressed by beta-islet cell of pancreas.		
Nestin	Pancreatic progenitor	Intracellular	Structural filament protein indicative of progenitor cell lines including pancreatic.	√	√
Pancreatic polypeptide	Pancreatic islet	Secreted	Expressed by γ -islet cell of pancreas.		
Somatostatin	Pancreatic islet	Secreted	Expressed by δ -islet cell of pancreas.		
STAT3	ESC	Intracellular	STAT family members are phosphorylated by the receptor associated kinases, and then form homo- or heterodimers that translocate to the cell nucleus where they act as transcription activators.		
TAT-SF1	ESC	Intracellular	General transcription elongation factor.		
TERT	ESC	Intracellular	Telomerase expression plays a role in cellular senescence, as it is normally repressed in postnatal somatic cells resulting in progressive shortening of telomeres.		
Pluripotent Stem Cells					
Alkaline phosphatase	ESC, ECC	Cell surface	Elevated expression of this enzyme is associated with undifferentiated pluripotent stem cell (PSC).	√	√
Alpha-fetoprotein (AFP)	Endoderm	Intracellular	Protein expressed during development of primitive endoderm; reflects endodermal differentiation.		
Bone morphogenetic protein-4	Mesoderm	Secreted	Growth and differentiation factor expressed during early mesoderm formation and differentiation.		
Brachyury	Mesoderm	Intracellular	Transcription factor important in the earliest phases of mesoderm formation and differentiation; used as the earliest indicator of mesoderm formation.		
CD9	ESC	Cell surface	CD9 is associated with leukemia inhibitory factor (LIF)-mediated maintenance of ES.	√	√
CD30	ESC, ECC	Cell surface	Surface receptor molecule found specifically on PSC.	√	√
CD56 (N-CAM)	Ectoderm	Cell surface	Cell-surface molecule that promotes cell-cell interaction; indicates primitive neuroectoderm formation.	√	√
CRABP2	Ectoderm, Neuronal Lineage	Intracellular	Function unknown. The inducibility of the CRABP2 gene suggests that this isoform is important in retinoic acid-mediated regulation of human skin growth and differentiation.		
Cripto (TDGF-1)	ESC, Cardiomyocyte	Cell surface	Gene for growth factor expressed by ES, primitive ectoderm, and developing cardiomyocyte.		
CRTR-1	ESC	Intracellular	CP2 related transcriptional repressor 1 (CRTR 1), is expressed in a developmentally regulated fashion <i>in vivo</i> and acts as a direct repressor of transcription.		
GATA-4 gene	Endoderm	Intracellular	Expression increases as ES differentiates into endoderm.		

Key: **BM** – Bone Marrow; **EB** – Embryoid Bodies; **MSC** – Mesenchymal Stem Cell; **ECC** – Embryonal Carcinoma Cell;

ESC – Embryonic Stem Cell; **HSC** – Hematopoietic Stem Cell; **Hu** – Human; **Ms** – Mouse; √ – Available from BD Biosciences

Stem Cells and Their Functionality *(continued)*

Stem Cell Markers *(continued)*

Marker Name	Cell Type	Localization	Significance	Hu	Ms
GCTM-2	ESC, ECC	Secreted	Synthesized by undifferentiated PSCs.		
Genesis	ESC, ECC	Intracellular	Transcription factor uniquely expressed by ES cells either in or during the undifferentiated state of PSC.		
Germ cell nuclear factor	ESC, ECC	Intracellular	Transcription factor expressed by PSC.		
Nanog	ESC	Intracellular	Recognizes the undifferentiated state of human and mouse stem cells. Nanog directs propagation of undifferentiated ES cells.		
Hepatocyte nuclear factor-4 (HNF-4)	Endoderm	Intracellular	Transcription factor expressed early in endoderm formation.		
Nucleostemin	ESC	Intracellular	Plays a role in tumor stem cell biology. Although the exact function is not yet known, it behaves like a molecular switch to control cell division.		
HLXB9 (homeo box HB9)	Ectoderm, Mesoderm	Intracellular	Mutation in the HLXB9 transcription factor causes an autosomal dominant form of sacral agenesis.		
Hsp27	ESC	Intracellular	Hsp27 has been identified as a switch between differentiation and apoptosis in murine embryonic stem cells.		
Nestin	Ectoderm, neural and pancreatic progenitor	Intracellular	Intermediate filaments within cells; characteristic of primitive neuroectoderm formation.	√	√
Oct-3	ESC, ECC	Intracellular	Member of the POU family of transcription factors; functions in pluripotent cells of early embryos and undifferentiated embryonal cell lines.		√
Podocalyxin	HSC, ESC	Cell surface	CD34-related marker of murine HSCs and ESCs.		
Pramel-4	ESC	Intracellular	Target gene of Oct-4 that is expressed in pluripotent stem cells.		
Pramel-5	ESC	Intracellular	Target gene of Oct-4 that is expressed in pluripotent stem cells.		
Rex-1	ESC	Intracellular	A marker gene for undifferentiated ES cells.		
Silylated keratan sulfate proteoglycan	ESC, ECC	Cell surface	Specific extracellular matrix-associated molecule is synthesized by undifferentiated PSC.		
Sox-2	ESC	Intracellular	SOX2 belongs to the SOX (SRY-like HMG box) family of transcription factors with diverse roles in development.		
Stage-specific embryonic antigen-1 (SSEA-1; CD15)	ESC, ECC	Cell surface	Glycoprotein specifically expressed in early embryonic development and by undifferentiated PSCs.	√	√
Stage-specific embryonic antigen-3 (SSEA-3)	ESC, ECC	Cell surface	Glycoprotein specifically expressed in early embryonic development and by undifferentiated PSCs.		
Stage-specific embryonic antigen-4 (SSEA-4)	ESC, ECC	Cell surface	Glycoprotein specifically expressed in early embryonic development and by undifferentiated PSCs.		
TAZ	MSC	Intracellular	Transcriptional coactivator with PPZ-binding motif that functions as a transcriptional modulator of MSC differentiation.		
Telomerase	ESC, ECC	Intracellular	An enzyme uniquely associated with immortal cell lines; useful for identifying undifferentiated PSCs.	√	√
Telomerase Reverse Transcriptase (TERT)	ESC	Intracellular	TERT plays a role in senescence, as it is normally repressed in postnatal somatic cells resulting in progressive shortening of telomeres.		
TRA-1-60	ESC	Cell surface	Pluripotency stem cell marker.		
TRA-1-81	ESC	Cell surface	Pluripotency stem cell marker.		
UTF1	ESC	Intracellular	Expressed mainly in ES cells. Oct-3/4 induces rapid proliferation and tumorigenic properties of ES cells through activation of the UTF1 gene.		
Vasa-1	ESC, germ cells	Intracellular	The human VASA gene is specifically expressed in late germ cell lineage development.		
Vimentin	Ectoderm, neural and pancreatic progenitor	Intracellular	Intermediate filaments within cells; characteristic of primitive neuroectoderm formation.	√	√
Skeletal Muscle/Cardiac/Smooth Muscle					
CD140α	NSC	Cell surface	Involved in cell proliferation, differentiation, and survival.		
HLA-ABC	All nucleated cells	Cell surface	Receptor for CD8/TCR, T/APC interaction.	√	
MyoD and Pax7	Myoblast, myocyte	Intracellular	Transcription factors that direct differentiation of myoblasts into mature myocytes.	√	√
Myogenin and MR4	Skeletal myocyte	Intracellular	Secondary transcription factors required for differentiation of myoblasts from muscle stem cells.	√	√
Myosin heavy chain	Cardiomyocyte	Intracellular	A component of structural and contractile protein found in cardio myocyte.		
Myosin light chain	Skeletal myocyte	Intracellular	A component of structural and contractile protein found in skeletal myocyte.		

Key: **BM** – Bone Marrow; **EB** – Embryoid Bodies; **MSC** – Mesenchymal Stem Cell; **ECC** – Embryonal Carcinoma Cell;

ESC – Embryonic Stem Cell; **HSC** – Hematopoietic Stem Cell; **Hu** – Human; **Ms** – Mouse; √ – Available from BD Biosciences

To investigate the biological properties of stem cells and tissue-specific progenitor cells *in vitro*, cell culture systems have been established to characterize the effects of key extracellular matrix (ECM) molecules, growth factors/cytokines, and other factors that influence stem/progenitor cell expansion and differentiation. For tissue engineering applications, stem cells have been used for *in vivo* studies to assess the production and functionality of key differentiated cell types in the context of tissue regeneration and repair.

Extracellular Matrix (ECM) Molecules, Growth Factors, and other Factors Used to Promote Stem/Progenitor Cell Proliferation and Differentiation *In Vitro*

Stem Cell Type	Progenitor or Differentiated Cell Type	ECM, Growth Factor/Cytokine, and/or Other Factor
Embryo		
Embryonic stem cell	Expansion of undifferentiated ESC ^{8,9,41-43}	BD Matrigel™ Matrix, BD™ PuraMatrix™ Peptide Hydrogel, Gelatin, LIF
	Pancreatic endocrine progenitor ⁸	bFGF
	Pancreatic islet ⁸	bFGF, nicotinamide
	Expansion of Nestin+ neural progenitors ⁹	Laminin, bFGF, poly-ornithine
	Neuron ^{9,10}	Laminin, poly-ornithine, RA
	Survival of ES-derived motor neurons ¹¹	BD Matrigel™ Matrix, BDNF, NT-3, CNTF, GDNF
	Glial progenitor cells ¹²	bFGF, PDGF-AA
	Adipocyte ¹³	RA
	Chondrocyte ¹⁴	BMP-2, BMP-4
	Osteoblast ¹⁵	RA, dexamethasone, ascorbate, β-glycerol phosphate
	Dendritic cells ¹⁶	GM-CSF, IL-3
Endothelial cells ¹⁷	Collagen IV, VEGF	
Bone Marrow		
Mesenchymal stem cell	Osteoblast ¹⁸⁻²⁰	BMP-2, bFGF, dexamethasone, ascorbate, β-glycerol phosphate
	Chondrocyte ^{19,21}	TGF-β3, dexamethasone
	Neuron, Glial cells ²²	EGF, BDNF, RA
	Adipocyte ^{19,20}	Dexamethasone, insulin, indomethacin, 1-methyl-3-isobutylxanthine
Multipotent adult progenitor cell	Expansion of undifferentiated MAPCs ²³	Fibronectin, EGF, PDGF-BB, 2% FBS
	Osteoblast ²³	Fibronectin, dexamethasone, ascorbate, β-glycerol phosphate
	Chondrocyte ²³	TGF-β1
	Endothelial cells ²⁴	Fibronectin, VEGF
	Hepatocyte-like cells ²⁵	BD Matrigel™ Matrix, FGF-4, HGF
Hematopoietic stem cell	Platelets, red/white blood cells ²⁶	IL-3, IL-6, G-CSF, erythropoietin, thrombopoietin
	Expansion of undifferentiated progenitors	IL-3, IL-6, Stem cell factor, Flt3 ligand, thrombopoietin
	Megakaryocytes, Platelets	IL-3, IL-6, GM-CSF, thrombopoietin
	Reticulocytes, Erythrocytes	IL-3, IL-6, IL-9, erythropoietin, GM-CSF
	B lymphocytes	CXCL12, IL-3, Flt3 ligand, IL-7
	Granulocytes	G-CSF, GM-CSF, IL-3, IL-6
	Monocytes, Macrophages, Dendritic cells	GM-CSF, M-CSF, TNF-α
	Mast Cells	IL-4, IL-6

Key: *BM* – Bone Marrow; *EB* – Embryoid Bodies; *MSC* – Mesenchymal Stem Cells; *ECC* – Embryonal Carcinoma Cells;

ESC – Embryonic Stem Cells; *HSC* – Hematopoietic Stem Cells; *Hu* – Human; *Ms* – Mouse; √ – Available from BD Biosciences

Color Key: ECM molecules are denoted in blue; growth factors/cytokines in red; other factors in black.


Stem Cells and Their Functionality (continued)

Extracellular Matrix (ECM) Molecules, Growth Factors, and Other Factors Used to Promote Stem/Progenitor Cell Proliferation and Differentiation *In Vitro* (continued)

Stem Cell Type	Progenitor or Differentiated Cell Type	ECM, Growth Factor/Cytokine, and/or Other Factor
Thymus		
Multipotential thymic precursor	T lymphocytes	IL-7
	NK cells	Flt3 ligand, IL-7, IL-15, IL-2, IL-12, IL-18
	Dendritic cells	GM-CSF, IL-3, IL-6, Flt3 ligand
Adult Tissue		
Liver	Expansion of undifferentiated oval cells ²⁷	Fibronectin, SCF, Flt-3 ligand, IL-3, LIF
	Hepatocyte ²⁸	Fibronectin, HGF, EGF
	Pancreatic islet ^{27,28}	Fibronectin, SCF, Flt-3 ligand, IL-3, nicotinamide
	Hepatocyte Progenitor ³⁹	BD™ PuraMatrix™ Peptide Hydrogel
Neural stem cell/ Neural progenitor cell	Expansion of NPCs ³⁹	βFGF, EGF, LIF
	Neuron, Glial cells ³⁰⁻³⁴	Fibronectin, Laminin, βFGF, EGF, PDGF-AA, PDGF-AB, PDGF-BB, NT-4, CNTF
	Neural, Glial Progenitors ⁴⁰	BD™ PuraMatrix™ Peptide Hydrogel
Adipose stem cell	Osteoblast ³⁵	Dexamethasone, ascorbate, β-glycerol phosphate
	Chondrocyte ^{35,36}	TGF-β1, insulin, ascorbate
	Adipocyte ³⁵	Dexamethasone, insulin, indomethacin, 1-methyl-3-isobutylxanthine
Muscle-derived stem cell	Expansion of MDSCs ³⁷	Collagen, IGF-1, EGF, SCF, FGF2
	Osteoblast ³⁸	Collagen, BMP-2

Key: ECM – extracellular matrix; EGF – epidermal growth factor; LIF – leukemia inhibitory factor; SCF – stem cell factor; HGF – hepatocyte growth factor; PDGF – platelet-derived growth factor; VEGF – vascular endothelial growth factor; BMP – bone morphogenetic protein; BDNF – brain derived neurotrophic factor; NT – neurotrophin; CNTF – ciliary neurotrophic factor; bFGF – basic fibroblast growth factor; TGF-β – transforming growth factor-beta; IL – interleukin; G-CSF – granulocyte-colony stimulating factor; GM-CSF – granulocyte-macrophage colony stimulating factor; IGF – insulin-like growth factor; RA – retinoic acid; FBS – fetal bovine serum.

Color Key: ECM molecules are denoted in blue; growth factors/cytokines in red; other factors in black.



Highlighted Products for Stem Cell Biology and Tissue Engineering Research

BD™ PuraMatrix™ Peptide Hydrogel

Cat. No. 354250


BD™ PuraMatrix™ Peptide Hydrogel is a synthetic matrix used to create defined, three-dimensional (3D) microenvironments for a variety of cell culture experiments. To achieve optimal cell growth and differentiation, it is necessary to determine the appropriate mixture of this material and bioactive molecules (eg, growth factors, extracellular matrix [ECM] proteins, and/or other molecules). BD PuraMatrix Peptide Hydrogel consists of standard amino acids (1% w/v) and 99% water. Under physiological conditions, the peptide component self-assembles into a 3D hydrogel that exhibits a nanometer scale fibrous structure. The hydrogel is readily formed in a culture dish, plate, or cell culture insert.

The resulting hydrogel has been shown to promote the differentiation of hepatocyte progenitor cells,³⁹ rat pheochromocytoma cells (PC12),⁴⁵ and hippocampal neurons.⁴⁰⁻⁴⁴ Studies have also demonstrated that BD PuraMatrix Peptide Hydrogel supports the attachment of a variety of primary (eg, neuronal, fibroblast, keratinocyte) and transformed (eg, MG-63, SH-SY5Y, HEK293, NIH3T3) cell types. Other potential applications include stem cell proliferation, tumor cell migration and invasion, angiogenesis assays, and *in vivo* model analyses of tissue regeneration. BD PuraMatrix Peptide Hydrogel is biocompatible, resorbable, and devoid of animal-derived material and pathogens. For *in vivo* studies in animals, the soluble material can be injected and will subsequently form a 3D hydrogel upon contact with the physiological environment.

To examine the potential utility of self-assembling peptides in the context of cardiac regeneration, recent work from the laboratory of Dr. Richard Lee demonstrated that endothelial cells improve cardiac myocyte survival and organization when these cell types were co-cultured in 3D using a peptide hydrogel (RAD16-II) closely related to BD PuraMatrix (RAD16-I).⁴⁵ In a subsequent study from the same group, RAD16-II injected into the myocardium of mice was found to promote vascular cell recruitment (endothelial and smooth muscle cells) into the microenvironment created by the hydrogel.⁴⁶

Differentiation of Hepatocyte Progenitor Cells

Rat hepatocyte progenitor cells (Lig-8) were encapsulated in BD PuraMatrix Peptide Hydrogel and cultured overnight in defined medium at 37°C. Samples were then used for bromodeoxyuridine (BrdU) uptake and *in situ* immunofluorescence analyses. As shown in the figure below, Lig-8 cells form spheroid colonies when cultured within the 3D hydrogel and express the hepatocyte markers CCAAT/enhancer binding protein (C/EBP) and cytochrome P450 1A1/1A2 (CYP1A1/1A2) in a manner that is independent of cellular mitotic activity.³⁹ Therefore, while some cells are proliferating, the entire colony exhibits differentiation potential.



Lig-8 cells cultured in BD™ PuraMatrix™ Peptide Hydrogel.

All cells in spheroid colonies, arrested or not, undergo differentiation.

Spheroids were isolated, transferred to adherent cultures, and incubated with BrdU for 24 hours.

- A. spheroid colony (phase contrast)
- B. same optical layer as A immunostained for C/EBP (red)
- C. same optical layer as A immunostained for BrdU (green)
- D. spheroid colony (phase contrast)
- E. same optical layer as D immunostained for CYP1A1/1A2 (red)
- F. same optical layer as D immunostained for BrdU

Data provided by 3D Matrix, Inc. and originally described in Semino, C.E., et al., *Differentiation* 71:262 (2003).


Features and Benefits

Purified Synthetic Peptide Composition (1% w/v)	Highly defined material that promotes cell attachment
3D Hydrogel Structure	Assembles into fibrous structure with average pore size of 50-200 nm
Easy Handling Injectable for <i>in vivo</i> Studies in Animals	Easily mixed with cells and/or bioactive molecules (eg, growth factors) prior to gelation; injectable for <i>in vivo</i> studies in animals
Transparent Hydrogel	Samples are readily visualized using standard staining methodologies and microscopy
Established Protocols	3D cell encapsulation cultures; Surface plating of adherent cells on microporous membrane inserts and microplates; cell recovery for sub-culturing or biochemical analyses; <i>in vivo</i> injection



Highlighted Products for Stem Cell Biology and Tissue Engineering Research *(continued)*

BD™ Three Dimensional Scaffolds


The BD™ Three Dimensional (3D) Scaffolds are ideal for *in vitro* and *in vivo* analyses of cell growth and differentiation. These scaffolds provide a 3D physical support matrix for *in vitro* cell culture as well as *in vivo* tissue regeneration.



BD™ 3D Calcium Phosphate Scaffold (Cat. No. 354617): Supports short- and long-term growth and differentiation of mesenchymal cells.



BD™ 3D Collagen Composite (Cat. No. 354613) (left) and BD™ OPLA® (Cat. No. 354614) (Open-Cell Poly-Lactic Acid) (right) Scaffolds: Support short- and long-term growth and differentiation of a variety of epithelial cells (eg, hepatocytes), neurons, endothelial cells, osteoblasts, chondrocytes, fibroblasts, and smooth muscle cells.



CD34 Human Hematopoietic Progenitor Cell Marker

The CD34 antigen is a single-chain transmembrane glycoprotein, MW 105 to 120 kilodaltons (kd). The antigen is associated with human hematopoietic progenitor cells and is a differentiation stage-specific leukocyte antigen. The CD34 antigen is present on immature hematopoietic precursor cells and all hematopoietic colony-forming cells in bone marrow and blood, including unipotent (CFU-GM, BFU-E) and pluripotent progenitors (CFU-GEMM, CFU-Mix, and CFU-Blast). Terminal deoxynucleotidyl transferase-positive B and T lymphoid precursors in normal bone marrow are CD34⁺. The CD34 antigen is present on early myeloid cells that express the CD33 antigen but lack the CD14 and CD15 antigens and on early erythroid cells that express the CD71 antigen and dimly express the CD45 antigen.

The CD34 antigen is also found on capillary endothelial cells and approximately 1% of human thymocytes. Normal peripheral blood lymphocytes, monocytes, granulocytes, and platelets do not express the CD34 antigen. CD34 antigen density is highest on early hematopoietic progenitor cells and decreases as cells mature. The antigen is absent on fully differentiated hematopoietic cells.

Uncommitted CD34⁺ progenitor cells are CD38⁻ and lack lineage-specific antigens such as CD71, CD33, CD10, and CD5, while CD34⁺ cells that are lineage-committed express the CD38 antigen in high density. Most CD34⁺ cells reciprocally express either the CD45RO or CD45RA antigens, with the CD45RO⁺ population being the more primitive. Approximately 60% of acute B lymphoid leukemias and acute myeloid leukemias (AMLs) and 1% to 5% of acute T lymphoid leukemias express the CD34 antigen. The antigen is not expressed on chronic lymphoid leukemias or lymphomas.

BD Biosciences offers a variety of different monoclonal antibodies against CD34. Please refer to the product list in the back of this brochure or visit us online.

Several multicolor combinations are available that include combinations of CD34 with other markers such as CD45, CD15, and HLA-DR.



Highlighted Products for Stem Cell Biology and Tissue Engineering Research *(continued)*

BD Biosciences offers products for *In Vitro* Diagnostic use, and is a leading provider of CD34 reagents for clinical flow cytometry applications.

The BD Procount™ progenitor cell enumeration system is a fully integrated clinical system for automated counting of CD34 positive cells. It combines powerful software, unique reagents, and BD Trucount™ technology to minimize variability in CD34 counts with an easy-to-use method. BD Procount software, when used with Procount reagents, identifies and enumerates CD34+ cells on a BD FACSCalibur™ flow cytometer.

BD Procount™ Progenitor Cell Enumeration Kit

Cat. No. 340498

The BD Procount™ progenitor cell enumeration kit is designed as a two-step, lyse-no-wash assay. The 25-test kit consists of a CD34 reagent, an isotype control reagent, and BD Trucount™ absolute count tubes to determine the absolute cell counts volumetrically. The BD Procount analysis strategy fulfills the need for a procedure with minimal sample handling and processing, combined with a consistent and standardized data analysis method. The BD Procount kit is intended for identifying and enumerating CD34+ cells. CD34+ cells may be present in peripheral blood and leukapheresis.

BD™ Stem Cell Control Kit

Cat. No. 340991

The BD™ Stem Cell Control Kit is intended as a complete, two-level process control for immunophenotyping and enumeration of leukocytes by flow cytometry. It is a control for antibody staining, red blood cell (RBC) lysis, instrument setup and performance, and data analysis.

Immunophenotyping by flow cytometry is a complex, multi-step process. The BD Stem Cell Control Kit is a stable control, with assigned values that can be used to monitor the immunophenotyping process for CD34+ cells.

BD FACSCalibur™ Flow Cytometry System*

The BD FACSCalibur™ system combines unique dual-laser technology, an automated sample loader option, and powerful software to provide the high throughput necessary to meet productivity requirements of clinical laboratories. The modularity and innovative technology of the BD FACSCalibur system also offers researchers the flexibility required for a variety of research applications.

- Three- or four-color fluorescence capability
- BD FACStation™ system data management system
- Flexible and modular for future upgrades

*Class I (1) Laser Product

Human Hematopoietic Stem Cell Research Reagents


APC Conjugated Mouse Anti-Human Hematopoietic Progenitor Cells (HPCs)

Cat. No. 557929

Monoclonal antibody BB9 recognizes a cell-surface glycoprotein of approximately 160 kDa, expressed on primitive hematopoietic progenitor cells. This antibody was generated using human bone marrow-derived stromal cells as immunogen. BB9 reacts with stromal cells, a minor subset of human bone marrow cells; it does not react with peripheral blood leukocytes.

Reacts with a subset of CD34⁺, CD90⁺ demonstrating low or absent expression of CD38, low retention of Rhodamine 123 and Ki-67. CD34⁺BB9⁺ cells are able to sustain hematopoiesis in pre-CFU culture stimulated with IL-3, IL-6, G-CSF and SCF. BB9 also reacts with a subset of CD34⁺ cells from mobilized peripheral blood.

The BB9 antibody is also available as a PE conjugate.



APC Hematopoietic Progenitor Cell, with isotype control.


Profile of anti-hematopoietic progenitor cell (BB9) reactivity on peripheral blood lymphocytes analyzed by flow cytometry.

CD117 (c-kit), R-Phycoerythrin (R-PE)-Conjugated Mouse Anti-Human Monoclonal Antibody

Cat. No. 555714

Monoclonal antibody YB5.B8 recognizes CD117 (c-kit), a 145 kDa cell-surface glycoprotein with tyrosine kinase activity. CD117 is present on hematopoietic progenitor cell subsets, thymocytes, mast cells, hepatocytes and histiocytes. CD117 serves as a cytokine receptor for steel factor (SLF), also known as stem cell factor (SCF) or mast cell growth factor (MGF). The interaction of c-kit and SLF is crucial to hematopoiesis, mast cell differentiation, melanogenesis, and germ cell development. The ability of YB5.B8 antibody to block the binding of c-kit ligand is still controversial.

The YB5.B8 antibody is also available as APC and PE-Cy 5 conjugate.




Log Fluorescence Intensity Profile of the erythroleukemia cell line, TF-1, analyzed by flow cytometry.

CD34, R-Phycoerythrin (R-PE)-Conjugated Mouse Anti-Human Monoclonal Antibody

Cat. No. 550619

Monoclonal antibody 563 recognizes the human form of the 105-120 kDa heavily O-glycosylated and N-glycosylated transmembrane glycoprotein, CD34, expressed on stem/progenitor hematopoietic cells, endothelial cells and some tissue. Clone 563 also cross reacts with a subset of peripheral blood mononuclear cells of both rhesus and cynomolgus macaque monkeys, but not baboon. The distribution of stemlike cells is similar to that observed with peripheral blood mononuclear cells from normal human donors.




Profile of CD34 Ab reactivity on peripheral blood (CD14-negative) stem cells from rhesus macaque (Macaca mulatta) analyzed by flow cytometry.

CDw338 (ABCG2) Mouse Anti-Human Monoclonal Antibody

Cat. No. 552823

Monoclonal antibody 5D3 reacts with ABCG2 (BCRP1), a multi-drug resistance protein that is a member of the ATP binding cassette (ABC) transporters. It is highly expressed on primitive stem cells as identified by the "side-population" (SP) phenotype. This SP phenotype is based on the efflux of fluorescent dyes such as Rhodamine 123 and Hoechst 33342. The expression of ABCG2 appears to be highly conserved as it has been identified in various species. Studies show that highly purified murine stem cells express BCRP1 mRNA and this expression declines sharply as the stem cells express CD34. The highest levels of BCRP1 mRNA expression have been seen in KDR⁺ human stem cells.



ABCG2 (BCRP1)-Purified Log Fluorescence Intensity.

Profile of anti-ABCG2 (clone 5D3) reactivity on MCF7-BCRP1 cells surface analyzed by flow cytometry. Second step staining with Cat. No. 550589.

For more information visit bdbiosciences.com

Alexa Fluor® Dyes


The Alexa Fluor® dye has become a widely used dye in stem cell research due to its brightness and photostability. BD Biosciences carries a wide selection of Alexa Fluor conjugated antibodies. Please visit our online catalog for a complete listing.

Alexa Fluor® 647 Mouse Anti-Human Cub Domain Containing Protein (CDCP1) Monoclonal Antibody

Cat. No. 558212

Monoclonal antibody CUB1 reacts with the Cub Domain Containing Protein 1 (CDCP1), a transmembrane protein described to contain three CUB domains in the extracellular region with a hexalysine stretch in the cytoplasmic region. It is expressed on a subset of CD34⁺ stem/progenitor cells in the bone marrow, cord blood and mobilized peripheral blood. Studies describe the detection of CDCP1 on CD34⁺/CD38⁻ bone marrow stem/progenitor cells, but not on mature peripheral blood leukocytes. CDCP1 is expressed on lung, colon and breast malignant cells and is coexpressed with CD34 and CD133 on myeloid leukemic cells. It may play a role of adhesion or interaction with the extracellular matrix. Reports describe transplanted human CDCP1⁺ cells into NOD/SCID mice to successfully engraft and undergo multi-lineage differentiation.

The CUB1 antibody is also available as a PE conjugate.



Alexa Fluor® CDCP1 and Isotype Control

Profile of CDCP1 (CUB1) reactivity on HT29 cells analyzed by flow cytometry.


Highlighted Products for Stem Cell Biology and Tissue Engineering Research (continued)

Mouse HSC Research Reagents

R-Phycoerythrin (R-PE)-Conjugated Rat Anti-Mouse CD135 (Flk-2/Fit3, Ly-72) Monoclonal Antibody

Cat. No. 553842

The A2F10 antibody reacts with Flk-2/Fit3 (Ly-72, CD135), a receptor protein tyrosine kinase closely related to c-kit, c-fms, and PDGF Receptor of the immunoglobulin superfamily. The Flt3 message is detected in hematopoietic stem cells and primitive progenitor cells in fetal liver, adult bone marrow, and fetal and adult thymus, as well as brain, placenta, and testis; but it is absent in more mature hematopoietic cells. In flow cytometric analysis, the A2F10 antibody recognizes Flt3-transfected Y3 cells (but not the parent cell line) and early B lymphoid lineage cells in juvenile and adult bone marrow. A role for CD135 in the regulation of hematopoiesis is suggested by the observations that soluble Flk-2/Fit3 ligand can both stimulate proliferation of stem cell-enriched fetal liver, fetal thymus, and adult bone marrow populations and enhance their responses to other growth factors *in vitro*.




Expression of CD135 on bone-marrow leukocytes. BALB/c bone-marrow cells were stained with FITC-conjugated anti-mouse CD45R/B220 and CD11b (Integrin α_M chain) monoclonal antibodies (Cat. no. 553087/553088 and 554982, respectively) and either PE-conjugated Rat IgG_{2a}, κ isotype control mAb R35-95 (Cat. no. 553930, Panel A) or PE-conjugated mAb A2F10.1 (Panel B).

APC Mouse Lineage Antibody Cocktail, with Isotype Control CD3e, CD11b, CD45R/B220, Erythroid Cells, and Ly-6G and Ly-6C,

Cat. No. 558074

The APC Mouse Lineage Antibody Cocktail has been designed to react with cells from the major hematopoietic lineages, such as T lymphocytes, B lymphocytes, monocytes/macrophages, NK cells, erythrocytes, and granulocytes. This pre-diluted Cocktail of five APC-conjugated antibodies is designed for the flow cytometric identification of hematopoietic progenitors in mouse bone marrow. The APC Mouse Lineage Isotype Control Cocktail contains equivalent concentrations of isotype-matched negative-control immunoglobulin. Additional fluorochrome-labeled reagents may be combined with the APC Mouse Lineage Antibody Cocktail, and the APC Mouse Lineage Isotype Control Cocktail, to further characterize hematopoietic progenitor subpopulations.




Identification of CD34⁺ and CD34⁻ subpopulations of hematopoietic progenitors. BALB/c bone marrow cells were treated with Mouse BD Fc Block™ purified anti-CD16/CD32 mAb 2.4G2 (Cat. No. 553141/553142) and stained with either FITC rat IgG_{2a} κ isotype control (not shown) or FITC anti-CD34 mAb RAM34 (Cat. No. 553733) and with either APC Mouse Lineage Isotype Control Cocktail (left panel) or APC Mouse Lineage Antibody Cocktail (right panel). Dead cells were excluded from analysis by staining with propidium iodide (Cat. No 556463).

R-Phycoerythrin (R-PE)-Conjugated Rat Anti-Mouse CD34 Monoclonal Antibody

Cat. No. 551387

The RAM34 antibody reacts with the CD34 glycoprotein on the surface of three independently derived mouse CD34-transfected cell lines. Normal thymocytes and splenocytes are negative for CD34 expression. In the bone marrow, 7-10% of cells are stained with RAM34 mAb, including most of the Ly-6A/E (Sca-1)⁺ CD90 (Thy-1) low Lineage Marker- hematopoietic stem cell-enriched subpopulation and myeloerythroid progenitors. CD34 is also expressed on a small percentage of fetal liver cells, including NK-cell progenitors. CD34 has been reported to be expressed on the endothelium of capillaries and, in this form, to function as a ligand for L-selectin. Consistent with this observation, RAM34 antibody stains endothelial cells in spleen, thymus, and postcapillary HEVs in the lymph nodes. It is reported that RAM34 antibody can be used to select CD34⁺ CD117 (c-Kit)⁺ Ly-6A/E (Sca-1)⁺ Lineage Marker- bone marrow-derived hematopoietic stem cells, capable of short-term multi-lineage reconstitution of lethally irradiated mice; while the CD34⁻ CD117⁺ Sca-1⁺ Lineage Marker- population contains self-renewing hematopoietic stem cells. Similarly, the bone marrow population with high dye efflux capacity and which is highly enriched for long-term reconstituting hematopoietic stem cells is CD34⁻ CD117 (c-Kit)⁺ Ly-6A/E (Sca-1)⁺ Lineage Marker.



The expression of CD34 on Lineage Marker-negative bone-marrow cells. C57BL/6 bone-marrow leukocytes were stained with PE-conjugated RAM34 mAb before (Panel A) and after (Panel B) depletion of the lineage-committed cells with the Mouse Lineage Panel No. 559971) followed by a magnetic beads separation method. The open histograms represent stained cells, while the shaded histograms represent unstained negative-control cells.

The RAM34 antibody is also available as FITC, APC, and Biotin conjugate.

Ly-6A/E (Sca-1) R-Phycoerythrin (R-PE)-Conjugated Rat Anti-Mouse Monoclonal Antibody

Cat. No. 553108

The D7 antibody reacts with Ly-6A.2 and Ly-6E.1, which are allelic members of the Ly-6 multigene family. Sca-1 (Ly-6A/E), a phosphatidylinositol-anchored protein of about 18 kDa, is expressed on the multipotent hematopoietic stem cells (HSC) in the bone marrow of mice with both Ly-6 haplotypes. In mice expressing the Ly-6.2 haplotype (eg, AKR, C57BL, C57BR, C57L, C58, DBA/2, PL, SJL, SWR, 129), Ly-6A/E is also expressed on distinct subpopulations of bone marrow and peripheral B lymphocytes and thymic and peripheral T lymphocytes.

The D7 antibody is also available as FITC, PE-Cy 7, and Biotin conjugate.

For more information visit bdbiosciences.com

Tools for Stem Cell Research

Instrumentation

Next Generation Flow Cytometry High-Speed Sorter

BD FACSAria™ Cell Sorter*


The BD FACSAria™ cell sorter sets a new standard for high-performance flow cytometry. Based on an entirely new design in instrumentation, the BD FACSAria cell sorter is the first benchtop sorter to incorporate a fixed alignment cuvette flow cell. The numerous technological advances embodied in the BD FACSAria cell sorter reduce the cost of owning a high-speed sorter and accelerate research by providing unparalleled, easy-to-use sorting and analysis.

BD FACSAria software provides pioneering instrument control for flexible sorting and data analysis. All instrument and sorting controls are embedded in the software to simplify instrument setup and operation. New features programmed into the software provide the complete solution for high-performance cell sorting.



A customized 355nm (UV) or 375nm (near UV) laser option for the BD FACSAria™ allows for precise sorting of side population stem cells labeled with Hoechst 33342. The collected cells can be cultured and used in further experiments.

*Class I (1) Laser Product




Tools for Stem Cell Research (continued)

Revolutionary Optics System of the BD FACSAria™ Flow Cytometer

High Sensitivity

This new optical system increases the quality and quantity of the information acquired from each sample.



Octagon-shaped collection optics.

For more information, visit bdbiosciences.com/FACSAria

Class I (1) Laser Product

BD FACSVantage™ SE System*

High-speed cell sorting with maximum flexibility

The BD FACSVantage™ SE flow cytometry system is the most flexible high-speed cell sorting instrument for the research laboratory. This instrument platform can accommodate any laser wavelength and a variety of nozzle sizes for maximum performance for any research application.

To meet the demands of today's research, we've improved the BD FACSVantage SE and its many options to allow for customization to meet any possible cell sorting need. This technology is accessible through a familiar user interface and efficient ergonomic design. The BD FACSVantage SE flow cytometry system is fully modular.

BD FACSVantage SE System

- Multicolor analysis and sorting
- BD FACStation™ data management system
- Flexible and modular for future upgrades
- Non-rectangular sort windows



Options

- Digital option
- BD TurboSort™ option
- BD MacroSort™ option
- BD CloneCyt™ Plus
- BD FACS™ Accudrop
- Pulse Processor

*Class I (1) Laser Product

BD Biosciences Fluorochrome Reference Chart

Typical Instrument Configuration

Instrument	Laser	Excitation Laser Line (nm)	Fluorescence Channel	Fluorochromes Provided by BD Biosciences
BD FACSCalibur™	Argon (L1)	488	FL1 Green	FITC Alexa Fluor® 488
			FL2 Yellow	PE
			FL3 Red	PE-Cy5* PerCP PerCP-Cy5.5 PE-Cy7
	Red Diode (L2)	635	FL4 Red	APC* Alexa Fluor® 647
BD FACSCanto™	Argon (L1)	488	Green	FITC
			Yellow	PE
			Red	PerCP PerCP-Cy5.5
	HeNe (L2)	633	Infra Red	PE-Cy7
			Red	APC*
		Infra Red	APC-Cy7	
BD FACSVantage™ SE (typical setup)	Argon (L1)	488	FL1 Green	FITC Alexa Fluor® 488
			FL2 Yellow	PE
			FL6 Red	PE-Texas Red® PE-Cy5* PerCP-Cy5.5
			FL3 Infra Red	PE-Cy7
	Krypton (L2)	407	FL4 (1) Blue	Alexa Fluor® 405 Pacific Blue®
	HeNe (L2 or L3)	633	FL4 (2) Red	APC* Alexa Fluor® 647
		FL5 Infra Red	APC-Cy7	
BD™ LSR II (typical setup)	Argon (L1)	488	Green	FITC Alexa Fluor® 488
			Yellow	PE
			Red	PE-Texas Red® PE-Cy5* PerCP PerCP-Cy5.5
	HeNe (L2)	633	Red	APC* Alexa Fluor® 647
			Infra Red	APC-Cy7
	UV (L3)	355	Violet	
	Violet (L4)	405	Blue	Alexa Fluor® 405 Pacific Blue®
			Green	AmCyan
		Blue	Alexa Fluor® 405 Pacific Blue®	
BD FACSAria™ (typical setup)	Argon (L1)	488	Green	FITC Alexa Fluor® 488
			Yellow	PE
			Red	PE-Texas Red® PE-Cy5* PerCP
			Far Red	PerCP-Cy5.5
	HeNe (L2)	633	Infra Red	PE-Cy7
			Red	APC* Alexa Fluor® 647
Violet (L3)	407	Infra Red	APC-Cy7	
		Green	AmCyan	
		Blue	Alexa Fluor® 405 Pacific Blue®	
BD FACSAria™ (typical setup)	Green Diode (L1)	532	Yellow	PE
			Far Red	PerCP-Cy5.5 PE-Cy7
	Red Diode (L2)	635	Red	APC* Alexa Fluor® 647
			Infra Red	APC-Cy7

* APC and PE-Cy5 may be used together on instruments with cross-beam compensation.

Note: Custom BD FACSAria™ instruments equipped with either 355nm or 375nm lasers can be ordered to enable stem cell applications using Hoescht or Dapi dyes.

Identification of a Primitive Progenitor Cell Subset Using Hoechst 33342

Characterization of hematopoietic stem cells and their subsets has traditionally been performed with cell surface staining using fluorescently labeled monoclonal antibodies. A cell marker of choice is CD34,⁴⁷ a transmembrane phosphoglycoprotein that is selectively expressed on vascular endothelial cells, tissue fibroblasts or stromal cells, and hematopoietic cells. CD34 has been used in different combinations with other cell surface markers and is the most prevalent marker currently used for identification of mouse and human progenitor cells. The cells identified by these methods have biological properties that include self renewal and differentiation into mature cell lineages.^{48,49}

Fluorescent probes such as rhodamine 123^{50,51} and pyronin Y⁵² have been shown to fractionate subpopulations of hematopoietic progenitor cells when used with other cell surface markers. Goodell^{53,54,55} and co-workers have described a method using dual wavelength flow cytometry to define a small (0.01 – 0.10%) population of primitive progenitor cells that selectively efflux the dye Hoechst 33342 (bis-benzimide) in mouse, primate, and human bone marrow, as well as in a variety of tissues. This population of cells is termed “side population” (SP) cells. These cells are enumerated by exciting the Hoechst-labeled cells with the 350 nm line of a UV laser and then simultaneously measuring the fluorescence using blue and far red emission filters. Of particular interest is that these cells express CD34-low to negative and are thought to be the most primitive of progenitor cells identified to date.




Figure 1. Optical diagram showing the position of the filter set that includes the following filters: DM610 SP, BP450/20, and EF675 LP.

Method



The SP hematopoietic stem cells are identified by the cells' ability to differentially efflux the Hoechst dye via a multi-drug-like transport mechanism. This is a dynamic process; therefore, care must be taken during the labeling protocol by judiciously observing the time and temperature requirements.⁵⁶

The BD Biosciences Stem Cell Side-Population Filter Set (Cat. No. 341063) provides operators of the BD FACSVantage™ SE (with or without the BD FACSDiva digital electronics option) with a convenient way to set up flow cytometers equipped with an ultraviolet laser in the second or third positions to perform these experiments. The DM610 SP used to distinguish the red from the blue fluorescence signals comes premounted on a mirror holder that is attached to the instrument's dichroic mirror holder. The BP450/20nm and EF675-LP filters are then installed in front of the PMT detectors as shown in **Figure 1**.

The BP450/20nm is used to discriminate the fluorescent blue signal and the EF675-LP red signal. The events are collected using linear amplification (**Figure 2**).

For best results:

1. Tune UV laser to 350 nm and 100 mW for defining this rare sub-population of cells.
2. Display data as correlated Hoechst blue and Hoechst red in linear mode.
3. Align cytometer to minimize the CVs in all fluorescence channels for this application.
4. Keep the sample at 4°C throughout the analysis or sort to minimize dye efflux.
5. Collect sorted S-P cells at 4°C in a small volume of DMEM containing 10% FBS and 10 mM HEPES (use a 1.5 ml Eppendorf centrifuge tube).
6. Maintain a low sample differential for maximum resolution.
7. Collect a large list-mode data file of 100,000 events for the analysis.
8. Gate out the dead cells using propidium iodide (PI).

The phenomenon of “red fluorescence” from Hoechst 33342 is not fully understood. It is thought that the way in which the dye binds to the chromatin conformation has a direct effect on the emission spectra.⁵⁷




Figure 2. Murine C57BL/6 bone marrow was stained using Hoechst 33342. A data set of 100,000 events was collected in list mode. The contour plot at 5% probability of ungated events shows the S-P population (R1) with frequencies that range from 0.02 to 0.08%.



BD FACSVantage SE

BD IMag™ Magnetic Cell Separation System


Mouse Hematopoietic Progenitor (Stem) Cell Enrichment Set – DM

Cat. No. 558451

The BD IMag™ Mouse Hematopoietic Progenitor Cell Enrichment Set – DM reacts with cells from the major hematopoietic cell lineages, such as T lymphocytes, B lymphocytes, monocytes/macrophages, granulocytes, and erythrocytes. The Biotinylated Mouse Lineage Depletion Cocktail contains biotinylated monoclonal antibodies to mouse CD3e (CD3ε chain), CD11b (Integrin α_M chain), CD45R/B220, Ly-6G and Ly-6C (Gr-1), and


TER-119/ Erythroid Cells (Ly-76). The BD IMag™ Streptavidin Particles Plus – DM are magnetic nanoparticles that have streptavidin covalently conjugated to their surfaces. This Set is designed for the immunomagnetic enrichment of hematopoietic progenitors from mouse bone marrow by depletion of cells committed to the T- and B-lymphocytic, myeloid (monocytic and granulocytic), and erythroid lineages.

The set contains sufficient reagents to label 1x10⁹ bone marrow cells.



Depletion of lineage-committed cells from mouse bone marrow. BALB/c bone-marrow cells were labeled with the BD IMag™ Mouse Hematopoietic Progenitor Enrichment Set – DM and separated on the BD IMagnet™ (Cat. no. 552311) according to the accompanying protocol. To demonstrate the efficiency of the depletion, unmanipulated bone marrow cells and the final depleted fraction were stained with APC-conjugated anti-mouse CD117 mAb 2B8 (Cat. No. 553356) to detect hematopoietic progenitors, and with PE-conjugated mAb TER-119 (Cat. No. 553673), PE-conjugated mAb RA3-6B2 (Cat. No. 553089/553090), and PE-conjugated mAb M1/70 (Cat. No. 557397/553311) to detect lineage-committed cells. The percentage of positive cells is indicated in each panel; placement of each marker is based on staining with the appropriate isotype control (data not shown). The final depleted fraction contains a greatly increased proportion of CD117+ cells and less than 5% of lineage-positive contaminants.

Depletion Flow Chart for Mouse Hematopoietic Progenitor (Stem) Cells



Imaging – Seeing is Believing

BD CARV II™ Confocal Imager

The BD CARV II™ confocal imager delivers high-resolution CCD confocal imaging in an easy-to-use and cost-effective optical package that fits on your existing microscope. High-speed multi-point confocal scanning, combined with high quantum efficiency CCD cameras, minimizes photobleaching and allows real-time imaging and recording at up to 100 fps. A long life arc source coupled to the instrument via an alignment-free light guide allows for full spectrum (360 nm - 700 nm) confocal imaging of virtually any fluorescent probe. Automation of internal multi-position excitation, dichroic, and emission filter wheels permits fast multi-dimensional imaging of up to five or more fluorescent probes in the same sample.




The BD CARV II™ system provides an easy-to-use confocal upgrade to your existing microscope without sacrificing image quality or three dimensional imaging capabilities.

BD Pathway™ Bioimager


The BD Pathway™ bioimager is an automated, confocal, real-time, single-cell kinetic and endpoint imaging system that has been integrated into a single, compact unit. The BD Pathway was designed to provide high-resolution, automated confocal imaging with sophisticated imaging software to assist in developing better cell-based assays. The system allows the user to explore biological events without many of the restrictions of conventional microscope-based high-content imaging systems.




The BD Pathway™ bioimager provides a wide range of imaging capabilities including kinetic, endpoint, and three dimensional confocal imaging. The environmentally controlled imaging chamber with liquid handling provides unsurpassed live cell kinetic imaging possibilities.




Merged montaged* image of hES cells labeled with SSEA-4 and Hoechst 33342 on BD Matrigel™ Matrix, 20X objective.



Merged montaged* transmitted light image of hES cells labeled with Oct-3, 20X objective.



Merged montaged* image of hES cells labeled with SSEA-1 and Hoechst 33342 on MEF feeder cells.



Merged montaged* image of live (blue - stained with Hoescht 33342) and dead (pink - stained with Propidium Iodide) hES cells.

* These are all montaged images: Stem cell colonies are large and at 20X objective only on a small fraction of the colony would be visible in the image. The precision of the BD Pathway™ bioimager stage allows users to take multiple images right next to each other (a montage) to piece together a large subject matter as if it were one image.

Product Listing

Reagents for Human Stem Cell Research

Antibodies to Stem Cell Growth Factors

DESCRIPTION	REACT	CLONE	ISOTYPE	APPS	REG	FORMAT	SIZE	CAT. NO.
CD271 (NGF Receptor)	Hu	C40-1457	Mouse IgG ₁ , κ	FCM, IF	RUO	Purified	0.1 mg	557194
				FCM	RUO	Biotin	100 tests	557195
				FCM	RUO	PE	100 tests	557196
basic FGF	Chick, Dog, Hu, Ms, Rat	3	Mouse IgG ₁	IF, WB	RUO	Purified	50 µg	610870
				IF, WB	RUO	Purified	150 µg	610871
	6	Mouse IgG ₁	WB	RUO	HRP	150 µg	610873	
			IF, IHC, IP, WB	RUO	Purified	50 µg	610072	
			IF, IHC, IP, WB	RUO	Purified	150 µg	610073	
EGF Receptor	Hu, Ms	13	Mouse IgG ₁	IF, IP, WB	RUO	Purified	50 µg	610016
				IF, IP, WB	RUO	Purified	150 µg	610017
				IF	RUO	FITC	50 µg	612554
	Hu	EGFR1	Mouse IgG _{2b} , κ	IF	RUO	FITC	150 µg	612555
				FCM, IF	RUO	Purified	0.1 mg	555996
				FCM	RUO	PE	100 Tests	555997
EGF Receptor, Activated Form	Hu	74	Mouse IgG ₁	IF, IP, WB	RUO	Purified	50 µg	610025
				IF, IP, WB	RUO	Purified	150 µg	610026
G-CSF	Hu	BVD11-37G10	Rat IgG _{2a}	ELISA Det.	RUO	Biotin	0.5 mg	554670
		BVD13-3A5	Rat IgG ₁	ELISA Cap., WB	RUO	Purified	1.0 mg	551342
GM-CSF	Hu	BVD2-21C11	Rat IgG _{2a}	IC/FCM Block, IP/WB	RUO	Purified	0.1 mg	554503
				ELISA Det.	RUO	Biotin	0.5 mg	554505
				IC/FCM	RUO	PE	0.1 mg	554507
		BVD2-23B6	Rat IgG _{2a}	Neu	RUO	NA/LE	0.5 mg	554501
				ELISA Cap., IP/WB	RUO	Purified	0.5 mg	554502
IL-3	Hu	BVD3-1F9	Rat IgG ₁	ELISA Det.	RUO	Biotin	0.5 mg	554674
				IC/FCM	RUO	PE	0.1 mg	554676
		BVD8-3G11	Rat IgG ₁	Neu	RUO	NA/LE	0.5 mg	554671
				ELISA Cap., WB	RUO	Purified	0.5 mg	554672
IL-6	Hu	AS12	Mouse IgG ₁ , κ	FCM, IC	RUO (GMP)	FITC	50 tests	340526
				FCM, IC	RUO (GMP)	PE	50 tests	340527
		MQ2-13A5	Rat IgG ₁	Neu	RUO	NA/LE	0.5 mg	554541
				ELISA Cap., IC/FCM Block, WB	RUO	Purified	0.5 mg	554543
				IC/FCM	RUO	FITC	0.1 mg	554544
		MQ2-39C3	Rat IgG _{2a}	IC/FCM	RUO	PE	0.1 mg	554545
				ELISA Det., WB	RUO	Biotin	0.5 mg	554546
		MQ2-6A3	Rat IgG _{2a} , κ	IC/FCM Block, ICC	RUO	Purified	0.25 mg	559068
				IC/FCM	RUO	FITC	0.1 mg	554696
				IC/FCM	RUO	PE	0.1 mg	554697
IC/FCM	RUO			PE	100 tests	559331		
Leukemia Inhibitory Factor (LIF, HILDA)	Hu	1F10	Mouse IgG ₁ , κ	FCM	RUO	APC	0.1 mg	558098
				FCM	RUO	PE	100 tests	558571

Reagents for Human Stem Cell Research

Antibodies to Stem Cell Markers

DESCRIPTION	REACT	CLONE	ISOTYPE	APPS	REG	FORMAT	SIZE	CAT. NO.
CD4	Hu	RPA-T4	Mouse IgG ₁ , κ	FA, FCM	RUO	NA/LE	0.5 mg	555343
				FCM, IHC(Fr)	RUO	Purified	0.1 mg	555344
				IHC(Fr), IHC(Zn)	RUO	Purified	1 ml	550369
				FCM	RUO	Biotin	100 tests	555345
				FCM	RUO	FITC	100 tests	555346
				FCM	RUO	PE	100 tests	555347
				FCM	RUO	APC	100 tests	555349
				FCM	RUO	Alexa Fluor® 488	100 tests	557695
				FCM	RUO	Alexa Fluor® 647	100 tests	557707
				FCM	RUO	Alexa Fluor® 700	0.1 mg	557922
				FCM	RUO	Pacific Blue®	0.1 mg	558116
				FCM	RUO	APC-Cy7	100 tests	557871
				FCM	RUO	PE-Cy5	100 tests	555348
		SK3	Mouse IgG ₁ , κ	FCM	RUO (GMP)	Purified	200 tests	346320
				FCM	RUO (GMP)	Biotin	100 tests	347321
				FCM	RUO (GMP)	PE	100 tests	347327
				FCM	RUO (GMP)	PerCP	100 tests	347324
				FCM	RUO (GMP)	PerCP-Cy5.5	50 tests	341654
				FCM	RUO (GMP)	APC	100 tests	340443
				FCM	RUO (GMP)	APC-Cy7	100 tests	341095
				FCM	RUO (GMP)	PE-Cy7	100 tests	348789
		SK3, SK4	Mouse IgG ₁ , κ	FCM	RUO (GMP)	FITC, FITC	100 tests	347413
CD4 v4	Hu	L120	Mouse IgG ₁ , κ	FCM	RUO (GMP)	FITC	50 tests	340422
				FCM	RUO (GMP)	PE	50 tests	340419
CD8	Hu	G42-8	Mouse IgG _{2b} , κ	FCM	RUO	Purified	0.1 mg	551346
				FCM	RUO	FITC	0.1 mg	551347
		HIT8a	Mouse IgG ₁ , κ	FCM	RUO	NA/LE	0.5 mg	555630
				FCM, IHC(Fr)	RUO	Purified	0.1 mg	555631
				IHC(Fr), IHC(Zn)	RUO	Purified	1 ml	550372
				FCM	RUO	FITC	100 tests	555634
				FCM	RUO	PE	100 tests	555635
				FCM	RUO	PE-Cy5	100 tests	555636
		RPA-T8	Mouse IgG ₁ , κ	FCM	RUO	NA/LE	0.5 mg	555363
				FCM, IHC(Fr)	RUO	Purified	0.1 mg	555364
				FCM	RUO	Biotin	100 tests	555365
				FCM	RUO	FITC	100 tests	555366
				FCM	RUO	PE	100 tests	555367
				FCM	RUO	APC	100 tests	555369
				FCM	RUO	Alexa Fluor® 488	100 tests	557696
				FCM	RUO	Alexa Fluor® 647	100 tests	557708
				FCM	RUO	Alexa Fluor® 700	0.1 mg	557945
				FCM	RUO	Pacific Blue®	0.1 mg	558207
				FCM	RUO	PE-Cy5	100 tests	555368
				FCM	RUO	PE-Cy7	100 tests	557746
				SK1	Mouse (BALB/c) IgG ₁ , κ	FCM	RUO (GMP)	Purified
		FCM	RUO (GMP)			Biotin	100 tests	347311
		FCM	RUO (GMP)			FITC	100 tests	347313
		FCM	RUO (GMP)			PerCP	100 tests	347314
		FCM	RUO (GMP)			PerCP-Cy5.5	50 tests	341051
		FCM	RUO (GMP)			APC	100 tests	340584
		FCM	RUO (GMP)			APC-Cy7	100 tests	348793
		FCM	RUO (GMP)			PE-Cy7	100 tests	335787
		FCM	RUO			APC-Cy7	100 tests	557834
		FCM	RUO (GMP)			Purified	100 tests	347350
		SK2	Mouse IgG _{2b} , κ	FCM	RUO (GMP)	Purified	100 tests	347350

Product Listing *(continued)*

Reagents for Human Stem Cell Research

Antibodies to Stem Cell Markers (continued)

DESCRIPTION	REACT	CLONE	ISOTYPE	APPS	REG	FORMAT	SIZE	CAT. NO.		
CD15 and SSEA-1	Hu ,Ms,Rat	HI98	Mouse IgM, κ	FCM, IHC(F), IHC(Fr)	RUO	Purified	0.1 mg	555400		
				IHC(F), IHC(Fr)	RUO	Purified	1 ml	550382		
				FCM	RUO	FITC	100 tests	555401		
				FCM	RUO	PE	100 tests	555402		
		MMA	Mouse IgM,κ	FCM	RUO	APC	100 tests	551376		
				FCM	RUO (GMP)	Purified	100 tests	347420		
		W6D3	Mouse IgG ₁ , κ	FCM	RUO (GMP)	FITC	100 tests	347423		
				FA, IHC(F)	RUO	Purified	0.2 mg	559045		
						FCM	RUO	Purified	0.1 mg	557895
						FCM	RUO	PE	100 tests	558084
CD24	Hu	ML5	Mouse IgG _{2b} , κ	FCM, IHC(Fr)	RUO	Purified	0.1 mg	555426		
				FCM	RUO	FITC	100 tests	555427		
				FCM	RUO	PE	100 tests	555428		
CD30 (Ki-1)	Hu	Ber-H83	Mouse IgG ₁ , κ	FCM	RUO (GMP)	FITC	50 tests	341644		
				FCM	RUO (GMP)	PE	50 tests	341645		
		Ber-H8	Mouse IgG ₁ , κ	FCM, IHC(F), IHC(Fr)	RUO	Purified	0.5 mg	555827		
				FCM	RUO	Biotin	0.5 mg	555828		
				FCM	RUO	FITC	100 tests	555829		
FCM	RUO	PE	100 tests	550041						
CD33	Hu	HIM3-4	Mouse IgG ₁ , κ	FCM	RUO	Purified	0.1 mg	555625		
				FCM	RUO	FITC	100 tests	555626		
		P67.6	Mouse IgG ₁ , κ	FCM	RUO (GMP)	Purified	100 tests	347780		
				FCM	RUO (GMP)	FITC	50 tests	340533		
				FCM	RUO (GMP)	PE	100 tests	347787		
				FCM	RUO (GMP)	PerCP-Cy5.5	50 tests	341650		
		WM53	Mouse IgG ₁ , κ	FCM	RUO (GMP)	APC	100 tests	340474		
				FCM	RUO (GMP)	PE-Cy7	100 tests	333946		
				FCM	RUO	Purified	0.1 mg	555449		
				FCM	RUO	PE	100 tests	555450		
				FCM	RUO	APC	100 tests	551378		
				FCM	RUO	PE-Cy5	100 tests	551377		
CD34 (gp 105-120)	Hu	563	Mouse IgG ₁ , κ	FCM	RUO	Purified	0.1 mg	550760		
				FCM	RUO	PE	100 tests	550761		
		581	Mouse IgG ₁ , κ	FCM, IHC(F), IHC(Fr)	RUO	Purified	0.1 mg	555820		
				IHC(F)*, IHC(Fr)	RUO	Purified	1 ml	550390		
				FCM	RUO	FITC	100 tests	555821		
				FCM	RUO	PE	100 tests	555822		
						FCM	RUO	APC	100 tests	555824
						FCM	RUO	PE-Cy5	100 tests	555823
		8G12	Mouse IgG ₁ , κ	FCM	RUO (GMP)	Purified	100 tests	348050		
				FCM	RUO (GMP)	FITC	100 tests	348053		
				FCM	RUO (GMP)	PE	100 tests	348057		
				FCM	RUO (GMP)	PerCP	50 tests	340430		
				FCM	RUO (GMP)	PerCP-Cy5.5	50 tests	347203		
				FCM	RUO (GMP)	APC	100 tests	340441		
				FCM	RUO (GMP)	PE-Cy7	100 tests	348791		
My10	Mouse IgG ₁ , κ	FCM	RUO (GMP)	Purified	100 tests	347660				

Reagents for Human Stem Cell Research

Antibodies to Stem Cell Markers (continued)

DESCRIPTION	REACT	CLONE	ISOTYPE	APPS	REG	FORMAT	SIZE	CAT. NO.
CD34	Cyno, Rhe, Hu	563	Ms IgG ₁ , κ	FCM	RUO	PE	50 tests	550619
CD34 (gp 105-120)	Ms	RAM34	Rat IgG _{2b} , κ	FCM	RUO	Purified	0.5 mg	553731
				IHC(Fr), IHC(Zn)	RUO	Purified	1 ml	550537
				FCM	RUO	Biotin	0.5 mg	553732
				FCM	RUO	FITC	0.5 mg	553733
				FCM	RUO	PE	0.2 mg	551387
CD34 (gp 105-120)	Dog	2E9	Ms IgG ₁ , κ	FCM	RUO	Biotin	0.1 mg	550427
CD34 (gp 105-120)	Dog	1H6	Ms IgG ₁ , κ	FCM	RUO	PE	0.1 mg	559369

CD34 Antibodies – Multicolor

DESCRIPTION	REACT	CLONE	ISOTYPE	APPS	REG	FORMAT	SIZE	CAT. NO.
CD15/CD34	Hu	MMA, 8G12	Ms IgM, κ, Ms IgG ₁ , κ	FCM	RUO	FITC, PE	50 tests	341128
CD45/CD34 combination	Hu	2D1, 8G12	Ms IgG ₁ , κ Ms IgG ₁ , κ	FCM	RUO	FITC, PE	50 tests	341071
HLA-DR/CD34	Hu	L243, 8G12	Ms IgG _{2b} , κ, Ms IgG ₁ , κ	FCM	RUO	FITC, PE	50 tests	341146

Product Listing *(continued)*

Reagents for Human Stem Cell Research

Antibodies to Stem Cell Markers (continued)

DESCRIPTION	REACT	CLONE	ISOTYPE	APPS	REG	FORMAT	SIZE	CAT. NO.
CD38	Hu	22	Mouse IgG ₁ , κ	IF, WB	RUO	Purified	50 µg	611114
				IF, WB	RUO	Purified	150 µg	611115
		HB7	Mouse IgG ₁ , κ	FCM	RUO (GMP)	Purified	100 tests	347680
				FCM	RUO (GMP)	FITC	50 tests	340927
				FCM	RUO (GMP)	PE	50 tests	342371
				FCM	RUO (GMP)	PE	100 tests	347687
	FCM			RUO (GMP)	APC	100 tests	340439	
	HIT2	Mouse IgG ₁ , κ	FCM, IHC(Fr)	RUO	Purified	0.1 mg	555458	
			FCM	RUO	FITC	100 tests	555459	
			FCM	RUO	PE	100 tests	555460	
			FCM	RUO	PerCP-Cy5.5	0.1 mg	551400	
			FCM	RUO	APC	100 tests	555462	
			FCM	RUO	PE-Cy5	100 tests	555461	
CD45 (Leukocyte Common Antigen, Ly-5)			Hu	2D1	Mouse IgG ₁ , κ	FCM	RUO (GMP)	Purified
	FCM	RUO (GMP)				FITC	100 tests	347463
	FCM	RUO (GMP)				PerCP	100 tests	347464
	FCM	RUO (GMP)				PerCP-Cy5.5	50 tests	340953
	FCM	RUO (GMP)				APC	100 tests	340943
	Hu, Ms, Rat	69	Mouse IgG ₁ , κ	FCM	RUO (GMP)	APC-Cy7	100 tests	348795
				FCM	RUO	APC-Cy7	100 tests	557833
				IF, IP, WB	RUO	Purified	50 µg	610265
	Hu	HI30	Mouse IgG ₁ , κ	IF, IP, WB	RUO	Purified	150 µg	610266
				FCM, IHC(F), IHC(Fr)	RUO	Purified	0.1 mg	555480
				FCM	RUO	Biotin	100 tests	555481
				FCM	RUO	FITC	100 tests	555482
				FCM	RUO	PE	100 tests	555483
Hu	HI30	Mouse IgG ₁ , κ	FCM	RUO	APC	100 tests	555485	
			FCM	RUO	PE-Cy5	100 tests	555484	
			FCM	RUO	PE-Cy7	100 tests	557748	
			FCM	RUO	PE-Cy7	100 tests	557748	
			FCM	RUO	PE-Cy7	100 tests	557748	
CD45/CD34 combination	Hu	2D1, 8G12		FCM	RUO (GMP)	FITC, PE	50 tests	341071
CD56 (N-CAM)	Hu	B159	Mouse IgG ₁ , κ	FCM	RUO	NA/LE	0.5 mg	555513
				FCM	RUO	Purified	0.1 mg	555514
				FCM	RUO	Biotin	100 tests	555515
				FCM	RUO	PE	100 tests	555516
				FCM	RUO	APC	100 tests	555518
				FCM	RUO	Alexa Fluor® 488	100 tests	557699
				FCM	RUO	Alexa Fluor® 647	100 tests	557711
				FCM	RUO	Alexa Fluor® 700	0.1 mg	557919
				FCM	RUO	PE-Cy7	100 tests	557747
				FCM	RUO	PE-Cy5	100 tests	555517
	Bab, Cyno, Hu, Rhe	MY31	Mouse IgG ₁ , κ	FCM	RUO (GMP)	Purified	100 tests	347740
				FCM	RUO (GMP)	PE	100 tests	347747
				FCM	RUO	Purified	0.2 mg	559049
	Hu	NCAM16.2	Mouse (BALB/c) IgG _{2b} , κ	FCM	RUO (GMP)	FITC	50 tests	340410
				FCM	RUO (GMP)	PE	50 tests	340363
FCM				RUO (GMP)	APC	100 tests	341025	
FCM				RUO (GMP)	PE-Cy7	100 tests	335791	
FCM				RUO	Purified	0.1 mg	558171	
Hu	p282 (H19)	Mouse IgG _{2b} , κ	ELISA, FA, FCM, IHC(Fr), IP	RUO	Purified	0.2 mg	559043	
			FCM, IF	RUO	Purified	0.1 mg	555761	
CD59	Hu	p282 (H19)	Mouse IgG _{2b} , κ	FCM	RUO	Biotin	100 tests	555762
				FCM	RUO	FITC	100 tests	555763
				FCM	RUO	PE	100 tests	555764
				FCM	RUO	PE	100 tests	555764

Reagents for Human Stem Cell Research

Antibodies to Stem Cell Markers (continued)

DESCRIPTION	REACT	CLONE	ISOTYPE	APPS	REG	FORMAT	SIZE	CAT. NO.
CD71 (Transferrin Receptor)	Hu	L01.1	Mouse IgG _{2b} , κ	FCM	RUO (GMP)	Purified	100 tests	347510
				FCM	RUO (GMP)	FITC	100 tests	347513
				FCM	RUO (GMP)	APC	100 tests	341028
		M-A712	Mouse IgG _{2b} , κ	FCM, IHC(Fr)	RUO	Purified	0.1 mg	555534
				FCM	RUO	Biotin	100 tests	555535
				FCM	RUO	FITC	100 tests	555536
				FCM	RUO	PE	100 tests	555537
				FCM	RUO	APC	100 tests	551374
FCM	RUO	PE-Cy5	100 tests	551143				
CD90 (Thy-1)	Hu, Pig, Rab	5.00E+10	Mouse IgG ₁ , κ	FCM	RUO	NA/LE	0.5 mg	555592
				FCM, IHC(Fr)	RUO	Purified	0.1 mg	555593
				IHC(Fr), IHC(Zn)	RUO	Purified	1 ml	550402
				FCM	RUO	Biotin	0.1 mg	555594
				FCM	RUO	FITC	0.1 mg	555595
				FCM	RUO	PE	0.1 mg	555596
				FCM	RUO	APC	0.1 mg	559869
				FCM	RUO	PE-Cy5	0.1 mg	555597
CDw93 (C1qRp)	Hu	R139	Mouse IgG _{2b} , κ	FA, FCM	RUO	NA/LE	0.5 mg	552954
				FCM, IF, WB	RUO	Purified	0.1 mg	551087
				FCM	RUO	FITC	0.1 mg	551531
		R3	Mouse IgM, κ	FCM	RUO	PE	0.1 mg	551509
				FCM, IP, WB	RUO	Purified	0.1 mg	551454
				FCM	RUO	Biotin	100 tests	552117
CD105 (Endoglin)	Hu	266	Mouse IgG ₁ , κ	IHC(Fr), IHC(Zn), IP	RUO	Purified	0.1 mg	555690
		35	Mouse IgG ₁ , κ	IF, WB	RUO	Purified	50 µg	611314
				IF, WB	RUO	Purified	150 µg	611315
CD109	Hu	TEA 2/16	Mouse IgG ₁ , κ	FCM, IHC(Fr)	RUO	Purified	0.1 mg	556039
				FCM	RUO	PE	100 tests	556040
CD117 (c-kit)	Hu	104D2	Mouse IgG ₁ , κ	FCM	RUO (GMP)	PE	50 tests	340529
				FCM	RUO (GMP)	PerCP-Cy5.5	50 tests	333944
				FCM	RUO (GMP)	APC	100 tests	341096
		28	Mouse IgG _{2b} , κ	WB	RUO	Purified	50 µg	612318
				WB	RUO	Purified	150 µg	612319
		Dog, Hu, Rab	YB5.B8	Mouse IgG ₁ , κ	FCM, IF	RUO	Purified	0.1 mg
	FCM				RUO	PE	0.1 mg	555714
	FCM				RUO	APC	0.1 mg	550412
FCM	RUO	PE-Cy5	0.1 mg	559879				
CD144 (VE-cadherin)	Hu	55-7H1	Mouse IgG ₁ , κ	IHC	RUO	Purified	0.1 mg	555661
CD184 (CXCR4, Fusin)	Hu	12G5	Mouse IgG _{2b} , κ	FA, FCM	RUO	NA/LE	0.5 mg	555971
				FCM	RUO	Purified	0.1 mg	555972
				FCM	RUO	Biotin	100 tests	555973
				FCM	RUO	PE	100 tests	555974
				FCM	RUO	APC	100 tests	555976
				FCM	RUO	PE-Cy5	100 tests	555975
		1D9	Rat IgG _{2b} , κ	FCM	RUO	Purified	0.1 mg	551413
				FCM	RUO	Biotin	100 tests	551970
				FCM	RUO	PE	0.2 mg	551510
CD318 (CDCP1)	Hu	CUB1	Mouse IgG _{2b} , κ	FCM	RUO	PE	0.1 mg	558209
				FCM	RUO	Alexa Fluor® 647	0.1 mg	558212
CDw338 (ABCG2)	Hu	5D3	Mouse IgG _{2b} , κ	FCM	RUO	Purified	0.1 mg	552823
γ Enolase (NSE)	Hu, Ms	NSE-P2	Mouse IgG ₁ , κ	IHC, WB	RUO	Purified	100 µg	557820
Disialoganglioside GD2	Hu	14.G2a	Mouse IgG _{2b} , κ	FCM, IF, IHC(Fr)	RUO	Purified	0.1 mg	554272
Disialoganglioside GD3	Hu	MB3.6	Mouse IgG ₃ , κ	FA, FCM, IF, IHC(Fr)	RUO	Purified	0.1 mg	554274

Product Listing *(continued)*

Reagents for Human Stem Cell Research

Antibodies to Stem Cell Markers (continued)

DESCRIPTION	REACT	CLONE	ISOTYPE	APPS	REG	FORMAT	SIZE	CAT. NO.
EGF Receptor	Hu, Ms	13	Mouse IgG ₁	IF, IP, WB	RUO	Purified	50 µg	610016
				IF, IP, WB	RUO	Purified	150 µg	610017
				IF	RUO	FITC	50 µg	612554
	Hu	EGFR1	Mouse IgG _{2b} , κ	FCM, IF	RUO	Purified	0.1 mg	555996
				FCM	RUO	PE	100 Tests	555997
				FCM	RUO	PE-Cy7	0.1 mg	558018
GAD65	Hu, Ms, Pig, Rat	GAD-6	Mouse IgG _{2b} , κ	IHC(Fr), IP, WB	RUO	Purified	0.1 mg	559931
HLA-DR	Dog, Hu	G46-6 (L243)	Mouse IgG _{2b} , κ	FCM	RUO	NA/LE	0.5 mg	555809
				FCM	RUO	Purified	0.1 mg	555810
				FCM	RUO	FITC	100 tests	555811
				FCM	RUO	PE	100 tests	555812
				FCM	RUO	APC	100 tests	559866
	Hu	L243	Mouse IgG _{2b} , κ	FCM	RUO (GMP)	Purified	100 tests	555813
				FCM	RUO (GMP)	Biotin	100 tests	347360
				FCM	RUO (GMP)	FITC	100 tests	347363
				FCM	RUO (GMP)	PE	100 tests	347367
				FCM	RUO (GMP)	PerCP	100 tests	347364
	Dog, Hu, Rab	TÜ36	Mouse IgG _{2b} , κ	FCM	RUO (GMP)	APC	100 tests	340549
				FCM	RUO	Purified	0.1 mg	555559
				FCM	RUO	PE	100 tests	555561
				FCM	RUO	APC	100 tests	559868
HLA-DR/CD34	Hu	L243, 8G12	FCM	FCM	RUO (GMP)	FITC, PE	50 tests	341146
				FCM	RUO (GMP)	PE	100 tests	551375
Leukemia Inhibitory Factor (LIF, HILDA)	Hu	1F10	Mouse IgG ₁ , κ	FCM	RUO	PE	100 tests	558571
Nestin	Hu, Ms, Rat	25	Mouse IgG ₁ , κ	WB	RUO	Purified	50 µg	611658
				WB	RUO	Purified	150 µg	611659
p63	Hu, Ms	4A4	Mouse IgG _{2b} , κ	IHC(F), IHC(Fr), WB	RUO	Purified	0.1 mg	559951
		5G8	Mouse IgG ₁ , κ	WB	RUO	Purified	0.1 mg	559952
		Y4A3	Mouse IgG _{2b} , κ	IHC(F), IHC(Fr), WB	RUO	Purified	0.1 mg	550025
RIP	Hu	G322-2	Mouse IgG ₁ , κ	IP, WB	RUO	Purified	50 µg	551041
				IP, WB	RUO	Purified	150 µg	551042
		Polyclonal	Rabbit	WB	RUO	Serum	0.1 ml	559689

Reagents for Human Stem Cell Research

Magnetic Cell Separation Reagents

DESCRIPTION	REACT	CLONE	ISOTYPE	APPS	REG	FORMAT	SIZE	CAT. NO.
CD4 Magnetic Particles - DM	Hu, Rhe	L200	Mouse IgG ₁ , κ	Sep	RUO	BD IMag-DM	1 x 10 ⁹ cells	557767
CD4 Magnetic Particles - MSC	Hu, Rhe	L200	Mouse IgG ₁ , κ	Sep	RUO	BD IMag-MS	1 x 10 ⁹ cells	557768
CD4 T Lymphocyte Enrichment Set - DM	Hu			Sep	RUO	BD IMag-DM, Biotin	1 x 10 ⁹ cells	557939
CD8 Magnetic Particles - DM	Hu, Rhe	SK1	Mouse (BALB/c) IgG ₁ , κ	Sep	RUO	BD IMag-DM	1 x 10 ⁹ cells	557766
CD8 Magnetic Particles - MSC	Hu, Rhe	SK1	Mouse (BALB/c) IgG ₁ , κ	Sep	RUO	BD IMag-MS	1 x 10 ⁹ cells	557765
CD8 T Lymphocyte Enrichment Set - DM	Hu			Sep	RUO	BD IMag-DM, Biotin	1 x 10 ⁹ cells	557941
CD56 Magnetic Particles - DM	Hu	NCAM16.2	Mouse (BALB/c) IgG _{2b} , κ	Sep	RUO	BD IMag-DM	1 x 10 ⁹ cells	557775

Reagents for Mouse Stem Cell Research

Antibodies to Stem Cell Growth Factors

DESCRIPTION	REACT	CLONE	ISOTYPE	APPS	REG	FORMAT	SIZE	CAT. NO.
basic FGF	Chick, Dog, Hu, Ms, Rat	3	Mouse IgG ₁	IF, WB	RUO	Purified	50 µg	610870
				IF, WB	RUO	Purified	150 µg	610871
				WB	RUO	HRP	150 µg	610873
		6	Mouse IgG ₁	IF, IHC, IP, WB	RUO	Purified	50 µg	610072
				IF, IHC, IP, WB	RUO	Purified	150 µg	610073
EGF Receptor	Hu, Ms	13	Mouse IgG ₁	IF, IP, WB	RUO	Purified	50 µg	610016
				IF, IP, WB	RUO	Purified	150 µg	610017
GM-CSF	Ms	MP1-22E9	Rat IgG _{2a}	Neu	RUO	NA/LE	0.5 mg	554403
				ELISA Cap., FCM, WB	RUO	Purified	0.5 mg	554404
				IC/FCM	RUO	PE	0.1 mg	554406
		MP1-31G6	Rat IgG ₁	ELISA Det., IP/WB	RUO	Biotin	0.5 mg	554407
IL-3	Ms	MP2-43D1	Rat IgG _{2a}	ELISA Det., WB	RUO	Biotin	0.5 mg	554384
				Neu	RUO	NA/LE	0.5 mg	554379
		MP2-8F8	Rat IgG ₁	ELISA Cap., IC/FCM, Block	RUO	Purified	0.5 mg	554381
				IC/FCM	RUO	PE	0.1 mg	554383
IL-6	Ms	MP5-20F3	Rat IgG ₁	Neu	RUO	NA/LE	0.5 mg	554398
				ELISA Cap., IC/FCM, Block, WB	RUO	Purified	0.5 mg	554400
				IC/FCM	RUO	PE	0.1 mg	554401
		MP5-32C11	Rat IgG _{2a}	ELISA Det., WB	RUO	Biotin	0.5 mg	554402

Reagents for Mouse Stem Cell Research

Antibodies to Stem Cell Markers

DESCRIPTION	REACT	CLONE	ISOTYPE	APPS	REG	FORMAT	SIZE	CAT. NO.		
CD4 (L3T4)	Ms	GK1.5	Rat (Lewis) IgG _{2b} , κ	Block, Cost, Cyt, FCM	RUO	NA/LE	0.5 mg	553726		
				FCM, IP	RUO	Purified	0.5 mg	553727		
				FCM	RUO	Biotin	0.5 mg	553728		
				FCM	RUO	FITC	0.1 mg	557307		
				FCM	RUO	FITC	0.5 mg	553729		
				FCM	RUO	PE	0.1 mg	557308		
				FCM	RUO	PE	0.2 mg	553730		
				FCM	RUO	APC-Cy7	0.1 mg	552051		
				H129.19	Rat (Lou/WS1) IgG _{2b} , κ	Block, Cyt, FCM	RUO	NA/LE	0.5 mg	553646
						FCM, IP	RUO	Purified	0.5 mg	553647
						IHC(Fr), IHC(Zn)	RUO	Purified	1 ml	550278
						FCM	RUO	Biotin	0.1 mg	553648
						FCM	RUO	Biotin	0.5 mg	553649
						FCM	RUO	FITC	0.1 mg	553650
						FCM	RUO	FITC	0.5 mg	553651
		FCM	RUO			PE	0.1 mg	553652		
		FCM	RUO			PE	0.2 mg	553653		
		RM4-4	Rat (SD) IgG _{2b} , κ	FCM	RUO	PE-Cy5	0.1 mg	553654		
				FCM	RUO	Purified	0.5 mg	553053		
				FCM	RUO	Biotin	0.1 mg	557443		
		RM4-5	Rat (DA) IgG _{2b} , κ	FCM	RUO	FITC	0.5 mg	553055		
				FCM, IP	RUO	Purified	0.5 mg	553043		
				IHC(Fr), IHC(Zn)	RUO	Purified	1 ml	550280		
				FCM, Sep	RUO	Biotin	0.1 mg	553044		
				FCM, Sep	RUO	Biotin	0.5 mg	553045		
				FCM	RUO	FITC	0.1 mg	553046		
				FCM	RUO	FITC	0.5 mg	553047		
FCM*	RUO			PE	0.1 mg	553048				
FCM*	RUO			PE	0.2 mg	553049				
FCM	RUO			PerCP	0.1 mg	553052				

Product Listing *(continued)*

Reagents for Mouse Stem Cell Research

Antibodies to Stem Cell Markers (continued)

DESCRIPTION	REACT	CLONE	ISOTYPE	APPS	REG	FORMAT	SIZE	CAT. NO.
				FCM	RUO	PerCP-Cy5.5	0.1 mg	550954
				FCM	RUO	APC	0.1 mg	553051
				FCM	RUO	Alexa Fluor® 488	0.1 mg	557667
				FCM	RUO	Alexa Fluor® 647	0.1 mg	557681
				FCM	RUO	Alexa Fluor® 700	0.1 mg	557956
				FCM	RUO	PE-Cy7	0.1 mg	552775
				FCM	RUO	PE-Cy5	0.1 mg	553050
				FCM	RUO	Pacific Blue®	0.1 mg	558107
CD24 (Heat Stable Antigen)	Ms	30-F1	Rat (LOU/Ws1/M) IgG _{2c} , κ	FCM	RUO	Purified	0.1 mg	558777
				FCM	RUO	Biotin	0.5 mg	555296
		J11d	Rat (Lewis) IgM, κ	FCM, IHC(Fr), WB	RUO	Purified	0.5 mg	553146
		M1/69	Rat (DA) IgG _{2b} , κ	Cyt, FCM, IHC(F), IHC(Fr), IHC(Zn), WB	RUO	Purified	0.1 mg	557436
				FCM	RUO	Biotin	0.5 mg	553260
				FCM	RUO	FITC	0.5 mg	553261
				FCM*	RUO	PE	0.1 mg	553262
CD30	Ms	mCD30.1	Armenian Hamster IgG ₁ , κ	FCM	RUO	Biotin	0.5 mg	553824
				FCM*	RUO	PE	0.2 mg	553825
CD34	Ms	RAM34	Rat IgG _{2a} , κ	FCM	RUO	Purified	0.5 mg	553731
				IHC(Fr), IHC(Zn)	RUO	Purified	1 ml	550537
				FCM	RUO	Biotin	0.5 mg	553732
				FCM	RUO	FITC	0.5 mg	553733
				FCM	RUO	PE	0.2 mg	551387
				FCM	RUO	APC	0.1 mg	558038
CD38	Ms	90	Rat IgG _{2a} , κ	FCM	RUO	Biotin	0.5 mg	553762
				FCM	RUO	FITC	0.1 mg	558813
				FCM*	RUO	PE	0.2 mg	553764
CD45 (Leukocyte Common Antigen)	Hu, Ms, Rat	69	Mouse IgG ₁ , κ	IF, IP, WB	RUO	Purified	150 µg	610266
CD45 (Leukocyte Common Antigen, Ly-5)	Ms	30-F11	Rat (LOU/Ws1/M) IgG _{2b} , κ	Cyt, FCM, IP	RUO	Purified	0.5 mg	553076
				IHC(F), IHC(Fr), IHC(Zn)	RUO	Purified	1 ml	550539
				FCM	RUO	Biotin	0.1 mg	553077
				FCM	RUO	Biotin	0.5 mg	553078
				FCM	RUO	FITC	0.1 mg	553079
				FCM	RUO	FITC	0.5 mg	553080
				FCM	RUO	PE	0.1 mg	553081
				FCM	RUO	PerCP	0.1 mg	557235
				FCM	RUO	PerCP-Cy5.5	0.1 mg	550994
				FCM	RUO	APC	0.1 mg	559864
				FCM	RUO	APC-Cy7	0.1 mg	557659
				FCM	RUO	PE-Cy5	0.1 mg	553082
				FCM	RUO	PE-Cy7	0.1 mg	552848
				CD56 (N-CAM)	Chick, Hu, Ms, Rat	12F11	Rat IgG _{2a}	IHC(Fr), WB
Hu, Ms, Rat	12F8	Rat IgM	IHC(Fr), WB		RUO	Purified	0.1 mg	556325
Ms, Rat	N-CAM 13	Mouse IgG _{2a}	IHC(Fr), WB		RUO	Purified	0.1 mg	556324
CD71 (Transferrin Receptor)	Ms	C2	Rat (Wistar/Furth) IgG ₁ , κ	FCM, IF, IHC(Fr)	RUO	Purified	0.5 mg	553264
				FCM	RUO	Biotin	0.1 mg	557416
				FCM	RUO	FITC	0.5 mg	553266
				FCM	RUO	PE	0.1 mg	553267
CD90 (Thy-1)	Ms	G7	Rat (Lewis) IgG _{2c} , κ	Cost, FCM	RUO	NA/LE	0.5 mg	553015
				FCM, IHC(Fr), IHC(Zn)	RUO	Purified	0.5 mg	553016
CD90.1 (Thy-1.1)	G Pig, Ms, Rab, Rat	OX-7	Mouse (BALB/c) IgG ₁ , κ	EM, FCM, IHC(Fr), IHC(Zn), IP, WB	RUO	Purified	0.5 mg	554895
				IHC(Fr), IHC(Zn)	RUO	Purified	1 ml	550571
				FCM	RUO	Biotin	0.5 mg	554896

Reagents for Mouse Stem Cell Research

Antibodies to Stem Cell Markers

DESCRIPTION	REACT	CLONE	ISOTYPE	APPS	REG	FORMAT	SIZE	CAT. NO.
				FCM	RUO	FITC	0.5 mg	554897
				FCM	RUO	PE	0.1 mg	551401
				FCM	RUO	PE	0.2 mg	554898
				FCM	RUO	PerCP	0.1 mg	557266
CD90.2 (Thy-1.2)	Ms	30-H12	Rat (LOU/Ws1/M) IgG _{2b} κ	Cost, Cyt, FCM	RUO	NA/LE	0.5 mg	553008
				Cyt, FCM	RUO	Purified	0.5 mg	553009
				FCM, Sep	RUO	Biotin	0.1 mg	553010
				FCM, Sep	RUO	Biotin	0.5 mg	553011
				FCM	RUO	FITC	0.1 mg	553012
				FCM	RUO	FITC	0.5 mg	553013
				FCM	RUO	PE	0.1 mg	553014
		53-2.1	Rat (LOU/Ws1/M) IgG _{2a} κ	Cyt, EM, FCM, IP	RUO	Purified	0.5 mg	553000
				IHC(Fr)	RUO	Purified	1 ml	550543
				FCM	RUO	Biotin	0.1 mg	553001
				FCM	RUO	Biotin	0.5 mg	553002
				FCM	RUO	FITC	0.1 mg	553003
				FCM	RUO	FITC	0.5 mg	553004
				FCM	RUO	PE	0.1 mg	553005
FCM	RUO	PE	0.2 mg	553006				
FCM	RUO	APC	0.1 mg	553007				
CD105 (Endoglin)	Ms	MJ7/18	Rat IgG _{2a} κ	IHC(Fr)	RUO	Purified	1 ml	550546
CD117 (c-kit)	Ms	2B8	Rat (Wistar) IgG _{2b} κ	FCM, IP	RUO	Purified	0.5 mg	553352
				FCM	RUO	Biotin	0.5 mg	553353
				FCM	RUO	FITC	0.5 mg	553354
				FCM	RUO	PE	0.2 mg	553355
				FCM	RUO	APC	0.1 mg	553356
				FCM	RUO	PE-Cy7	0.1 mg	558163
		ACK45	Rat IgG _{2b} κ	Block, FCM	RUO	NA/LE	0.5 mg	553867
				FCM	RUO	Purified	0.5 mg	553868
				FCM	RUO	PE	0.2 mg	553869
				FCM	RUO	PE	0.2 mg	553869
CD127 (IL-7 Receptor α chain)	Ms	B12-1	Rat (Wistar) IgG _{2a} λ	FCM	RUO	Biotin	0.5 mg	555288
				Block, FCM	RUO	NA/LE	0.5 mg	550426
		SB/199	Rat IgG _{2b} κ	FCM	RUO	PE	0.2 mg	552543
				FCM	RUO	PE	0.2 mg	552543
CD144 (VE-Cadherin)	Ms	11D4.1	Rat (Lewis) IgG _{2a} κ	FCM, IP	RUO	Purified	0.5 mg	555289
				IHC(Fr)	RUO	Purified	1 ml	550548
CD153 (CD30 Ligand)	Ms	RM153	Rat (Sprague-Dawley) IgG _{2b} κ	FCM	RUO	PE	0.1 mg	559232
CD184 (CXCR4)	Ms	2B11/CXCR4	Rat IgG _{2b}	FCM	RUO	Purified	0.1 mg	551852
				FCM	RUO	Biotin	0.1 mg	551968
				FCM	RUO	FITC	0.1 mg	551967
				FCM	RUO	PE	0.1 mg	551966
Early B Lineage	Ms	493	Rat (Lewis) IgG _{2a} κ	FCM, IP	RUO	Purified	0.1 mg	550433
				FCM	RUO	Biotin	0.1 mg	550434
		AA4.1	Rat (Sprague-Dawley) IgG _{2b} κ	FCM	RUO	FITC	0.5 mg	559156
				FCM	RUO	PE	0.1 mg	558039
GFAP	Bov, Chick, Dog, G Pig, Hu, Ms, Pig, Rab, Rat, Sheep	2.00E+01	Mouse IgG _{2b}	IF, IHC(F), WB	RUO	Purified	0.5 mg	556329
		4A11	Mouse IgG _{2b}	IF, IHC(F), WB	RUO	Purified	0.5 mg	556327
Lineage Panel	Ms	145-2C11, M1/70, RA3-6B2, RB6-8C5, TER-119		FCM, Sep	RUO	Biotin	1,000 tests	559971
Lineage Cocktail, with Isotype Control	Ms	145-2C11, M1/70, RA3-6B2, RB6-8C5, TER-119		FCM	RUO	APC	100 tests	558074

Product Listing *(continued)*

Reagents for Mouse Stem Cell Research

Antibodies to Stem Cell Markers (continued)

DESCRIPTION	REACT	CLONE	ISOTYPE	APPS	REG	FORMAT	SIZE	CAT. NO.	
Ly-6A/E (Sca-1)	Ms	D7	Rat (Lewis) IgG _{2a} , κ	FCM, IHC(Fr), IP, WB	RUO	Purified	0.1 mg	557403	
				FCM	RUO	Biotin	0.1 mg	557404	
				FCM	RUO	FITC	0.1 mg	557405	
				FCM	RUO	PE	0.2 mg	553108	
			E13-161.7	Rat (Wistar) IgG _{2a} , κ	Cyt, FCM, IHC(Fr), IP	RUO	PE-Cy7	0.1 mg	558162
					FCM	RUO	Biotin	0.5 mg	553333
					FCM	RUO	FITC	0.5 mg	553335
					FCM	RUO	PE	0.2 mg	553336
Notch1	Hu, Ms	mN1A	Mouse IgG ₁ , κ	ELISA, IC/FCM, IHC, WB	RUO	Purified	0.1 mg	552466	
				IC/FCM	RUO	PE	0.1 mg	552768	
3-Oct	Ms	40	Mouse IgG ₁	WB	RUO	Purified	50 µg	611202	
				WB	RUO	Purified	150 µg	611203	
Tyrosine Hydroxylase	Ms, Rat	45	Mouse IgG ₁	WB	RUO	Purified	50 µg	612300	
				WB	RUO	Purified	150 µg	612301	

Reagents for Mouse Stem Cell Research

Magnetic Cell Separation Reagents

DESCRIPTION	REACT	CLONE	ISOTYPE	APPS	REG	FORMAT	SIZE	CAT. NO.
Hematopoietic Progenitor Cell Enrichment Set - DM	Ms			Sep	RUO	BD IMag-DM, Biotin	2 x 10 ⁹ cells	558451

Reagents for Mouse Stem Cell Research

Stem Cell Growth Factors

DESCRIPTION	REACT	CLONE	ISOTYPE	APPS	REG	FORMAT	SIZE	CAT. NO.
GM-CSF, Recombinant mouse	Ms			ELISA Std., FA, IC/FCM Block	RUO		10 µg	554586
IL-3, Recombinant mouse	Ms			ELISA Std., FA, IC/FCM Block	RUO		10 µg	554579
IL-6, Recombinant mouse	Ms			ELISA Std., FA, IC/FCM Block	RUO		5 µg	554582
IL-7, Recombinant mouse	Ms			FA	RUO		2 µg	554583

Other Reagents for Stem Cell Research

Extracellular Matrix Molecules and Growth Factors

DESCRIPTION	SIZE	CAT. NO.
Lymphokine Products		
<i>Colony-Stimulating Factors</i>		
GM-CSF, human recombinant	1 µg	354048
Interleukins		
IL-1β, human recombinant	2 µg	354042
IL-2, human recombinant	10,000 µg	354043
IL-2, human recombinant	50,000 µg	356043
IL-2, mouse recombinant	10,000 µg	356078
IL-2, mouse recombinant	25,000 µg	354078
IL-2, rat natural	4,000 µg	354110
IL-3, mouse recombinant	10 µg	354058
IL-4, human recombinant	5 µg	354068
Stem Cell Factor		
SCF, human recombinant	10 µg	354105
Tumor Necrosis Factor		
TNF-β, human recombinant	10 µg	354066
TNF-β, human recombinant	5x10 µg	356066
Growth Factor Products		
<i>Epidermal Growth Factors</i>		
EGF, mouse natural, (culture grade)	100 µg	354001
EGF, mouse natural, (culture grade) (1 mg)	10x100 µg	356001
EGF, mouse natural, (receptor grade)	100 µg	354010
EGF, mouse natural, (receptor grade) (0.5 mg)	5x100 µg	356010
EGF, human recombinant	100 µg	354052
EGF, human recombinant (1 mg)	10x100 µg	356052
<i>Fibroblast Growth Factors</i>		
bFGF, bovine natural	10 µg	356037
bFGF, human recombinant	10 µg	354060
bFGF, human recombinant (50 µg)	5x10 µg	356060
bFGF, human recombinant (100 µg)	10x10 µg	356061
<i>Hepatocyte Growth Factor</i>		
HGF/SF, human recombinant	5 µg	354103
<i>Insulin-Like Growth Factors</i>		
IGF-1, human recombinant (culture grade)	10 µg	354037
<i>Nerve Growth Factors</i>		
2.5S NGF, mouse natural	10 µg	354005
2.5S NGF, mouse natural	100 µg	356004
2.5S NGF, mouse natural (1 mg)	2x500 µg	356005
7S NGF, mouse natural	100 µg	354009
<i>Platelet-Derived Growth Factors</i>		
PDGF-BB, human recombinant	10 µg	354051
PDGF-BB, human recombinant (0.1 mg)	10x10 µg	356051
<i>Transforming Growth Factors</i>		
TGF-β, human natural	1 µg	354039
TGF-β, human natural (5 µg)	5x1 µg	356039
TGF-β, human natural (10 µg)	5x2 µg	356040
<i>Vascular Endothelial Growth Factor</i>		
VEGF, human recombinant	10 µg	354107
Extracellular Matrix Products		
<i>BD Matrigel™ Basement Membrane Matrix</i>		
BD Matrigel Matrix	5 ml	356234
BD Matrigel Matrix	10 ml	354234
BD Matrigel Matrix (5x10 ml)	50 ml	356235
BD Matrigel Matrix Phenol-Red Free	10 ml	356237


DESCRIPTION	SIZE	CAT. NO.
Extracellular Matrix Products (continued)		
<i>BD Matrigel™ Matrix</i>		
BD Matrigel Matrix High Concentration (HC)	10ml	354248
BD Matrigel Matrix High Concentration (HC) Phenol-Red Free	10ml	354262
BD Matrigel Matrix High Concentration (HC) Growth Factor Reduced	10ml	354263
Growth Factor Reduced (GFR) BD Matrigel Matrix	5 ml	356230
GFR BD Matrigel Matrix	10 ml	354230
GFR BD Matrigel Matrix Phenol-Red Free	10 ml	356231
Dispase	100 ml	354235
BD Cell Recovery Solution	100 ml	354253
Extracellular Matrix, human	1 mg	354237
<i>BD Cell-Tak™ Cell and Tissue Adhesive</i>		
BD Cell-Tak Cell and Tissue Adhesive	1 mg	354240
BD Cell-Tak Cell and Tissue Adhesive	5 mg	354241
BD Cell-Tak Cell and Tissue Adhesive (2x5 mg)	10 mg	354242
<i>Collagens</i>		
Collagen I, bovine	30 mg	354231
Collagen I, human	0.25 mg	354243
Collagen I, rat tail	100 mg	354236
Collagen I, rat tail (10x100 mg)	1 g	356236
Collagen I, human recombinant	250 µg	354254
Collagen II, bovine	5 mg	354257
Collagen III, human	0.25 mg	354244
Collagen I High concentration rat tail	100mg	354249
Collagen IV, human	0.25 mg	354245
Collagen IV, mouse	1 mg	354233
Collagen IV, mouse (10x1 mg)	10 mg	356233
Collagen V, human	0.25 mg	354246
Collagen VI, Human	500 µg	354261
<i>Fibronectin</i>		
Fibronectin, human	1 mg	354008
Fibronectin, human	5 mg	356008
Fibronectin, human (5x5 mg)	25 mg	356009
<i>Laminin</i>		
Laminin, mouse	1 mg	354232
Ultra-pure Laminin, mouse	1 mg	354239
Laminin/Entactin High Concentration	500 µg	354259
<i>Poly-D-Lysine</i>		
Poly-D-Lysine	20 mg	354210
<i>BD™ PuraMatrix™ Peptide Hydrogel</i>		
BD™ PuraMatrix™ Peptide Hydrogel	5 ml	354250
<i>Vitronectin</i>		
Vitronectin, human	0.25 mg	354238
<i>BD™ 3D Scaffold products</i>		
BD 3D Collagen Composite scaffold	24 scaffolds in one BD Falcon™ 48-well plate	354613
BD 3D OPLA® scaffold	24 scaffolds in one BD Falcon™ 48-well plate	354614
BD 3D Calcium Phosphate scaffold	24 scaffolds in one BD Falcon™ 48-well plate	354617

Other Reagents for Stem Cell Research *(continued)*

IVD Reagents

For In Vitro Diagnostic Use.

DESCRIPTION	REACT	APPS	REG	SIZE	CAT. NO.
BD Procount Progenitor Cell Enumeration Kit	Hu	FCM	IVD	25 tests	340498
BD Stem Cell Control Kit (Bi-Level Control)	Hu	FCM	IVD	15 tests	340991
BD Procount v2.1 (OS9 and G4 computers or earlier only) for Macintosh		FCM	IVD		341158



References

1. Prockop, D.J., Marrow stromal cells as stem cells for nonhematopoietic tissues. *Science* 276:71-74 (1997).
2. Orkin, S.H., Diversification of hematopoietic stem cells to specific lineages. *Nat. Rev. Genet.* 1:57-64 (2000).
3. Odorico, J.S., Kaufman, D.S., Thomson, J.A., Multilineage differentiation from human embryonic stem cell lines. *Stem Cells* 19:193-204 (2001).
4. Blau, H.M., Brazelton, T.R., Weimann, J.M., The evolving concept of a stem cell: entity or function? *Cell* 105:829-841 (2001).
5. Watt, F.M., Hogan, B.L., Out of Eden: stem cells and their niches. *Science* 287:1427-1430 (2000).
6. Watt, F.M., Stem cell fate and patterning in mammalian epidermis. *Curr. Opin. Genet. Dev.* 11:410-417 (2001).
7. Panchision, D.M., McKay, R.D., The control of neural stem cells by morphogenic signals. *Curr. Opin. Genet. Dev.* 12:478-487 (2002).
8. Lumelsky, N., Blondel, O., Laeng, P., Velasco, I., Ravin, R., McKay, R., Differentiation of embryonic stem cells to insulin-secreting structures similar to pancreatic islets. *Science* 292:1389-1394 (2001).
9. Lee, S.H., Lumelsky, N., Studer, L., Auerbach, J.M., McKay, R.D., Efficient generation of midbrain and hindbrain neurons from mouse embryonic stem cells. *Nat. Biotechnol.* 18:675-679 (2000).
10. Bain, G., Kitchens, D., Yao, M., Huettner, J.E., Gottlieb, D.I., Embryonic stem cells express neuronal properties *in vitro*. *Dev. Biol.* 168, 342-357 (1995).
11. Wichterle H., Lieberam I., Porter J.A., Jessell T.M., Directed differentiation of embryonic stem cells into motor neurons. *Cell* 110:385-397 (2002).
12. Brustle, O., Jones, K.N., Learish, R.D., Karram, K., Choudhary, K., Wiestler, O.D., Duncan, I.D., McKay, R.D., Embryonic stem cell-derived glial precursors: a source of myelinating transplants. *Science* 285: 754-756 (1999).
13. Dani, C., Smith, A.G., Dessolin, S., Leroy, P., Staccini L., Villageois, P., Darimont, C., Ailhaud, G., Differentiation of embryonic stem cells into adipocytes *in vitro*. *J. Cell Sci.* 110:1279-1285 (1997).
14. Kramer, J., Hegert, C., Guan, K., Wobus, A.M., Muller, P.K., Rohwedel, J., Embryonic stem cell-derived chondrogenic differentiation *in vitro*: activation by BMP-2 and BMP-4. *Mech. Dev.* 92:193-205 (2000).
15. Buttery, L.D., Bourne, S., Xynos, J.D., Wood, H., Hughes, F.J., Hughes, S.P., Episkopou, V., Polak, J.M., Differentiation of osteoblasts and *in vitro* bone formation from murine embryonic stem cells. *Tissue Eng.* 7:89-99 (2001).
16. Fairchild, P.J., Brook, F.A., Gardner, R.L., Graca, L., Strong, V., Tone, Y., Tone, M., Nolan, K.F., Waldmann, H., Directed differentiation of dendritic cells from mouse embryonic stem cells. *Curr. Biol.* 10:1515-1518 (2000).
17. Yamashita, J., Itoh, H., Hirashima, M., Ogawa, M., Nishikawa, S., Yurugi, T., Naito, M., Nakao, K., Nishikawa, S., Flk1-positive cells derived from embryonic stem cells serve as vascular progenitors. *Nature* 408:92-96 (2000).
18. Hanada, K., Dennis, J.E., Caplan, A.L., Stimulatory effects of basic fibroblast growth factor and bone morphogenetic protein-2 on osteogenic differentiation of rat bone marrow-derived mesenchymal stem cells. *J. Bone Miner. Res.* 12:1606-1614 (1997).
19. Pittenger, M.F., Mackay, A.M., Beck, S.C., Jaiswal, R.K., Douglas, R., Mosca, J.D., Moorman, M.A., Simonetti, D.W., Craig, S., Marshak, D.R., Multilineage potential of adult human mesenchymal stem cells. *Science* 284:143-147 (1999).
20. Colter, D.C., Sekiya, I., Prockop, D.J., Identification of a subpopulation of rapidly self-renewing and multipotential adult stem cells in colonies of human marrow stromal cells. *Proc. Natl. Acad. Sci. U S A* 98:7841-7845 (2001).
21. Mackay, A.M., Beck, S.C., Murphy, J.M., Barry, F.P., Chichester, C.O., Pittenger, M.F., Chondrogenic differentiation of cultured human mesenchymal stem cells from marrow. *Tissue Eng.* 4:415-428 (1998).
22. Sanchez-Ramos, J., Song, S., Cardozo-Pelaez, F., Hazzi, C., Stedeford, T., Willing, A., Freeman, T.B., Saporta, S., Janssen, W., Patel, N., Cooper, D.R., Sanberg, P.R., Adult bone marrow stromal cells differentiate into neural cells *in vitro*. *Exp. Neurol.* 164:247-256 (2000).
23. Reyes, M., Lund, T., Lenvik, T., Aguiar, D., Koodie, L., Verfaillie, C.M., Purification and *ex vivo* expansion of postnatal human marrow mesodermal progenitor cells. *Blood* 98:2615-2625 (2001).
24. Reyes, M., Dudek, A., Jahagirdar, B., Koodie, L., Marker, P.H., Verfaillie, C.M., Origin of endothelial progenitors in human postnatal bone marrow. *J. Clin. Invest.* 109:337-346 (2002).
25. Schwartz, R.E., Reyes, M., Koodie, L., Jiang, Y., Blackstad, M., Lund, T., Lenvik, T., Johnson, S., Hu, W.S., Verfaillie, C.M., Multipotent adult progenitor cells from bone marrow differentiate into functional hepatocyte-like cells. *J. Clin. Invest.* 109:1291-1302 (2002).
26. Morrison, S.J., Uchida, N., Weissman, I.L., The biology of hematopoietic stem cells. *Annu. Rev. Cell Dev. Biol.* 11:35-71 (1995).
27. Yang, L., Li, S., Hatch, H., Ahrens, K., Cornelius, J.G., Petersen, B.E., Peck, A.B., *In vitro* trans-differentiation of adult hepatic stem cells into pancreatic endocrine hormone-producing cells. *Proc. Natl. Acad. Sci. USA* 99:8078-8083 (2002).
28. Petersen, B.E., personal communication.
29. Carpenter, M.K., Cui, X., Hu, Z.Y., Jackson, J., Sherman, S., Seiger, A., Wahlberg, L.U., *In vitro* expansion of a multipotent population of human neural progenitor cells. *Exp. Neurol.* 158:265-278 (1999).
30. Palmer, T.D., Schwartz, P.H., Taupin, P., Kaspar, B., Stein, S.A., Gage, F.H., Progenitor cells from human brain after death. *Nature* 411:42-43 (2001).
31. Reynolds, B.A., Weiss, S., Clonal and population analyses demonstrate that an EGF-responsive mammalian embryonic CNS precursor is a stem cell. *Dev. Biol.* 175:1-13 (1996).
32. Qian, X., Shen, Q., Goderie, S.K., He, W., Capela, A., Davis, A.A., Temple, S., Timing of CNS cell generation: a programmed sequence of neuron and glial cell production from isolated murine cortical stem cells. *Neuron*. 28:69-80 (2000).
33. Caldwell, M.A., He, X., Wilkie, N., Pollack, S., Marshall, G., Wafford, K.A., Svendsen, C.N., Growth factors regulate the survival and fate of cells derived from human neurospheres. *Nat. Biotechnol.* 19:475-479 (2001).
34. Johe, K.K., Hazel, T.G., Muller, T., Dugich-Djordjevic, M.M., McKay, R.D., Single factors direct the differentiation of stem cells from the fetal and adult central nervous system. *Genes Dev.* 10:3129-3140 (1996).
35. Zuk, P.A., Zhu, M., Mizuno, H., Huang, J., Futrell, J.W., Katz, A.J., Benhaim, P., Lorenz, H.P., Hedrick, M.H., Multilineage cells from human adipose tissue: implications for cell-based therapies. *Tissue Eng.* 7:211-228 (2001).
36. Erickson, G.R., Gimble, J.M., Franklin, D.M., Rice, H.E., Awad, H., Guilak, F., Chondrogenic potential of adipose tissue-derived stromal cells *in vitro* and *in vivo*. *Biochem. Biophys. Res. Commun.* 290:763-769 (2002).
37. Deasy, B.M., Qu-Peterson, Z., Greenberger, J.S., Huard, J., Mechanisms of muscle stem cell expansion with cytokines. *Stem Cells* 20:50-60 (2002).
38. Lee, J.Y., Qu-Petersen, Z., Cao, B., Kimura, S., Jankowski, R., Cummins, J., Usas, A., Gates, C., Robbins, P., Wernig, A., Huard, J., Clonal isolation of muscle-derived cells capable of enhancing muscle regeneration and bone healing. *J. Cell Biol.* 150:1085-1100 (2000).
39. Semino, C.E., et al., Functional differentiation of hepatocyte-like spheroid structures from putative liver progenitor cells in threedimensionalpeptide scaffolds. *Differentiation* 71:262 (2003).
40. Semino, C.E., et al., Entrapment of migrating hippocampal neural cells in 3D peptide nanofiber scaffold. *Tissue Engineering* 10:643(2004).
41. Xu, C. et al *Nature Biotechnology*, Vol 19, pp 971-974. (2001)
42. Zwaka, T and Thomson, J.A. *Nature Biotechnology*, Vol 21, pp 319-321. (2003)
43. Unpublished results
44. Holmes, TC, et al. *PNAS USA* 97:6728 (2000).
45. Narmoneva, DA, et al. *Circulation* 110:962 (2004).
46. Davis, M.E., et al. *Circulation* 111:442 (2005).
47. Civin, C., Straus, L.C., Brovall, C., Fackler, M.J., Schwartz, J.F., and Shaper, J.H., *J. Immunol.* 133, 157 (1984).
48. Baum, C.M., Weissman, I.L., Tsukamoto, A., Buckle, A., Peault, B., *Proc. Natl. Acad. Sci. USA* 89:2804 (1992).
49. Spangrude, G.J., Heimfeld, S., and Weissman, I., *Science* 241, 58 (1988).
50. Srour, E.F., Leemhuis, T., Brandt, J.E., van Beisen, K., and Hoffman, R., *Cytometry* 12:179 (1991).
51. Uchida, N., Combs, J., Chen, S., Zanjani, E., Hoffman, R., Tsukamoto, A., *Blood* 88(4):1297 (1996).
52. Ladd, A.C., Pyatt, R., Gothot, A., Rice, S., McMahl, J., Traycoff, C.M., and Srour, E., *Blood* 90(2):658 (1997).
53. Goodell, M.A., Rosenzweig, M., Kim, H., Marks, D.F., DeMaria, M., Paradis, G., Grupp, S.A., Sief, C.A., Mulligan, R.C., and Johnson, R.P., *Nat. Med.* 3(12): 1337 (1997).
54. Storms, R.W., Goodell, M.A., Fisher, A., Mulligan, R.C., and Smith, C., *Blood* 96(6):2125 (2000).
55. Wulf, G.G., Wang, R., Kuehnle, I., Weider, D., Marini, F., Brenner, M., Andreef, M., and Goodell, M.A., *Blood* 98(4):1166 (2001).
56. Protocol from: www.bcm.tmc.edu/genetherapy/goodell
57. Watson, J.V., Hakeff, A., Chambers, S.H., and Smith, P.J., *Cytometry* 6:310 (1985).

Additional Resources

bdbiosciences.com/stemcellsource

Listed below are additional resources that can be ordered free of charge or easily downloaded from our website at bdbiosciences.com



BD Procount™
Setting the standard for progenitor cell enumeration. Those of you interested in a standardized, accurate assay approach to counting CD34+ cells will find the BD Procount white paper reports an effective and more convenient method for progenitor cell enumeration.



BD Oncomark™
Multicolor reagents designed for your clinical testing.



BD Oncomark™ Poster
The evolution of blood cells.




BD IMag™ Cell Separation System
The BD IMag™ Cell Separation System is based on a simple yet highly effective direct magnet technology that allows for the rapid positive and/or negative selection (enrichment) of cell populations without the use of magnetic separation columns.




BD™ PuraMatrix™ Peptide Hydrogel
BD™ PuraMatrix™ Peptide Hydrogel, a novel synthetic matrix ideal for creating optimized 3D cell culture environments.




Tools for Stem Cell Biology and Tissue Engineering Research
BD™ Extracellular Matrix Proteins, Growth Factors, and Three Dimensional Scaffolds.



Driving Cellular Communication
BD™ Extracellular Matrix Proteins.



BD BioCoat™ Cellware
BD Falcon™ Tissue Culture Plates, Dishes, Flasks, Coverslips, and Culture Slides are available pre-coated with a variety of ECM proteins and attachment factors.



BD Falcon™ Cell Culture Products
BD Falcon™ Cell Culture Plates, Dishes, Flasks, Roller Bottles and CultureSlides.

Regional Offices

Asia Pacific
BD Singapore
Tel 65.6861.0633
Fax 65.6860.1593

Australia/New Zealand
Australia
Toll Free: 1800 656 100
Tel 61.2.8875.7000
Fax 61.2.8875.7200
bd_anz@bd.com

New Zealand
Toll Free: 0800 572.468
Tel 64.9.574.2468
Fax 64.9.574.2469
bd_anz@bd.com

Europe
Belgium
Tel 32.53.720.550
Fax 32.53.720.549
customer_service_bdbelgium@europe.bd.com

Canada
BD Biosciences
Toll Free 888.259.0187
Tel 905.542.8028
Fax 888.229.9918
canada@bd.com

Japan
Nippon Becton Dickinson
Toll Free 0120.8555.90
Tel 81.24.593.5405
Fax 81.24.593.5761

United States
BD Biosciences
Toll Free 877.232.8995
Customer/Technical Service
Tel 301.340.9775
Bioimaging Systems
Tel 301.340.9775
Discovery Labware
Tel 978.901.7491
Immuncytometry Systems
Tel 800.325.9637
Pharming
Tel 800.325.9637
bdbiosciences.com

Local Offices and Distributors

Argentina/Paraguay/Uruguay
Toll Free: 0800.444.5523
Tel 54.11.4551.7100
Fax 54.11.4551.7400
crc_argentina@bd.com

Austria
Tel 43.1.706.36.60.20
Tel 43.1.706.36.60.11
customerservice.bdb.at@europe.bd.com

Brazil
Unidade Biosciences
Tel 55.11.5185.9995
Fax 55.11.5185.9895
biosciences@bd.com.br

Central America/Caribbean
Tel 502.2.440.4195
Fax 502.2.440.3570

Chile
Tel 56.2.482.7800 x16
Fax 56.2.482.7802

China
Tel 8610.6418.1608
Fax 8610.6418.1610

Colombia
Tel 57.1.572.4060 x244 / 258
Fax 57.1.572.3650 / 4250

Denmark
Tel 45.43.43.45.66
Fax 45.43.43.41.66
bdbnordic@europe.bd.com

East Africa
Tel 254.2.341157
Fax 254.2.341161
bd@africaonline.co.ke

Eastern Europe
Tel 49.6221.305.161
Fax 49.6221.305.418
bdb.ema@europe.bd.com

Egypt
Tel 202.268.0181
Fax 202.266.7562

Finland
Tel 358.9.88.70.7832
Fax 358.9.88.70.7817
bdbnordic@europe.bd.com

France
Tel 33.4.76.68.36.09
Fax 33.4.76.68.35.06
customerservice.bdb.france@europe.bd.com

Germany
Tel 49.6221.305.551
Fax 49.6221.303.609
customerservice.bdb.de@europe.bd.com

Greece
Tel 30.1.940.77.41
Fax 30.1.940.77.40

Hong Kong
Tel 852.2575.8668
Fax 852.2803.5320

Hungary
Tel 36.1.345.7090
Fax 36.1.345.7093
bdb.ema@europe.bd.com

India
Tel 91.124.238.3566.70
Fax 91.124.238.3225

Indonesia
Tel 62.21.577.1920
Fax 62.21.577.1925

Italy
Tel 39.02.48.240.1
Fax 39.02.48.20.33.36

Korea
Tel 822.3404.3700
Fax 822.557.4048

Malaysia
Tel 603.7725.5517
Fax 603.7725.4772

Mexico
Tel 52.55.5999.8296
Fax 52.55.5999.8288
bioscienciasmexico@bd.com

Middle East
Tel 971.4.337.95.25
Fax 971.4.337.95.51
bdb.ema@europe.bd.com

The Netherlands
Tel 31.20.582.94.20
Fax 31.20.582.94.21
customer_service_bdholland@europe.bd.com

North Africa
Tel 33.4.76.68.35.03
Fax 33.4.76.68.35.44
bdbiosciences_maghreb@europe.bd.com

Norway
Tel 47.73.59.12.00
Fax 47.73.59.12.01
norge@europe.bd.com

Pakistan
Tel 92.42.5718051
Fax 92.42.5718056
bdb_pak@bd.com

Peru/Bolivia/Ecuador
Tel 51.1.712.0160
Fax 51.1.712.0165

Philippines
Tel 632.815.8981
Fax 632.815.6644

Poland
Tel 48.22.651.75.88
Fax 48.22.651.75.89
bdb.ema@europe.bd.com

Portugal
Enzifarma
Tel 351.21.421.93.30
Fax 351.21.421.93.39

South Africa
Tel 27.11.603.2620
Fax 27.11.804.0544 / 46
bdb.ema@europe.bd.com

Spain
Tel 34.902.27.17.27
Fax 34.91.848.81.04

Sweden
Tel 46.8.775.51.10
Fax 46.8.775.51.11
bdbnordic@europe.bd.com

Switzerland
Tel 41.61.485.22.22
Fax 41.61.485.22.00
customerservice.bdb.ch@europe.bd.com

Taiwan
Tel 8862.2722.5660
Fax 8862.2725.1768

Thailand
Tel 662.646.1800
Fax 662.646.1801

Turkey
Tel 90.212.328.2720
Fax 90.212.328.2730
bdb.ema@europe.bd.com

United Kingdom & Ireland
Tel 44.1865.78.16.88
Fax 44.1865.78.16.27
BDUK_CustomerService@europe.bd.com

Venezuela
Tel 58.212.241.3412 x248
Fax 58.212.241.7389

West Africa
Sobidis
Tel 225.20.33.40.32
Fax 225.20.33.40.28

