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1. Steinman RM. The dendritic cell system and its role in immunogenicity. *Annu Rev Immunol.* 1991;9:271-296.
2. Hart D, Fabre J. Demonstration and characterization of Ia-positive dendritic cells in the interstitial connective tissues of rat heart and other tissues, but not brain. *J Exp Med.* 1981;153:347-361.
3. Schuler G, Steinman R. Murine epidermal Langerhans cells mature into potent immunostimulatory dendritic cells in vitro. *J Exp Med.* 1985;161:526-546.
4. Austyn JM, Hankins DF, Larsen CP, Morris PJ, Rao AS, Roake JA. Isolation and characterization of dendritic cells from mouse heart and kidney. *J Immunol.* 1994;152:2401-2410.
5. Fossum S. The life history of dendritic leukocytes (DL). *Curr Top Pathol.* 1989;79:101-124.
6. Austyn JM. New insights into the mobilization and phagocytic activity of dendritic cells [comment]. *J Exp Med.* 1996;183:1287-1292.
7. Peters J, Gieseler R, Thiele B, Steinbach F. Dendritic cells: from ontogenetic orphans to myelomonocytic descendants. *Immunol Today.* 1996;17:27-278.
8. Freudenthal PS, Steinman RM. The distinct surface of human blood dendritic cells, as observed after an improved isolation method. *Proc Natl Acad Sci USA.* 1990;87:7698-7702.
9. Cameron PU, Lowe MG, Sotzik F, Coughlan AF, Crowe SM, Shortman K. The interaction of macrophage and non-macrophage tropic isolates of HIV-1 with thymic and tonsillar dendritic cells in vitro. *J Exp Med.* 1996;183:1851-1856.
10. Caux C, Dezutter DC, Schmitt D, Banchereau J. GM-CSF and TNF-alpha cooperate in the generation of dendritic Langerhans cells. *Nature.* 1992;360:258-261.
11. Romani N, Gruner S, Brang D, et al. Proliferating dendritic cell progenitors in human blood. *J Exp Med.* 1994;180:83-93.
12. Zhou LJ, Tedder TF. Human blood dendritic cells selectively express CD83, a member of the immunoglobulin superfamily. *J Immunol.* 1995;154:3821-3835.
13. Zhou LJ, Tedder TF. CD14<sup>+</sup> blood monocytes can differentiate into functionally mature CD83<sup>+</sup> dendritic cells. *Proc Natl Acad Sci USA.* 1996;93:2588-2592.
14. Hock BD, Starling GC, Daniel PB, Hart DN. Characterization of CMRF-44, a novel monoclonal antibody to an activation antigen expressed by the allostimulatory cells within peripheral blood, including dendritic cells. *Immunology.* 1994;83:573-581.
15. Fearnley DB, McLellan AD, Mannering SI, Hock BD, Hart DNJ. Isolation of human blood dendritic cells using the CMRF-44 monoclonal antibody: implications for studies on antigen-presenting cell function and immunotherapy. *Blood.* 1997;89:3708-3716.
16. Olweus J, BitMansour A, Warnke R, et al. Dendritic cell ontogeny: a human dendritic cell lineage of myeloid origin. *Proc Natl Acad Sci USA.* 1997;94:12551-12556.
17. Hart DNJ, McKenzie JL. Isolation and characterization of human tonsil dendritic cells. *J Exp Med.* 1988;168:157-170.
18. Grouard G, Risoan M, Filgueira L, Durand I, Banchereau J, Liu J. The enigmatic plasmacytoid T cells develop into dendritic cells with interleukin-3 (IL-3) and CD40 ligand. *J Exp Med.* 1997;185:1101-1111.
19. O'Doherty U, Peng M, Gezelter S, et al. Human blood contains two subsets of dendritic cells, one immunologically mature and the other immature. *Immunology.* 1994;82:487-493.
20. Thomas R, Lipsky PE. Human peripheral blood dendritic cell subsets. Isolation and characterization of precursor and mature antigen-presenting cells. *J Immunol.* 1994;153:4016-4027.
21. Reid CD, Stackpoole A, Meager A, Tikerpae J. Interactions of tumor necrosis factor with granulocyte-macrophage colony-stimulating factor and other cytokines in the regulation of dendritic cell growth in vitro from early bipotent CD34<sup>+</sup> progenitors in human bone marrow. *J Immunol.* 1992;149:2681-2688.
22. Young JW, Szabolcs P, Moore MA. Identification of dendritic cell colony-forming units among normal human CD34<sup>+</sup> bone marrow progenitors that are expanded by c-kit-ligand and yield pure dendritic cell colonies in the presence of granulocyte/macrophage colony-stimulating factor and tumor necrosis factor alpha. *J Exp Med.* 1995;182:1111-1120.
23. Grouard G, Durand I, Filgueira L, Banchereau J, Liu Y. Dendritic cells capable of stimulating T cells in germinal centres. *Nature.* 1996;384:364-367.
24. Hsu FJ, Benike C, Fagnoni F, et al. Vaccination of patients with B-cell lymphoma using autologous antigen-pulsed dendritic cells. *Nat Med.* 1996;2:52-58.
25. Grabbe S, Beissert S, Schwarz T, Granstein RD. Dendritic cells as initiators of tumor immune responses: a possible strategy for tumor immunotherapy? *Immunol Today.* 1995;16:117-121.
26. Girolomoni G, Ricciardi-Castagnoli P. Dendritic cells hold promise for immunotherapy. *Immunol Today.* 1997;18:102-104.
27. Agis H, Fureder W, Bankl HC, et al. Comparative immunophenotypic analysis of human mast cells, blood basophils and monocytes. *Immunology.* 1996;87:535-543.

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28. Agis H, Beil WJ, Bankl HC, et al. Mast cell-lineage versus basophil lineage involvement in myeloproliferative and myelodysplastic syndromes: diagnostic role of cell-immunophenotyping. *Leukemia and Lymphoma*. 1996;22:187-204.
29. Valent P, Besemer J, Muhm M, Majdic O, Lechner K, Bettelheim P. Interleukin 3 activates human blood basophils via high-affinity binding sites. *Proc Natl Acad Sci USA*. 1989;86:5542-5546.
30. Valent P. Immunophenotypic characterization of human basophils and mast cells. *Chem Immunol*. 1995;61:34-48.
31. Stain C, Stockinger H, Scharf M, et al. Human blood basophils display a unique phenotype including activation linked membrane structures. *Blood*. 1987;70:1872-1879.

## Contributors

Kerstin Willmann and Johanna Olweus  
Becton Dickinson Immunocytometry Systems, San Jose, CA USA

FACSComp, AutoCOMP, PAINT-A-GATE<sup>PRO</sup>, Attractors, CellQuest はBecton, Dickinson and Companyの商標、FACS, Vacutainer, FALCONは登録商標です。

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**BD Biosciences**  
Clontech  
Discovery Labware  
Immunocytometry Systems  
Pharmingen