BD® CBA Analysis Software 簡易マニュアル V1.0.2

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Step 1. Choose Samples and Standards

You may drag and move files between sections

Home

Next Step

Samples

0 Files, Sorted by Timestamp

Select or Drop
FCS Sample Files Here

Max file size: 100 MB, File type: FCS

Standard

0 Files, Sorted by Timestamp

Select or Drop
FCS Standard Files Here

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Step 1. Choose Samples and Standards

You may drag and move files between sections

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Create New

Search...

Name	Model	Type	Description	Analytes	Detail	Edit	Remove
Human Th1/Th2 Cytokine Kit II	Human	1D	BD® CBA				
Test2	Human	1D					

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①-1 2D: Human Soluble Protein - 54 beads

NameModel①-2 TypeDescription

Human Soluble ProteinHuman2D BD CBA Flex™ Sets

Debris Parameters

②X FSC-AY SSC-A

Cluster Parameters

③X SelectY Select

Reporter Parameter

④ Select

⑤Analytes54 Total Beads

Human IgE

Human MIG

Human IFN-γ

Human IL-13

Human IL-12p70

Human IL-12/IL-23p40

Human MIP-1β

⑥Save OnlySave and Select

Select from BD Catalog (optional)...

Select from BD Catalog (optional)...

1D Kits and Flex Sets

1D: Human Anaphylatoxin Kit - 3 beads

1D: Human Chemokine Kit - 5 beads

1D: Human Th1/Th2 Cytokine Kit - 6 beads

1D: Human Th1/Th2 Cytokine Kit II - 6 beads

1D: Human Inflammatory Cytokine Kit - 6 beads

1D: Human Th1/Th2/Th17 Cytokine Kit - 7 beads

1D: Mouse Inflammation Kit - 5 beads

1D: Mouse Inflammation Kit II - 6 beads

1D: Mouse Th1/Th2/Th17 Cytokine Kit - 7 beads

1D: Mouse Immunoglobulin Isotyping Kit - 7 beads

1D: Non-Human Primate Th1/Th2 Cytokine Kit - 6 beads

2D Kits and Flex Sets

2D: Human Immunoglobulin - 6 beads

2D: Human Enhanced Sensitivity (ES) - 11 beads

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Step 4. Gating Your Beads
Assign or verify gated beads for each file

①

Tube_001.fcs
1 of 15 Files

②

Apply Gates To All Files

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ID	Name
A1	Human IFN-γ
A2	Human TNF
A3	Human IL-10
A4	Human IL-6
A5	Human IL-4
A6	Human IL-2

Flex setの場合

Step 4. Gating Your Beads
Assign or verify gated beads for each file

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● Tube_001.fcs
Collision detected between A6 and A5
Collision detected between A5 and A6

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© 2021 BD and EyeCue Lab. All rights reserved. Unless otherwise noted, BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.

8, Standardsカーブの確認(Step5)

Step 5. Standards

Set your regression type for each file

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9

© 2021 BD and EyeCue Lab. All rights reserved. Unless otherwise noted, BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.

9, Report/結果(Step6)

Step 6. Your Report

Review and export your experimental data

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
FlowGrid Clustering Algorithm

For debris gating and bead population clustering:

- Ye, X., Ho, J. Ultrafast clustering of single-cell flow cytometry data using FlowGrid. *BMC SystBiol* 13, 35 (2019). <https://doi.org/10.1186/s12918-019-0690-2>

Levenberg-Marquardt Curve Fitting Algorithm

For estimating Non-Linear Regressions

- <https://github.com/mljs/levenberg-marquardt> - which is based on:
 - Brown, K.M., Dennis, J.E. Derivative free analogues of the Levenberg-Marquardt and Gauss algorithms for nonlinear least squares approximation. *Numer. Math.* 18, 289–297 (1971). <https://doi.org/10.1007/BF01404679>
 - Gavin, H. (2013). *The Levenberg-Marquardt method for nonlinear least squares curve-fitting problems* c ©. <https://people.duke.edu/~hpgavin/ce281/lm.pdf>

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