



Reaction Phenotyping Methods using Recombinant Enzymes and HLM

Christopher Patten, PhD

BD Biosciences

July 24, 2012

Today's Presentation

- Introduction
 - Drug Metabolism overview
 - *In Vitro* Model Systems
- Reaction Phenotyping

Goal of *In Vitro* ADME Testing: Identify Ideal Drug Candidates early in Development Process

Ideal Drug from ADME/Drug Development View

- **Administration is oral and the drug is easily absorbed**
 - Good permeability and aqueous solubility (class I compounds)
- **Metabolically stable (but not too stable)**
 - Swallow the pill once a day
- **Predictable Metabolism**
 - Linear drug metabolism kinetics
- **Balanced Clearance**
 - Renal or biliary secretion of parent drug
 - Metabolism to limited number of inactive products
 - Metabolism by several P450s (>2)
 - Metabolism should not depend on Polymorphic P450s
- **Not an Inhibitor or Inducer of ADME Enzymes**
 - P450, UGT, and MDR1 (P-gp)
- **Small First-Pass Effect (liver or gut)**
- **Wide Therapeutic Index**



Phase I and Phase II Biotransformation

Metabolism/Biotransformation is divided into two groups

Phase I: Addition or unmasking of functional, polar moiety (CYP, FMO, MAO, AO)

- Oxidation and/or Reduction
- Hydrolysis
- Most typical is Hydroxyl group created or exposed


Phase II: Conjugation with small, endogenous substance, often takes advantage of functional group added in Phase I

End Result: increase polarity and aqueous solubility of drug which facilitates elimination from the body



Phase I and Phase II Metabolism

- Can be active or inactive at target site
- Toxic
- Mutagen or Carcinogen



- Expose functional group that can be conjugate
- Small increase in hydrophilicity

- Large increase in hydrophilicity
- Conjugates are generally inactive

Sites of Drug Metabolism

Primary site of Metabolism is Liver


- Extra-Hepatic metabolism can be important
 - “Portals of entry” tissues often have significant drug metabolism capability
 - GI (small intestine) is major site of metabolism next to liver
 - Other sites include: Lung, skin, nasal mucosa, kidney
 - CYPs profile in extra-hepatic tissues is often different than liver

First Pass Effect

- Drug absorbed in small intestine and transported to liver via portal vein
- Extensive metabolism in Liver and/or intestine
- Limited systemic availability



Bioavailability for Oral Medication (F)




P450s Important for Drug Metabolism

- >60% of drugs on market are metabolized by P450s
- P450 pathway is often rate limiting pathway
- Families 1-3 are the Drug Metabolism CYPs
- “Big 7” CYPs: 1A2, 2B6, 2C8, 2C9, 2C19, 2D6, 3A4
- Second tier CYPs
 - 2A6, 2E1, 4F(2/3/12), 2J2, 4A11



UDP-Glucuronosyltransferases (UGTs)

- Most important Phase II drug metabolizing enzyme
- UGTs conjugate glucuronic acid to lipophilic substrates to more water-soluble metabolites, glucuronides, to facilitate excretion
- Glucuronidation rxns:
 - O-glucuronidation
 - N-glucuronidation
 - Acyl-glucuronidation (NSAID)
- UGT Isoforms
 - Hepatic: UGT1A1, UGT1A3, UGT1A4, UGT1A6, UGT1A9, UGT2B7, UGT2B15 and UGT2B10
 - GI specific: UGT1A7, UGT1A8 and UGT1A10



Model Systems for Predicting Metabolic Pathways

In Vitro Systems

I. Hepatocytes: Prepared from fresh human livers (organ donors)

- Gold-Standard for DM Studies
 - Contain all the enzymes/transporters and co-factors for drug metabolism
 - Metabolic stability (Screening for long half-life drugs)
 - Metabolite profiling (structures of metabolites)
 - Liver toxicity studies
 - Enzyme induction studies (P450 induction)
 - *In vitro* / *In vivo* scaling

Not used for:

- Enzyme Mapping/Reaction Phenotyping
- DDI

II. Liver Slices

- Similar to hepatocytes in that they contain the full complement of hepatic DMEs
 - Harder to prepare than other systems
 - Seldom used for ADME studies



In Vitro Systems

III. Liver Microsomes

- Contain all P450s, FMOs, and UGTs
- Easy to prepare and can be stored for long periods (-80°C)
- Withstand several Freeze/Thaw cycles
- Can make Donor Pool; 20 to 150 donors (average patient in population)
- **BD UltraPool™ HLM 150**
 - 150 donors
 - **Average patient**
 - **Very low lot to lot variability**
 - Equal gender ratio
 - Adult donors only (no pediatrics)
 - Pool contains equal amounts of microsomal protein from each donor
 - Pool quality livers not used (based on P450 spectra; low P420 content)
- Uses: Drug half-life (in vitro scaling), DDI, metabolite profiles, enzyme mapping, mechanistic studies

IV. Liver S9

- Same as microsomes, but contains cytosolic enzymes (SULT, GST, AO/XO, ADHs, NATs)
- Same advantages and uses as HLM
- P450 activity ~five-fold lower vs HLM



In Vitro Systems

V. Recombinant Enzymes

- Single DME expressed in a cell line
 - Baculovirus-insect cells (BD Supersomes™)
 - Yeast
 - E. coli
 - Mammalian cells
- Co-expressed with P450 OR and in some cases with b5
- Reaction phenotyping (EM), DDI, mechanistic studies
- Can be useful for studying highly stable compounds (more active than HLM)
- Requires RAF or ISEF for comparing to HLM activity



Reaction Phenotyping Methods



Reaction Phenotyping

- Pharmaceutical companies are required to characterize all the routes of clearance (drug elimination)
- Why? Because impairment of a route of elimination pathway (by a DDI, polymorphism or disease) can elevate drug levels and cause toxicity
- Generally, pharmaceutical companies are looking for drugs with multiple routes of elimination
 - If any one route is impaired, the others can compensate
- Reaction phenotyping measures the proportion of metabolism elimination which is carried out by the different enzymes
- Requires a sophisticated experimental approach
- Three basic approaches; approach used will depend on development stage





Reduce chance of Drug Candidate becoming a Victim – Increase the Number of Routes of Clearance



One (bad)



Two (better)



Many (ideal)



Integrated Approach to Reaction Phenotyping

Which Enzyme is Most Important (Principle P450)

cDNA Expressed Systems
(recombinant enzyme
Panel)

Specific Chemical or
Antibody Inhibitors with
Pooled HLM

Correlation Analysis with
Single Donor HLM Panel
(12 single donors)

- Incubations with panel of recombinant CYPs to identify important P450s. Use Relative Activity Factors (RAFTs) approach and/or relative hepatic abundance of the enzymes to determine relative importance when multiple CYPs involved. Scaled to predict *in vivo* clearance.
- The effect of co-incubated CYP-selective chemical or monoclonal antibody inhibitors on rates of metabolism in HLM can be used to identify primary DMEs
- A correlation of rate of metabolism can be made with a panel of HLM donors ($n \geq 10$) that have been phenotyped for the major DMEs
- Each can have own limitations – combined approach is typically employed



Reaction Phenotyping in Drug Development Stages

- Reaction phenotyping conducted at multiple stages during drug development process
- Approaches will vary depending on stage (and specific laboratory preferences)
- **Discovery Stage**
 - Loss of parent assay method at a low drug concentration ($\sim 1 \mu\text{M}$); radiolabelled compounds not available for metabolite quantitation
 - Typically involves high-throughput systems (multi-well formats)
 - Generally use single method Rx Phenotyping approach; limit testing to major CYPs
 - Pooled HLM with selected chemical/antibody inhibitors, or
 - cDNA expressed enzyme panel (in some labs cDNA panel is first tier approach)
- **Early Development Stage**
 - Expand studies with HLM and cDNA-expressed CYPs; use comprehensive panels of inhibitors and/or cDNA-expressed enzymes. Run tests at low drug concentration (single concentration).
 - Measure metabolite formation; identify major metabolites; use radiolabelled test compounds
 - Identify major CYPs for each major metabolite
 - Determine if any Polymorphic CYPs are involved
- **Full Development Stage**
 - Determine full kinetics in HLM (K_m and V_{max}); metabolite formation with radiolabelled compound
 - Consider non-specific binding in microsomes (fu)
 - Detailed Rx Phenotyping studies, focus on major CYPs (>25% of clearance pathway), include all 3 methods
 - Inhibition in HLM
 - cDNA-expressed enzyme panel
 - Correlation analysis with single donor HLM panel



Common First Step: Determine CYP/non-CYP Involvement

- Flavin Containing Monooxygenase (FMO): microsomal enzyme, dependant on NADPH and O₂ for activity (same as CYP)
 - Both CYP and FMO catalyze hetero-atom oxidation (N- and S-oxidation)
 - C-oxidations carried out by CYP, not FMO
- Aldehyde Oxidase: cytosolic and do not require co-factor
 - Can carry out similar reactions as CYPs
 - Measure activity in HLM vs cytosol and/or S9 (activity in absence of co-factor, and cytosolic > LM)
- Monoamine Oxidase (MAO-A/B): both forms abundant in liver mitochondria
 - No co-factor requirements, both inhibited by Pargyline

Methods to rule out non-CYP oxidative pathways

- 1-Aminobenzotriazole (ABT) is a general CYP inhibitor that can be used to distinguish CYP from non-CYP pathways
 - ~1 mM ABT pre-incubated with HLM, with and without NADPH
 - CYP activity is decreased by 1-ABT in the presence of NADPH
 - Inhibition is not equal among CYPs (most potent for CYP3A4, least for CYP2C9)
 - 1-ABT can be used to inactivate CYPs in hepatocytes
- FMO vs CYP
 - FMO are heat labile in absence of NADPH: pre-heat HLM to 45°C for 1 minute with/without NADPH
 - In absence of NADPH heat will inactivate FMO, but not CYP
 - Non-Ionic detergent inactivation of CYPs: CYPs are sensitive to non-ionic detergents (1% Triton X-100), while FMOs are resistant
 - FMOs not inhibited by general CYP inhibitors: e.g. 1-ABT (methimazole inhibits FMO, but also several CYPs)



Establishing Linear Assay Conditions

- Determine linear assay conditions for time and protein concentration in HLM
 - Establish using pooled HLM
- Test linearity using multiple time points (e.g. 5 to 60 min) and protein concentrations (e.g. 0.2, 0.4 and 0.8 mg/mL) at multiple substrate concentrations (e.g. 1, 10 and 100 μ M)
- Hold one parameter constant (e.g. protein) and vary the other (time), and test at each substrate concentration
- Repeat the process holding the other parameter constant
- The linear assay conditions established for pooled HLM will be used for all future studies involving HLM: Kinetic analysis and phenotyping with HLM (inhibition and correlation analysis)



Substrate Concentration Considerations

- In Discovery Stage kinetic parameters are not known; typically start with low substrate concentration and measure loss of parent
- Kinetic parameters (K_m and V_{max}) for HLM (pooled) are typically determined in development. Measure metabolite(s) formation with radiolabeled test compounds.
- Assay needs to be carried out under linear conditions for both time and protein concentration
 - Substrate utilization <10% if possible (can be difficult at low $[S]$), not to exceed 20%.
 - May have to re-adjust time or protein to achieve low substrate utilization
 - Adjust for non-specific binding (determine free-fraction); extensive binding (lipophilic compounds) will decrease effective concentration and increase K_m
- Typically use 12 to 16 substrate concentrations: $0.1 \times K_m$ to $10 \times K_m$
 - Initial K_m can be estimated from linearity studies
- Plot data to obtain kinetic parameters: K_m & V_{max}
 - Determine if single or multiple K_m values
 - Simple Michaelis-Menten or Allosteric Kinetics (S_{50})
 - Intrinsic Clearance (Cl_{int}) = V_{max}/K_m (mL/min/mg HLM)



P450 Enzyme Kinetic Plots

Michaelis-Menten Kinetics


$$v = \frac{V_{\max} * S}{K_m + S}$$

Linear Transformation

Rectangular Hyperbola



Direct Plot: V vs [S]

$V = V_{\max}$




K_m : substrate concentration that gives $\frac{1}{2} V_{\max}$

Eadie-Hofstee Plot: $v = V_{\max} - K_m (v/[S])$
Exaggerates deviations from the Michaelis-Menten kinetics



Rx Phen at two [S]; high and low



BD

Reaction Phenotyping with Chemical and Antibody Inhibitors (“knock-out” method)

- Enzyme source: pooled HLM (represent average patient), e.g. BD UltraPool™ HLM 150
- Substrate concentration: pharmacologically relevant [S]
 - May require 2 concentrations (high and low) if multiple CYPs involved (high/low K_m isoforms)
 - If *in vivo* concentration is not known, then $[S] < K_m$
 - Common [S] in absence of kinetic or *in vivo* data is 1 μM (typically below K_m , first-order reaction kinetics)
- CYPs can be inhibited by common solvents used to dissolve test compounds and inhibitors
 - Keep solvent concentrations to a minimum: Methanol (<1%), ACN (<2%) and DMSO (<0.2%)
 - Include solvent negative controls for chemical inhibitors
- Negative controls for antibodies: Pre-immune serum or ascites, or irrelevant antibody
 - Antibodies are seldom purified and contain albumin; may result in non-specific binding of test compound or other artifacts



Chemical Inhibitors

- Chemical inhibitors are typically inexpensive and easy to use
- Selective, potent inhibitors not available for all CYP isoforms (e.g. inhibitors to distinguish CYP3A4 & 3A5 involvement)
- Can have a narrow range where potency and specificity are at a maximum
 - E.g. Ketoconazole, CYP3A4 inhibitor, can inhibit other CYPs at low micromolar concentrations: CYP4F2/12/3B, 2C8/9, 2B6, 2J2, 1A1 and 1B1
 - Azamulin shown to be more selective vs KTZ (Stresser et.al.)
- Mechanism of inhibition needs to be considered (competitive, non competitive, time dependent [increased potency over time])
 - Competitive inhibitor potency depends on [S]
 - TDI requires pre-incubation with NADPH
- Chemical inhibitors can be subject to metabolism/depletion in HLM; reduced potency over time
- Inhibitor potency may be substrate-dependent (IC_{50} varies across substrates)
- Inhibitor of one CYP may activate another; e.g. α -naphthoflavone potent inhibitor of CYP1A2, but can activate CYP3A4



CYP Inhibitors (Zhang, Expert Opinion, Review, 2007)

| CYP form(s) | Inhibitor | Inhibitor concentrations |
|---------------------|-----------------------------------------|-------------------------------|
| CYP1A2 | Furafylline α -Naphthoflavone | 10 – 30 μ M* 1 μ M |
| CYP2A6 | Methoxsalen | 1 μ M |
| CYP2B6 | ThioTEPA | 50 μ M |
| CYP2C8 | Montelukast | 0.1 μ M† |
| CYP2C9 | Sulfaphenazole | 10 μ M |
| CYP2C19 | Benzylrivanol | 1 μ M |
| CYP2D6 | Quinidine | < 2 μ M |
| CYP2D6 and CYP3A4/5 | | 10 μ M |
| CYP2E1 | Diethyldithiocarbamate | 50 μ M |
| CYP3A4/5 and others | Ketoconazole | 1 and 10 μ M |
| CYP3A4/5 | Troleandomycin | 50 μ M* |

*Inhibition potency increases with preincubation in the presence of NADPH (time-dependent inhibitors). A 15- to 30-min preincubation is recommended.

†Montelukast is a potent CYP2C8 inhibitor. However, protein binding is a significant factor and concentrations of the inhibitor should be adjusted for liver microsomal protein content.



Inhibitory Antibodies

- Commonly used in Reaction Phenotyping Studies
- Complete panel to cover all isoforms not available
- Can be more specific than chemical inhibitors (i.e. monoclonals and peptide antibodies)
- Often more costly vs chemical inhibitor
- Mechanism of inhibition is non-competitive; lowers V_{max} (Qin Mei et al, JPET, 1999 & DMD, 2002)
 - IH-Abs inhibit independent of test compound concentration (chemical inhibitor potency/selectivity can be effected by substrate and substrate-concentration)
- Inhibitory UGT, FMO antibodies are not available (Western blot antibodies only)



Inhibitory Antibodies (continued)

- Majority of commercially available inhibitory CYP antibodies are monoclonals (NIH); single epitope
- Not all IH-MABs are specific (cross-reactivity)
- Incomplete inhibition is somewhat common
 - MAB-CYP/Substrate complex remains productive
 - Can make data interpretation difficult in cases where multiple CYPs are involved
- Peptide antibodies can be made highly specific (pre-determined epitope), however they are typically non-inhibitory...and when they are not all rabbits will produce inhibitory antibodies (need to inject multiple rabbits and screen serum)




Antibody Inhibition Assay Basics

- Antibody should be specific and potent
 - >90% inhibition in recombinant and HLM, otherwise results will be difficult to interpret if multiple CYPs involved
 - Typically the potency is greater in recombinant systems vs HLM
- Should run a titration curve to insure maximal inhibition has been achieved (partial inhibition of a CYP makes interpretation difficult)
- In Theory, if one CYP is involved, ~90% inhibition should be achieved...if only 50% is achieved, then one or more additional CYPs are involved
 - Best case scenario: % Inhibition should add up to ~100%
- A good practice is to conduct initial study in pooled HLM, and follow up with single donors high/low in the CYP of interest (confirmation studies)




CYP2C19 and 2C8 Inhibitory Antibody

Specificity of Inhibition by IH-MAB-2C19




- Both show high specificity
- >90% inhibition of activity in HLM

Specificity of Inhibition by IH-MAB-2C8



CYP1A2 – Inhibitory MAB

Specificity of Inhibition by MAB-1A2





- High specificity for 1A2
- Incomplete inhibition of 1A2 in HLM (~80%)

Antibody vs Chemical Inhibition

(Qin Mei et al., JPET, 1999, vol. 291)

- CYP3A4 contributes only 52 to 73% to Diazepam N-demethylation
- Mab-3A4a (A) detects the involvement of other CYPs
- Ketoconazole selectivity (B) is shown to be concentration dependent (over predicts the contribution of CYP3A4)



BD Gentest™ Chemical/Antibody Reagents for CYP Reaction Phenotyping Applications

| CYP | Probe Substrate | Specific Inhibitor |
|---------|-------------------------|------------------------------------------------------|
| CYP1A2 | Phenacetin | Furafylline ^{BD} |
| CYP2B6 | Bupropion ^{BD} | CYP2B6 Mab ^{BD} |
| CYP2A6 | Coumarin | CYP2A6 Mab ^{BD} |
| CYP2C8 | Amodiaquine | Montelukast |
| CYP2C9 | Diclofenac | Sulfaphenazole ^{BD} |
| CYP2C9 | Tolbutamide | Sulfaphenazole ^{BD} |
| CYP2C19 | (S)-Mephenytoin | Benzylrinvanol ^{BD} |
| CYP2D6 | Dextromethorphan | Quinidine |
| CYP2D6 | Bufuralol | Quinidine |
| CYP3A4 | Nifedipine | Ketoconazole ^{BD} or Azamulin ^{BD} |
| CYP3A4 | Testosterone | Ketoconazole ^{BD} or Azamulin ^{BD} |
| CYP3A4 | Midazolam ^{BD} | Ketoconazole ^{BD} or Azamulin ^{BD} |



Correlation Approach

- Takes advantage of wide inter-individual variability of CYP expression levels between donors (>100-fold variability in CYP activity between donors)
- Most direct method – not prone to experimental artifacts from other Rx Phenotyping systems (e.g. chemical/antibody specificity)
- Requires panel of characterized single donor HLM samples; at least 10 single donors
- Panel should have low internal correlation between any two CYPs
- Correlate rates of metabolism of NCE vs CYP-probe substrate activity in panel and/or CYP abundance in panel - unknown & reference activity/abundance
- Labor intensive, typically carried out in development stage. Used to confirm results from HLM inhibition & recombinant CYP studies
- Works best when one enzyme is involved per metabolite
 - Conduct analysis with CYP inhibitors when multiple CYPs suspected
- Need to examine the graph
 - Correlation dependent on single point outlier
 - Check if regression line passes through zero
 - Positive y-axis intercept and high correlation could mean another enzyme is involved




Low Correlation Among Enzyme Activities

| | CYP1A2 | CYP2C9 | CYP2C19 | CYP2B6 | CYP2D6 | CYP3A4 |
|---------|---------|---------|---------|---------|---------|---------------|
| CYP1A2 | | 0.2446 | 0.0685 | 0.0025 | -0.3785 | 0.2627 |
| CYP2C9 | 0.2446 | | -0.3924 | 0.3876 | 0.1798 | 0.4786 |
| CYP2C19 | 0.0685 | -0.3924 | | -0.1748 | -0.4706 | 0.0769 |
| CYP2B6 | 0.0025 | 0.3876 | -0.1748 | | 0.1844 | 0.4683 |
| CYP2D6 | -0.3785 | 0.1798 | -0.4706 | 0.1844 | | -0.3600 |
| CYP3A4 | 0.2627 | 0.4786 | 0.0769 | 0.4683 | -0.3600 | |


| UGT-Isoforms | UGT1A1 | UGT1A4 | UGT1A9 |
|--------------|--------|--------|--------|
| UGT1A1 | | -0.252 | 0.317 |
| UGT1A4 | | | 0.135 |
| UGT1A9 | | | |

Correlation Analysis Examples

CYP2C9 and UGT1A1



CYP2C9 Substrates



UGT1A1 Substrates



Reaction Phenotyping using Recombinant CYPs

- cDNA-Expressed P450s: Single CYP isoform expressed in a cell-based system (BD Supersomes™)
 - Co-expressed with P450 Oxidoreductase (OR), with or without cytochrome b5
 - b5 can stimulate, inhibit or have no effect depending on CYP and substrate
- Incubate test compound with a panel of cDNA-expressed enzymes
- Major advantage to cDNA approach is that all the major CYPs and UGT cDNA-expressed enzymes are commercially available (BD Supersomes™), unlike chemical or antibody inhibitors
 - Less common CYPs - CYP4F (2, 3A, 3B, 12), 2J2, 1A1, 1B1, 2C18 - are available as BD Supersomes
 - Easily determine all the important CYPs, regardless of number of isoforms involved
- Accessory proteins (OR, b5), and CYPs, are generally “over-expressed” relative to HLM, producing higher activity in recombinant system vs the same CYP in native HLM environment
 - High activity is an advantage when testing slowly metabolized compounds; HLM may not be activity enough to generate metabolites or show loss of parent
- Activity needs to be normalized to correct for the higher activity vs HLM



BD Supersome™ Characterization

- Gene sequence is perfect match with sequence published at U.S. National Library of Medicine
- Sold as insect cell membrane preparations (CYPs are membrane bound, same as in HLM; retain native membrane structure)
- Co-expressed with P450 OR and in some cases with Cytochrome b5
 - Cytochrome b5 can increase, decrease, or have no effect on activity depending on CYP and test compound
 - In most cases b5 increases Vmax (CYP3A4), but can also effect Km (e.g. CYP2E1; b5 decreases Km)
- Substrate specificity consistent with literature reports
- Michaelis-Menten constant (Km) is generally consistent with value observed in pooled HLM
 - Insures the membrane structure of the Supersome is the same as in native tissue (liver or intestine)
 - CYPs, UGTs and MAOs are all membrane bound proteins...NAT is our only soluble Supersome
 - Km: substrate concentration that gives ½ maximal activity (it is a “constant” for a given enzyme)

| Enzyme | Substrate | Km Value (uM) | |
|---------|-----------------|---------------|-----------------|
| | | Supersomes | Tissue Fraction |
| CYP3A4 | Terfenadine | 1.9 | 2.0 |
| CYP2C9 | Diclofenac | 4.2 | 2.9 |
| UGT1A1 | Estradiol | 30 | 33 |
| UGT1A4 | Trifluoperazine | 100 | 90 |
| UGT-1A9 | Propofol | 21 | 19 |
| MAO-A | Serotonin | 86 | 91 |


*Provides assurance that it
is the right protein
behaving the right way*



Involvement of Non-Traditional CYPs (CYP4F Metabolism of FTY-720, Yi *et.al.* 2011)

Correlation analysis of fingolimod metabolism rates with marker enzyme activities in a panel of HLM from individual donors

| Enzyme | Marker Reaction | Correlation (R value) |
|----------|--------------------------------------------------|-----------------------|
| CYP1A2 | Phenacetin <i>O</i> -deethylation | 0.08 |
| CYP2A6 | Coumarin 7-hydroxylation | 0.2 |
| CYP2B6 | (<i>S</i>)-Mephenytoin <i>N</i> -demethylation | 0.25 |
| CYP2C8 | Paclitaxel 6 α -hydroxylation | 0.18 |
| CYP2C9 | Diclofenac 4'-hydroxylation | 0.4 |
| CYP2C19 | (<i>S</i>)-Mephenytoin 4'-hydroxylation | 0.03 |
| CYP2D6 | Bufuralol 1'-hydroxylation | 0.39 |
| CYP2E1 | Chlorzoxazone 6-hydroxylation | 0.32 |
| CYP3A4/5 | Testosterone 6 β -hydroxylation | 0.32 |
| CYP4A | Lauric acid 12-hydroxylation | 0.32 |
| FMO | Methyl <i>p</i> -tolyl sulfide oxidation | 0.4 |




- Correlation study was negative for major CYPs
- Available chemical and antibody inhibitors showed no effect
- cDNA panel study: BD Supersomes™, CYP4F2 and 4F3B, were CYPs showing significant activity for the major metabolite (M12)
- CYP4F2 and HLM showed similar kinetics for FTY-720
 - $K_m \sim 100 \mu\text{M}$, showed substrate inhibition in both systems



BD Supersomes™ versus HLM

- High Expression Level of OR, b5 and CYP
- High Activity (typically several fold higher than HLM)
- Scaling Methods: Normalize rCYP data to accurately predict CYP activity in HLM environment; assign relative contribution when multiple CYPs shown to metabolize test compound
 - Relative abundance approach
 - Multiply BD Superome activity for NCE x CYP abundance in HLM
 - RAF
 - ISEF



When more than one CYP involved, which is most important?

“Relative Activity Factor” Method – RAF (Crespi, C *et.al.*)

- **RAF: Relates the activity of cDNA-expressed enzyme to the activity of the enzyme in its native environment - HLM.**
- **RAF = HLM activity of probe substrate/rCYP activity of probe substrate**
 - $(\text{pmol/min} \cdot \text{mg}) / (\text{pmol/min} \cdot \text{pmol}) = \underline{\text{pmol CYP/mg}}$
- Example (RAF): Recombinant enzyme panel shows that two P450 enzymes metabolize NCE (3A4 and 2C9)
 - rCYP3A4 Activity for NCE = 100 pmol/min/pmol
 - rCYP2C9 Activity for NCE = 20 pmol/min/pmol
 - rCYP3A4 Activity for Probe (Testosterone) = 200 pmol/min/pmol
 - rCYP2C9 Activity for Probe (Diclofenac) = 40 pmol/min/pmol
 - 3A4 Testosterone 6 β -hydroxylase in HLM (probe) = 3000 pmol/min/mg
 - 2C9 Diclofenac 4-hydroxylase in HLM (probe) = 1000 pmol/min/mg
- **RAF 3A4 = 3A4 activity for probe in HLM/r3A4 activity for probe = 15 pmol/mg**
 - Rate of 3A4 NCE activity In HLM = RAF3A4 x r3A4 NCE activity = 1500 pmol/min.mg
- **RAF 2C9 = 2C9 activity for probe in HLM/r2C9 activity for probe = 25 pmol/mg**
 - Rate of 2C9 NCE activity In HLM = RAF2C9 x r2C9 NCE activity = 500 pmol/min.mg



RAF (continued)

- **Predicated Total CYP NCE Activity in HLM:**
 - Sum of calculated NCE activity by all P450s involved:
 $1500 \text{ (3A4)} + 500 \text{ (2C9)} = 2000 \text{ pmol/min*mg}$
- **Percent Contribution of CYP3A4 and CYP2C9 to NCE Activity in HLM:**
 - CYP3A4: $1500/2000 \times 100 = 75\%$
 - CYP2C9: $500/2000 \times 100 = 25\%$
- **Confirm RAF result with P450 specific inhibitory antibodies or chemical inhibitors:**
 - CYP3A4: 1 μM Ketoconazole = ~75% Inhibition in HLM
 - CYP2C9: 2 μM Sulfaphenazole = ~25% Inhibition in HLM



RAF (continued): I-S-E-F

- ISEF is “Inter System Extrapolation Factor” (simCYP)

$$\text{ISEF} = \frac{V_{\max} \text{ or } CL \text{ (probe-HLM pool)}}{V_{\max} \text{ or } CL \text{ (probe-rCYP)} \times \text{CYP abundance in HLM pool (WB)}}$$

- ISEF has no units (cancel out): turnover# of CYP in HLM/ turnover # of rCYP
- Takes CYP abundance in HLM into account
- The ISEF is used to estimate NCE Clearance using rCYP (Supersomes) as the enzyme source. Can also estimate DDI severity in different populations
- ISEF allows for population variability to be brought into the model (ethnicity , age or disease state).
- Intrinsic Clearance of NCE for a specific CYP (mL/min*mg):

$$CL_{\text{int-NCE}} = [CL\text{-ISEF}_{\text{CYPx}} \times (V_{\max\text{-NCE}} [rCYP]/K_m\text{-NCE} [rCYP])] \times \text{CYP-abundance}_{(HLM)}$$


- *CYP-abundance* is where CYP population variability is brought into model
- ISEF should be consistent across all probe substrates for given CYP (work in progress at Simcyp)
- Scale Cl_{int} to *in vivo* and hepatic clearance



Qualification of ISEF Method

Chen Y. DMD, 2011

- Studied 20 compounds (10 marketed drugs, 10 in-house compounds)
- Determined ISEF for major CYPs (CYP1A2, 2C8, 2C9, 2C19, 2D6, 3A4)
- Used multiple probe substrates (2C9, 2D6, 3A4)
 - ISEF varies per probe substrate: but were within 2-fold of one another
 - Take average when using multiple probes
- Methods: loss of parent (CI), metabolite formation, V_{max}/K_m (CI)
 - ISEF determined from each method were similar (~within 2-fold)
- Compared ISEF method vs HLM (inhibition) for both Rx Phenotyping and Cl_{int}
 - Tight match between % contribution (ISEF vs HLM)
 - Tight match for Cl_{int} (within 2-fold); CI-ISEF showed lowest fold error



Conclusion

- FDA requires pathways of elimination identified ($\geq 25\%$ may require in vivo study)
- HLM and recombinant P450s are useful (predictive) models for reaction phenotyping
- An integrated approach involving the 3 methods discussed (or at least 2) is recommended
 - Correlation analysis is being used less; labor intense and interpretation issues when multiple CYPs involved
 - “Knock-out” in HLM and rCYP methods currently the most common
- All major (and most minor) recombinant CYPs available for Rx Phen testing; data should be normalized using either RAF or ISEF
- Antibody and chemical inhibitors are available for major CYPs, but not for less common CYP
 - Room for improvement in chemical/antibody specificity for major CYP inhibitors (more reliable data interpretation)
- Tools for mapping non-CYP metabolic enzymes and transporter pathways are still somewhat limited
- General trend is that compounds are becoming more stable and involve more non-CYP pathways for metabolism



Contact Us



Questions?

Contact Information:

Christopher Patten, PhD

e-mail: chris_patten@bd.com

Technical Support:

In the U.S.

tel: 877.232.8995 or or 978-901-7491

e-mail: labware@bd.com

Outside the U.S.

Contact your local distributor or visit
bdbiosciences.com/offices to locate
your nearest BD Biosciences office.

For research use only. Not intended for use in diagnostic or therapeutic procedures.
BD, BD Logo, and all other trademarks are the property of Becton, Dickinson and Company. ©2009 BD



Supplemental Slides

- BD Supersome™ Portfolio
- Scaling Factors
- Scaling Intrinsic Clearance
- Azamulin Specificity vs Ketoconazole

Human BD Supersomes™ Product Portfolio

| BD Supersome | Cytochrome b5 |
|--------------|-------------------|
| CYP1A1 | |
| CYP1A2 | |
| CYP1B1 | |
| CYP2A6 | (b ₅) |
| CYP2A13 | |
| CYP2B6 | (b ₅) |
| CYP2C8 | (b ₅) |
| CYP2C9 | (b ₅) |
| CYP2C9*2 | |
| CYP2C9*3 | |
| CYP2C18 | |
| CYP2C19 | (b ₅) |
| CYP2D6 | |
| CYP2E1 | (b ₅) |
| CYP2J2 | (b ₅) |
| CYP3A4 | (b ₅) |
| CYP3A5 | (b ₅) |
| CYP3A7 | (b ₅) |
| CYP4A11 | |
| CYP4F2 | (b ₅) |
| CYP4F3A | (b ₅) |
| CYP4F3B | (b ₅) |
| CYP4F12 | (b ₅) |
| CYP19 | |
| CYP51 | Custom |
| Total | 23 |

| Non-CYP BD Superomes |
|-------------------------------------------------------------------------------------------|
| CES1b (Liver) |
| CES1c (Liver) |
| CES2 (Intestine) |
| |
| UGT1A1, 1A3, 1A4, 1A6, 1A7, 1A8, 1A9, 1A10, 2B4, 2B7, 2B15, 2B17 (2B10, 2B28 - Dec. 2011) |
| |
| FMO1, FMO3, FMO5 |
| MAO-A, MAO-B |
| NAT1, NAT2 |



Hepatic Scaling Factors for Whole Liver Intrinsic Clearance

- Microsomal SF (MSF)
 - 40 mg microsomal protein/gr liver (Hakooz, 2006)
 - BD UltraPool™ HLM 150: 43 mg microsome/gr liver
- Hepatocellularity SF
 - 99 million cells/gr liver (Barter, 2007)
- Liver Weight (LW): 1400 grams
- Body Weight (BW): 70 kg
- 20 grams liver/kg BW

- $CL_{U_{int-in vivo}} = (CL_{U_{int}} \times MSF \times LW) / BW$
- Units = $\mu\text{L}/\text{min} \cdot \text{kg BW}$ (or $\text{mL}/\text{min} \cdot \text{kg BW}$)

$$CL_{U_{int-in vivo}} = CL_{U_{int}} (V_{max}/K_m) \times 40\text{mg microsome/gr liver} \times 20 \text{ gr liver/kg BW}$$



Scaling Intrinsic Clearance to *In Vivo* Hepatic Clearance

Vmax/Km or Half-Life
(hepatocyte/microsomes/rCYPs)




$CL\text{-}int_{in\ vitro}$



Scaling factors

$CL\text{-}int'_{in\ vivo}$



Models of
hepatic clearance

$CL\text{-}hep$

In Vivo Clearance

Scaling Factors

- 43 mg mic./gram liver
- 20 g liver/kg b.w.

Well-Stirred Model:

$$CL_h = \frac{f_{ub} * CL_{int}'_{in\ vivo} * Q_H}{f_{ub} * CL_{int}'_{in\ vivo} + Q_H}$$

Parallel Tube Model:

$$CL_h = Q_H - Q_H * \exp(-f_{ub} * CL_{int}'_{in\ vivo}/Q_H)$$

Dispersion Model:

$$CL_h = Q_H * \frac{1 - 4a}{(1 + a)^2 * \exp[(a-1)/2Dn] - (1-a)^2 * \exp[-(a+1)/2Dn]}$$

$$a = (1 + 4Rn * Dn)^{1/2}; Dn = 0.17 \text{ (dispersion no.);}$$

$$Rn = f_{ub} * CL_{int}'_{in\ vivo}/Q_H \text{ (efficiency no.)}$$

(f_{ub} = fraction unbound; Q_H = hepatic blood-flow)



Azamulin vs Ketoconazole as a selective CYP3A4 inhibitor (Stresser et.al.)

| | | | IC ₅₀ (μM) | IC ₅₀ (μM) | IC ₅₀ ratio |
|----------|----------------|------------------|-----------------------|-----------------------|------------------------|
| Enzyme | Enzyme source | Substrate | KTZ | AZA | AZA/KTZ |
| CYP1A1 | BD Supersomes™ | BzRes | 0.017 | > 50 | >2900 |
| CYP1A2 | BD Supersomes | CEC | > 50 | > 200 | - |
| CYP1A2 | HLM | Phenacetin | 55 | > 100 | > 1.8 |
| CYP1B1 | BD Supersomes | BzRes | 5.6 | > 50 | > 8.9 |
| CYP2A6 | BD Supersomes | Coumarin | > 50 | > 50 | - |
| CYP2B6 | BD Supersomes | EFC | 6.3 | > 50 | > 7.9 |
| CYP2B6 | HLM | Bupropion | ND | >50 | - |
| CYP2C8 | BD Supersomes | DBF | 4.4 | 41 | 9.3 |
| CYP2C8 | HLM | Paclitaxel | 8.6 | 73 | 8.5 |
| CYP2C9 | BD Supersomes | MFC | 5.1 | 57 | 11.2 |
| CYP2C9 | HLM | Diclofenac | 9 | > 100 | > 11 |
| CYP2C19 | BD Supersomes | CEC | 13 | 47 | 3.6 |
| CYP2C19 | HLM | S-mephenytoin | 57 | > 100 | > 1.8 |
| CYP2D6 | BD Supersomes | AMMC | ND | > 200 | - |
| CYP2D6 | HLM | Bufuralol | 67 | > 100 | > 1.5 |
| CYP2E1 | BD Supersomes | MFC | ND | > 50 | - |
| CYP2E1 | HLM | p-Nitrophenol | > 50 | > 50 | - |
| CYP2J2 | BD Supersomes | Terfenadine | 4.6 | 6.6 | 1.4 |
| CYP3A4 | BD Supersomes | BFC | 0.023 | 0.026 | 1.1 |
| CYP3A4 | BD Supersomes | Terfenadine | < 0.1 ^e | 0.12 | 1.6 |
| CYP3A4/5 | HLM | Testosterone | 0.06 | 0.24 | 4.0 |
| CYP3A4/5 | HLM | Midazolam | 0.04 | 0.15 | 3.7 |
| CYP3A5 | BD Supersomes | BFC | 0.11 | 0.38 | 3.5 |
| CYP3A7 | BD Supersomes | BFC | 0.16 | 0.34 | 2.1 |
| CYP4A11 | HLM | Lauric acid | > 100 | > 100 | - |
| CYP4F2 | BD Supersomes | LTB ₄ | < 1 ^f | 46 | 57 |
| CYP4F3b | BD Supersomes | LTB ₄ | 22 | >100 | > 4.5 |
| CYP4F12 | BD Supersomes | Terfenadine | 0.74 | 29 | 39 |

- Azamulin (AZA) is a Mechanism-Based Inhibitor of CYP3A4
- Shows enhanced specificity towards CYP3A4 vs KTZ
- KTZ can inhibit several CYPs at low micromolar concentrations (4F, 2C8/9, 2J2, 2B6, 1A1, 1B1)
- AZA can inhibit 2J2 at micromolar concentration; primarily extra-hepatic CYP (intestine)

