

The Length of Preincubation Times in Abbreviated Cytochrome P450 Time-dependent Inhibition Studies: One Size Fits All?

David M. Stresser, Ph.D.

BD Biosciences Discovery Labware




In the Context of Drug-Drug Interactions, there are Generally Two Types of CYP Inhibition

- **Reversible Inhibition**

- IC_{50} unchanged with incubation time
- Most drugs are these
 - Competitive, noncompetitive, “mixed”

- **Time-Dependent Inhibition (TDI)**


- Changes in IC_{50} with incubation time
- Irreversible - covalently bound
 - Mechanism based inactivation
- Quasi-Irreversible
 - Metabolite-intermediate complex [MIC]
- Reversible
 - Metabolite more inhibitory than parent



Consequences of TDI

- When TDI is the mode of inhibition, the inhibitory interaction will generally be greater over time following multiple dosing and be longer lasting after discontinuation
- Compared to reversible inhibition, prediction of inhibition response *in vivo* is more challenging

These Factors Drive Testing Behaviors



Percentage of respondents who test new chemical entities for their ability to cause time-dependent inhibition from a survey of 17 companies (2/3 large pharma) in March to April 2008 [Grimm SW et al (2009), Drug Metab. Dispos. 37:1355]



Goals of Abbreviated Testing for TDI

- Detect TDI in a robust but simplified test
 - Use information to assess need for further investigation, such as determination of K_I/k_{inact} or testing in additional systems
- Avoid false negatives
 - TDI missed in assay, found later unexpectedly, after significant resources consumed – to be avoided
- Avoid false positives
 - TDI found in assay, but proved not to be inactivator upon subsequent testing
 - Chasing issues unnecessarily. A nuisance, but tolerable

Regulatory Guidance - USFDA

“Time-dependent inhibition should be examined in standard in vitro screening protocols. A 30-minute pre-incubation of a potential inhibitor before the addition of substrate is recommended”

- *FDA DRAFT Guidance for Industry – Drug Interaction Studies (Sept. 2006)*

“TDI should be studied in standard in vitro screening protocols by pre-incubating the drug... before the addition of a substrate. Any time-dependent loss of initial product formation rate may indicate time-dependent inhibition, and definitive in vitro studies to obtain TDI parameters (i.e., k_{inact} and K_I). Details of this tiered approach were proposed by the PhRMA Drug Metabolism Technical Group (Grimm et al. 2009).”

- *FDA DRAFT Guidance for Industry – Drug Interaction Studies (Feb. 2012)*



Regulatory Guidance - EMEA

- *“...If the inhibition is enhanced by pre-incubations, time-dependent inhibition may be present. In this situation k_{inact} and K_I should be determined”*

- Draft EMEA Guideline on the Investigation of Drug Interactions (April, 2010)

Time-dependent Inhibition Testing Practices

- A typical initial test is an “IC₅₀ shift” assay
 - 47% of companies in 2009 PhRMA survey¹ use this model as initial TDI assessment
- 35% measure percent decline or inactivation rate at a single inhibitor concentration




¹Grimm et al, 2009



For Respondents Conducting the IC₅₀ Shift...

- About half of respondents perform IC₅₀ shift experiments by conducting the pre-incubation at a higher concentration of I and E and then diluting into incubation with probe substrate (like a typical k_{inact} experiment)
- ...while the other half conduct the activity incubation by adding the probe substrate with no dilution step

Dilution Method for IC₅₀ Shift



Preincubation

10X Test article, n = 7, 0.5 log spacing*

10X Microsomal protein

1X NADPH regenerating system

Probe Substrate not present

30 min Preincubation time

Secondary Incubation different vessel

1X Test article, n = 7, 0.5 log spacing

1X Microsomal protein


1X NADPH regenerating system

Probe Substrate present

5 min incubation time

*This fold dilution is a typical value

Non-dilution Method for IC₅₀ Shift



Preincubation

1X Test article, n = 7, 0.5 log spacing

1 X Microsomal protein

1X NADPH regenerating system

Probe Substrate not present

30 min Preincubation time

Secondary Incubation same vessel

1X Test article, n = 7, 0.5 log spacing

1X Microsomal protein


1X NADPH regenerating system

Probe Substrate present

5 min Incubation time

IC₅₀ Shift Defined

- Simple ratio of IC₅₀ values
 - “Plus NADPH” IC₅₀ value as denominator
 - “Minus NADPH” IC₅₀ value as numerator (or direct IC₅₀ value)
- If the ratio is > than cut-off value, this indicates TDI
- A recommended cut-off value is 1.5 to 2-fold (Grimm et al, 2009)
- Note: Comparison of percent inhibition data is also important, particularly when inhibition or solubility limitations preclude calculation of an IC₅₀ value




How Long to Preincubate?

- A 30 min preincubation period is often used
- Why?
 - There is a general assumption that this is adequate to permit detection of TDI
 - It was recommended in 2006 draft FDA guidance

Comparison of 10 and 30 Min Preincubation

- Results from our laboratory demonstrate that comparing pre-incubation times of 10 to 30 minutes distinguish rapid from slow acting time-dependent inhibitors¹.
- We found significant decreases in shifted IC_{50} values in the interval 10 to 30 minutes for some, but not all compounds, suggesting the 30 min period is arbitrary and possibly too short



¹ - Perloff et al (2009) Xenobiotica, 39:99-112; US patent 7,968,314



Longer Incubation Times to Enhance Sensitivity in Detecting TDI


Case Study with Diltiazem



Sequential Metabolism of Diltiazem is Responsible for Time-Dependent Inhibition of CYP3A



- N-desmethyl metabolite is a more potent inactivator than parent diltiazem
- Further oxidation to N-hydroxydesmethyl diltiazem leads to MIC and TDI

	K_{inact}	K_i
N-desmethyl diltiazem	0.047	1.1
Diltiazem	0.012	0.48



Adapted from Zhao et al (2007) Sequential Metabolism Is Responsible for Diltiazem-Induced Time-Dependent Loss of CYP3A. Drug Metab Dispos 35:704.

Sequential Metabolism Leading to TDI







- Conditions that promote sequential metabolism are expected to drive MIC formation
- This would include higher [I], higher protein (because $v \sim E$), and longer preincubation times

Effect of Increasing Preincubation Times in the Dilution and Non-dilution Method IC50 Shift Assay for Diltiazem for CYP3A4 with Testosterone and Midazolam as Substrates


Substrate	Precincubation time	Exp ¹	Dilution method			Non-dilution method		
			IC50 (- NADPH)	IC50 (+NADPH)	Shift	IC50 (- NADPH)	IC50 (+NADPH)	Shift
Midazolam	30 min	1	109	10	11	100	45	2.2
Testosterone	30 min	1	151	5.7	27	74	29	2.5
Midazolam	90 min	1	107	0.02	5711	98	13	7.5
Testosterone	90 min	1	197	0.04	4545	183	13	14
Midazolam	3 min	2	82	56	1.5	117	123	0.9
Midazolam	10 min	2	108	46	2.4	146	98	1.5
Midazolam	30 min	2	123	20	6.1	151	80	1.9
Midazolam	3 min	3	87	59	1.5	111	96	1.2
Midazolam	10 min	3	74	34	2.2	132	109	1.2
Midazolam	30 min	3	92	13	6.8	150	74	2.0
Midazolam	90 min	3	154	0.10	1493	234	41	5.7

¹ - Experiment number

Effect of Extended Incubation Times in the Dilution and Non-dilution Method IC50 Shift Assay - Diltiazem



Metabolism of 100 μ M Diltiazem in 0.1 or 1 mg/mL HLM




Confirmation that the longer incubation times as well as higher protein concentrations are driving metabolism




Other Examples

Inhibitor	Precincubation time	Exp ¹	Dilution method			Non-dilution method		
			IC50 (-NADPH)	IC50 (+NADPH)	Shift	IC50 (-NADPH)	IC50 (+NADPH)	Shift
Amiodarone	10 min	1	49	19	2.6	46	40	1.2
Amiodarone	30 min	1	51	6.7	7.5	61	22	2.7
Amiodarone	90 min	1	21	0.93	23	55	4.5	12
Amiodarone	10 min	2	81	34	2.4	40	25	1.6
Amiodarone	30 min	2	28	11	2.5	41	15	2.7
Amiodarone	90 min	2	13	0.8	16	30	4.0	7.5
Verapamil	3 min	1	26	15	1.7	41	31	1.3
Verapamil	10 min	1	29	7.5	3.9	40	20	2.0
Verapamil	30 min	1	26	0.54	48	37	10	3.6
Verapamil	3 min	2	21	11	1.9	29	21	1.4
Verapamil	10 min	2	18	3.7	4.8	28	12	2.2
Verapamil	30 min	2	23	0.32	70	30	4.1	7.2
Verapamil	90 min	2	28	0.03	962	25	1.2	21.0
Ketoconazole	3 min	1	0.0092	0.0113	0.8	0.0153	0.0144	1.1
Ketoconazole	10 min	1	0.0100	0.0114	0.9	0.0099	0.0107	0.9
Ketoconazole	30 min	1	0.0096	0.0126	0.8	0.0088	0.0113	0.8

¹ - Experiment number



Amiodarone




Verapamil

Negative control



AMG487

- A potent and selective CXCR3 antagonist, displayed a dose- and time-dependent reduction in oral clearance in a Phase I multiple dose clinical study.
- One explanation for this observation is time-dependent inhibition of CYP3A4, the enzyme primarily responsible for AMG487 metabolism
- The major phenol metabolite (M2), but not parent AMG487 shows TDI of CYP3A4 *in vitro*




Tonn GR, et al. (2009) An inhibitory metabolite leads to dose- and time-dependent pharmacokinetics of (R)-N-{1-[3-(4-ethoxyphenyl)-4-oxo-3,4-dihydro-pyrido[2,3-d]pyrimidin-2-yl]-ethyl}-N-pyridin-3-yl-methyl-2-(4-trifluoromethoxy-phenyl)-acetamide (AMG 487) in human subjects after multiple dosing Drug Metab Dispos 37:502-513.




AMG487 and M2 with 30 Min Preincubation



Midazolam



Testosterone



M2



Hypothesis: Driving metabolism with higher protein and longer preincubation time will permit detection of TDI of AMG487



AMG487 - Dilution and Substrate Addition Method (30 and 90 min Preincubation)

30 Minute Preincubation		Substrate Addition method			Dilution method		
Substrate	Inhibitor	IC50 (μM) (-NADPH)	IC50 (μM) (+NADPH)	IC50 Shift	IC50 (μM) (-NADPH)	IC50 (μM) (+NADPH)	IC50 Shift
Midazolam	AMG487	8.3	5.7	1.5	7	5.2	1.4
Testosterone	AMG487	20	22	0.9	26	25	1.1

90 Minute Preincubation		Substrate Addition method			Dilution method		
Substrate	Inhibitor	IC50 (μM) (-NADPH)	IC50 (μM) (+NADPH)	IC50 Shift	IC50 (μM) (-NADPH)	IC50 (μM) (+NADPH)	IC50 Shift
Midazolam	AMG487	9.2	2.9	3.2	8.0	2.7	2.9
Testosterone	AMG487	27	30	1.1	30	21	1.4

- Longer preincubation gave higher shifts
- No difference between dilution and non-dilution method

Henne et al (2012) Sequential Metabolism of AMG 487, a Novel CXCR3 Antagonist, Results in Formation of Quinone Reactive Metabolites that Covalently Modify CYP3A4 Cys239 and Cause Time-Dependent Inhibition of the Enzyme. Drug Metab Dispos. Manuscript in review

Summary and Conclusions

- The time-course of the “IC₅₀ shift” was examined and extended beyond the conventional 30 min preincubation period to examine the effect on assay sensitivity.
- All compounds tested exhibited increases in IC₅₀ shift when incubating for 90 min compared to 30 min.
- The increase occurred in both the dilution and non-dilution methods, but the increase was much larger for diltiazem and verapamil in the dilution method.

Summary and Conclusions

- In all cases, the increase in shift was attributable to a lowering of the +NADPH IC_{50} values.
- The finding of larger shifts at 90 min in both the dilution and non-dilution methods demonstrate enzyme inactivation is not complete at the conventional 30 min period
- The data suggest broad application of extended preincubation times represents a simple method to enhance sensitivity of the IC_{50} shift assay in detection of TDI and may improve opportunities to de-risk compounds

Acknowledgments

- Elke Perloff
- Andrew Mason
- Shangara Dehal
- Andy Blanchard
- Thuy Ho
- Charles Crespi
- Nathalie Boily
- Deqing Xiao
- Enne Akor
- Eric Gangl
- Kirk Henne (Amgen)





Questions?

Contact information

David M. Stresser, Ph.D.

e-mail: david_stresser@bd.com

or admetox@bd.com

Our services website: Click [here](http://www.bdbiosciences.com/research/admetox/services/index.jsp)

<http://www.bdbiosciences.com/research/admetox/services/index.jsp>

For research use only. Not intended for use in diagnostic or therapeutic procedures.
BD, BD Logo, and all other trademarks are property of Becton, Dickinson and Company. ©2010 BD

