

Helping all people
live healthy lives

Abstract

Overview


Objectives

This poster presents a compilation of our observations to date, and summarizes strategies for mitigating preanalytical variability.

Experimental Designs

- 
- BD™ P100*
- *For Research Use Only. Not For Diagnostic Use.
• Patented.

Procedures for peptide analysis



Results


```
graph TD; A[Blood Collection Process] --> B[Patient Information]; A --> C[Collection Device]; A --> D[Blood Derivative and Processing]; B --> B1[Gender]; B --> B2[Age]; B --> B3[Diet]; B --> B4[Genetics]; B --> B5[Medical Background]; B --> B6[Health Background]; B --> B7[Special conditions pregnancy, pre/post-menopausal, medications]; B --> B8[Social history alcohol intake, smoking status]; C --> C1[Tube or bag?]; C --> C2[Glass or plastic?]; C --> C3[Gel separator?]; C --> C4[Non-gel separator?]; C --> C5[Other tube additives?]; D --> D1[Needle gauge]; D --> D2[Blood collection set details of]; D --> D3[Phlebotomy]; D --> D4[Tourniquet technique?]; D --> D5[Patient position: seated/standing/lying]; D --> D6[First tube vs. last tube or track tube order]; D --> D7[Blood source: venipuncture or from existing line?];
```

Blood Collection Process

- **Patient Information**
 - Gender
 - Age
 - Diet
 - Genetics
 - Medical Background
 - Health Background
 - Special conditions (pregnancy, pre/post-menopausal, medications)
 - Social history (alcohol intake, smoking status)
- **Collection Device**
 - Tube or bag?
 - Glass or plastic?
 - Gel separator?
 - Non-gel separator?
 - Other tube additives?
- **Blood Derivative and Processing**
 - Needle gauge
 - Blood collection set (details of)
 - Phlebotomy
 - Tourniquet technique?
 - Patient position: seated/standing/lying
 - First tube vs. last tube (or track tube order)
 - Blood source: venipuncture or from existing line?

Adapted from (1)


Cell separation method affects plasma proteome content



- | Centrifuge time, min | 10 | 20 | 30 | 10 | 20 | 30 | 10 | 20 | 30 |
|----------------------|------|----|----|------|----|----|------|----|----|
| Centrifuge RCF | 1100 | | | 1600 | | | 2500 | | |


- The nature of the collection device itself can alter the sample protein content.
- This data represents another example of the need to carefully plan and track all aspects of blood sample acquisition for proteomics studies.

Variable peptide content in typical blood samples



- Serum and plasma were drawn from the same individual. Data shown for fresh "time zero" samples
- Serum contains significantly different peptides from plasma:
 - More peptides from 2000 to 3000 m/z.
 - Highly abundant peptides from 1000 to 2000 m/z.
 - The extra peptides are generated due to the *ex vivo* proteolysis during the clotting process^{2,3}.
- Differences between the three plasma samples are also observed.
- Since the nature of the blood sample imparts easily detected differences in the resulting peptide content, the choice of sample type is non-trivial, should be kept consistent throughout studies, and careful record-keeping of the nature of samples will be essential as part of later data analysis.


Protease inhibitors in plasma samples improve stability over time




- Variation in peptide content as a function of room-temperature incubation of EDTA plasma is evident. Similar extents of change are evident in citrate and heparin plasma and in serum³, suggesting a general instability of human plasma due to intrinsic proteolysis.
- The time-dependent variations of peptide content are reduced in P100 samples.
- P100 sample is reasonably stable for > 2 hours after collection, suggesting improved handling conditions including:
 - Use of protease-inhibited plasma, such as provided by BD P100.
 - Analysis of the sample within a 2-hour timeframe, or
 - Freezing the collected sample within 2 hours.

Stabilization by Protease Inhibitors: a “deeper look” with LC/MALDI

EDTA




P100




- LONALDI methods yield a virtual plot (HPLC retention time vs. m/z), which shows significantly more peptides than direct MALDI shown in Figure 4. Data shown for one individual, with blood drawn into EDTA and P100 tubes. Peptides were prepared from "time zero" plasma, or from aliquots allowed to incubate at room temperature for one hour.
- Consistent with direct MALDI analyses, the EDTA samples have more peptides at "time zero" than P100. In addition, the P100 sample is noticeably more stable over time than the EDTA sample.
- Protease inhibitors provide stability over time to plasma peptide content.

Example with a specific peptide: bradykinin

- Bradykinin is a 9 residue peptide, shown to be a marker of cardiac trauma. Presumably at low levels in healthy blood samples, it has been described as being artifactually high, to the point of uselessness, in drawn blood samples, due to activity of kallikrein cleaving high molecular weight kininogen⁴.



- EDTA and P100 samples from one healthy subject. Peptides from "time zero" plasma.
- Magnified view of MALDI spectra show that there are fewer low mass peptides in P100 than in EDTA sample, suggesting that P100 inhibits the generation of these peptides.
- Of particular note is bradykinin (red arrow), which is virtually absent in P100 samples, but easily detected in EDTA plasma.



- EDTA and P100 samples from a healthy subject, peptides isolated after room temperature plasma incubation times indicated. MALDI spectra at the same reverse-phase retention time are shown, providing a very clear view of bradykinin.
- LC/MALDI spectra as a function of time demonstrate that bradykinin concentration remains nearly undetectable in P100 plasma for at least 120 minutes.
- By comparison, EDTA samples maintain approximately 10x higher bradykinin levels from the earliest time points.
- These data suggest that P100 maintains a low, *in vitro* level for at least 2 hours, possibly enabling better utility for bradykinin in drawn blood samples.

Conclusions

- Data demonstrate that many preanalytical variables (how blood samples are collected and handled) result in measurable changes in subsequent proteomics analyses.
- Reproducibility of data relies on reducing preanalytical variables, and many details of blood sample collection are important, often more so than researchers recognize.
- Every potential source of variation should be tracked, and minimized when possible.

- Choice of specific blood sample type directly influences observed proteins and peptides, both in “fresh” samples and also over time.

- Variations in residual cellular contamination appear to alter detectable protein content, with the least amount of degradation apparently associated with lowest cell counts.

- Inclusion of protease inhibitors in the plasma sample provide increased stability, with demonstrable improvement both at “time zero” and over at least 2 hours of room temperature incubation.

References

- (1) Rai, A. J. et al. *Proteomics* (2005) 5:3262-3277.
- (2) Villanueva, J. et al. *J. Clin. Invest.* (2006) 116:271-284.
- (3) Yi, J. et al. *J. Proteome Res.* (2007), in press.
- (4) Campbell, D. J. *Braz. J. Med. Biol. Res.* (2000) 33:665-677.