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
1. Abrams C, Shattil SJ. Immunological detection of activated platelets in clinical disorders. *Thrombosis and Haemostasis*. 1991;65:467-473.
2. Shattil SJ, Cunningham M, Hoxie JA. Detection of activated platelets in whole blood using activation-dependent monoclonal antibodies and flow cytometry. *Blood*. 1987;70:307-315.
3. Abrams CS, Ellison N, Budzinski AZ, Shattil SJ. Direct detection of activated platelets and platelet-derived microparticles in humans. *Blood*. 1990;75:128-138.
4. Jennings LK, Ashmun RA, Wang WF, Dockter ME. Analysis of human platelet glycoproteins IIb-IIIa and Glanzmann's thrombasthenia in whole blood by flow cytometry. *Blood*. 1986;68:173-179
5. Jackson AL, Warner NL. Preparation, staining and analysis by flow cytometry of peripheral blood leukocytes. In: Rose NR, Friedman H, Fahey JL, eds. *Manual of Clinical Laboratory Immunology*. 3rd ed. Washington, DC: American Society for Microbiology; 1986:226-235
6. Shattil SJ, Motulsky HJ, Insel PA, Flaherty L, Brass LF. Expression of fibrinogen receptors during activation and subsequent desensitization of human platelets by epinephrine. *Blood*. 1986;68:1224-1231.
7. Jackson AL, Warner NL. Preparation, staining, and analysis by flow cytometry of peripheral blood leukocytes. In: Rose NR, Friedman H, Fahey JL, eds. *Manual of Clinical Laboratory Immunology*, 3rd ed. Washington, DC: American Society for Microbiology; 1986:226-235.
8. Slockblower J, Belgeri K, Bruck E, et al. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture (H3-A2)*. Villanova, PA: The National Committee for Clinical Laboratory Standards, 1984.
9. Fijnheer R, Modderman PW, Veldman H, et al. Detection of platelet activation with monoclonal antibodies and flow cytometry. *Transfusion*. 1990;30:20-25.

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
Quadrant Statistics

X Parameter: FL1-H PAC-1 (Log)
Y Parameter: FL2-H CD62P PE (Log)
Quad Location: 15, 18



Quad	Events	% Gated	Y Mean	Y Geo Mean
UL	370	2.24	79.38	54.23
UR	3243	19.60	71.34	50.65
LL	5290	31.98	4.45	3.28
LR	7639	46.18	5.44	3.78

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