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1. Froland, S.S., aig., B., BerdalP. (1971) Surface-bound immunoglobulin as a Marker for B Lymphocytes. *Nature* 234, 251.
2. Warner NL, (1974) Membrane Immunoglobulins and antigen receptors on B & T lymphocytes *Advances Immunol.* 19, 67.
3. Kohler G., Milstein C. (1975) Continuous Cultures of Fused Cells Secreting Antibody of Predefined Specificity. *Nature* 256, 495.
4. Cantor, H., Boyse, E.A. (1975) Functional Subclasses of T Lymphocytes Bearing Different Ly Antigens. I. The Generation of Functionally Distinct T Cell Subclasses in a Differentiative Process Independent of Antigen. *J Exp Med* 141, 1376.
5. Meeker, T.C., Miller, R.A., Link, M.P., Bindl, J. Warneke, R, Levy, R. (1984) A Unique B Lymphocyte Antigen Defined by a Monoclonal Antibody. *Hybridoma* 3, 305.
6. Moldenhauer G, Dorken B, Schwartz R, Pezzutto A, Knops J, Hammerling GJ. Analysis of ten B lymphocyte-specific workshop monoclonal antibodies. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human B Lymphocytes, Volume II*. New York: Springer-Verlag; 1986:61-67.
7. Nadler L. B Cell/Leukemia Panel Workshop: Summary and comments. In: Reinherz E. Haynes B, Nadler L, Bernstein I, eds. *Leukocyte Typing II: Human B Lymphocytes, Volume II*. New York: Springer-Verlag; 1986:3-43.
8. Reinhert, E.L., Haynes, B.F., Nadler, L.M., Bernstein, L.D., (eds). (1985) In *Leucocyte Typing II, vol. 2: Human B Lymphocytes* (Springer-Verlag, New York).
9. Loken, M.R., Shah, V.O., Dattilio, K.L., Civin, C.I. (1987) Flow Cytometric Analysis of Human Bone Marrow II. Normal B Lymphocyte Development. *Blood* 70, 1316
10. Ryan, D., Kossover, S., Mitchel, S., Frantz, C., Hennessy, L., Cohen, H. (1986) Subpopulations of Common Acute Lymphoblastic Leukemia Antigenpositive Lymphoid cells in Normal Bone Marrow identified by Hematopoietic differentiation antigens. *Blood* 68, 417.
11. Duen B, Moller P, Pezzutto A, Schwartz-Albiez R, Molderhauer G. B-cell antigens: CD20. In: Knapp W, Duen B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. Oxford: Oxford University Press; 1989:46-48.
12. Hultin LE, Hausner MA, Hultin PM, Giorgi JV. CD20 (pan-B cell) antigen is expressed at a low level on a subpopulation of human T lymphocytes. *Cytometry*. 1993;14:196-204.
13. Lampson LA, Levy R. Two populations of Ia-like molecules on a human B cell line. *J Immunol.* 1980;125:393.
14. Brodsky FM. A matrix approach to human class II histocompatibility antigens: Reactions of four monoclonal antibodies with the products of nine haplotypes. *Immunogenetics*. 1984;19:179.
15. Robbins PA, Evans EL, Ding AH, Warner NL, Brodsky FM. Monoclonal antibodies that distinguish between class II antigens (HLA-DP, DQ, and DR) in 14 haplotypes. *Human Immunol.* 1987;18:301.
16. Thompson BE, Wagner DK, Nelson DL, Sullivan JL. Activated lymphocytes during acute Epstein-Barr virus infection. *J Immunol.* 1987;139:3802-3807.
17. Levacher M, Tallet S, Dazza M, Dournon E, Rouveix B, Pocardalo J. T activation marker evaluation in ARC patients treated with AZT. Comparison with CD4+ lymphocyte count in non-progressors and progressors towards AIDS. *Clin Exp Immunol.* 1990;81:177-182.
18. Stites D. Casavant C, McHugh T, et al. Flow cytometric analysis of lymphocyte phenotypes in IDS using monoclonal

- antibodies and simultaneous dual immunofluorescence. Clin Immunol Immunopathol. 1986;38:161-177.
19. Terstappen LWMM, Hollander Z, Meiners H, Loken M. Quantitative comparison of myeloid antigens on five lineages of mature peripheral blood cells. J Leuk Biol. 1990;48:138-148.
 20. Warnke R, Levy R. Detection of T and B cell antigens with hibridoma monoclonal antibodies: A biotin-avidin-horseradish peroxidase method. J Histochem cytochem. 1980;28:771.
 21. Warnke R, Miller R, Grogan T, Pederson M, Dilley J, Levy R. Immunologic phenotype in 30 patients with diffuse large-cell lymphoma. N Eng J Med. 1980;303:293.
 22. Engleman EG, Warnke R, Fox RI, Levy R. Studies of a human T lymphocyte antigen recognized by a monoclonal antibody. Proc Natl Acad Sci USA. 1981;78:1791-1795.
 23. Zipf TF, Fox RI, Dilley J, Levy R. Definition of the high risk ALL patient by immunologic phenotyping with monoclonal antibodies. Cancer Research. 1981;41:4786.

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