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1. Fitzgerald-Bocarsly P, Herberman R, Hercend T, et al. A definition of natural killer cells, In: Ades E, Lopez C, eds. Natural Killer Cells and Host Defense. Basel: Karger; 1989:1.
2. Jackson A. Basic phenotyping of lymphocytes: Selection and testing of reagents and interpretation of data. Clin Immunol Newslett. 1990;10:43-55.
3. Kidd P, Vogt R. Report of the Workshop on the evaluation of T-cell subsets during HIV infection and AIDS. Clin Immunol Immunopathol. 1989;52:3-9.
4. Landay A, Muirhead K. Procedural guidelines for performing immunophenotyping by flow cytometry. Clin Immunol Immunopathol. 1989;52:48-60.
5. Schmidt R. Monoclonal antibodies for diagnosis of immunodeficiencies. Blut. 1989;59:200-206.
6. Nicholson J. Use of flow cytometry in the evaluation and diagnosis of primary and secondary immunodeficiency diseases. Arch Pathol Lab Med. 1989;113:598-605.
7. Foucar K, Gueken Jr. Clinical application of immunologic techniques to the diagnosis of lymphoproliferative and immunodeficiency disorders. Lab Medicine. 1982;13:403.
8. Cohen S, Weetman A. Activated interstitial and intraepithelial thyroid lymphocytes in autoimmune thyroid disease. Acta Endocr. 1988;119:161-166.
9. Smolen J, Chused T, Leiserson W, Reeves J, Alling D, Steinberg A. Heterogeneity of

immunoregulatory T-cell subsets in systemic lupus erythematosus. Am J Med. 1982;72:783-790.

10. Giorgi J, Hultin L. Lymphocyte subset alterations and immunophenotyping by flow cytometry in HIV disease. Clin Immunol Newslett. 1990;10:55-62.
11. Antel J, Bania M, Noronha A, Neely S. Defective suppressor cell function mediated by T8+ cell lines from patients with progressive multiple sclerosis. J Immunol. 1986;137:3436-3439.
12. Gratama J, Naipal A, Olijans P, et al. T lymphocyte repopulation and differentiation after bone marrow transplantation: Early shifts in the ratio between T4+ and T8+ T lymphocytes correlate with the occurrence of acute graft-versus-host disease. Blood. 1984;63(6):1416-1423.
13. Bishop G, Hall B, Duggin G, Horvath J, Sheil A, Tiller D. Immunopathology of renal allograft rejection analyzed with monoclonal antibodies to mononuclear cell markers. Kidney Internat. 1986;29:708-717.
14. Wolde-Mariam W, Peter J. Recent diagnostic advances in cellular immunology. Diagnost Med. 1984;7:25-32.
15. Landay A, Ohlsson-Wilhelm B, Giorgi J. Application of flow cytometry to the study of HIV infection. AIDS. 1990;4:479-497.
16. Cobbold S, Hale G, Waldmann H. Non-lineage, LFA-1 family, and leucocyte common antigens: New and previously defined clusters. In: McMichael A, ed. Leucocyte Typing III: White Cell Differentiation Antigens. Oxford: Oxford University Press; 1986:788-803.
17. Beverley P. Production and use of monoclonal antibodies in transplantation immunology. In: Touraine J, Trager J, Betuel H, eds. Transplantation and Clinical Immunology XI. Amsterdam: Excerpta Medica; 1980:87-94.
18. Rowe D, Beverley P. Characterization of breast cancer infiltrates using monoclonal antibodies to human leucocyte antigens. Br J Cancer. 1984;49:149-159.
19. Bernard A, Boumsell L, Hill C. Joint report of the First International Workshop on Human Leucocyte Differentiation Antigens by the investigators of the participating laboratories: T2 protocol. In: Bernard A, Boumsell L, Dausset J, Milstein C, Schlossman S, eds. Leucocyte Typing. Berlin: Springer-Verlag; 1984:25-60.
20. Dimitriu-Bona A, Burmester G, Waters S, Winchester R. Human mononuclear phagocyte differentiation antigens. I. Patterns of antigenic expression on the surface of human monocytes and macrophages defined by monoclonal antibodies. J Immunol 1983;130(1):145-152.
21. Herrmann F, Komischke B, Odenwald E, Ludwig W. Use of monoclonal antibodies as a diagnostic tool in human leukemia. I. Acute myeloid leukemia and acute phase of chronic myeloid leukemia. Blut. 1983;47:157-163.
22. Schwitzer R. Cluster report: CD45/CD45R. In: Knapp W, Duken B, Gilks W, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. Oxford University Press; 1989:628-634.
23. Goyert S, Tesio L, Ashman L, et al. Report on the CD14 Cluster Workshop. In: Knapp W, Duken B, Gilks W, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. Oxford: Oxford University Press; 1989:789-794.
24. Bernstein I, Self S. Joint report of the myeloid section of the Second International Workshop on Human Leucocyte Differentiation Antigens. In: Reinherz E, Haynes b, Nadler L, Bernstein I, eds. Leucocyte Typing II: Human Myeloid and Hematopoietic Cells. New York: Springer-Verlag; 1984;3:1-25.
25. Jayaram Y, Hogg N. Surface expression of CD14 molecules on human neutrophils. In: Knapp W, Duken B, Gilks W, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. Oxford University Press; 1989:796-797.
26. Haynes B. Summary of T cell studies performed during the Second International Workshop and Conference on Human Leucocyte Differentiation Antigens. In: Reinherz E, Haynes B, Nadler L, Bernstein I, eds. Leucocyte Typing II: Human T Lymphocytes. New York: Springer-Verlag; 1986:1:3-30.
27. Ledbetter J, Evans R, Lipinski M, Cunningham-Rundles C, Good R, Herzenberg L. Evolutionary conservation of surface molecules that distinguish T lymphocyte helper/inducer and T cytotoxic/suppressor subpopulations in mouse and man. J Exp med. 1981;153 (February):310-323.
28. Kan E, Wang C, Wang L, Evans R. Non-covalently bonded subunits of 22 and 28 kd are rapidly internalized by T cells reacted with anti-Leu-4 antibody. J Immunol. 1983;131(2):536-539.
29. Knowles R. Immunochemical analysis of the T cell-specific antigens. In: Reinherz E, Haynes B, Nadler L, Bernstein I, eds. Leucocyte Typing II: Human T Lymphocytes. New York: Springer-Verlag; 1986:259-288.
30. Nadler L, B cell/leukemia panel workshop: Summary and comments. In: Reinherz E, Haynes B, Nadler L, Bernstein I, eds. Leucocyte Typing II: Human B Lymphocytes. New York: Springer-Verlag; 1986:2:3-43.
31. van Dongen J, Krissansen G, Wolvers-Tettero I, et al. Cytoplasmic expression of the CD3 antigen as a diagnostic marker for immature T-cell malignancies. Blood 1988;71(3):603-612.
32. Brenner M, Groh V, Porcelli F, Hochstenbach H, et al. Structure and distribution of the human gd T-cell receptor. In: Knapp W, Duken B, Gilks W, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. Oxford: Oxford University Press; 1989:1049-1053.
33. Duken B, Moller P, Pezzutto A, Schwartz-Albiez R, Molderhauer G. B cell antigens: CD19. In: Knapp W, Duken B, Gilks W, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. Oxford: Oxford University Press; 1986:34-36.
34. Clevers H, Alarcón B, Wileman T, Terhorst C. The T cell receptor/CD3 complex: a dynamic protein ensemble. Annual Rev Immunol. 1988;(6):629.
35. Evans R, Wall D, Platsoucas C, et al. Thymus-dependent membrane antigens in man: Inhibition of cell-mediated lympholysis by monoclonal antibodies to the TH2 antigen. Proc Natl Acad Sci USA. 1981;78(1):544-548.
36. Maddon P, Dalglish A, McDougal J, Clapham P, Weiss R, Axel R. The T4 gene encodes the AIDS virus receptor and is expressed in the immune system and the brain. Cell. 1986;47:333-348.
37. Dalglish A, Beverly P, Clapham P, Crawford D, Greaves M, Weiss R. The CD4 (T4) antigen is an essential component of the receptor for the AIDS virus. Nature; 1984;312(December):763-767.
38. Kurlle R. Cluster report: CD3, In: Knapp W, Duken B, Gilks W, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. Oxford: Oxford University Press; 1989:290-293.
39. Wood F, Warner N, Warnke R. Anti-Leu-3/T4 antibodies react with cells of monocyte/macrophage and Langerhans lineage. J Immunol. 1983;131(1):212-216.
40. Lanier L, Le A, Phillips J, Warner N, Babcock G. Subpopulations of human natural killer cells defined by expression of the Leu-7 (HNK-1) and Leu-11 (NK-15) antigens. J Immunol. 1983; 131(4):1789-1796.
41. Anderson P, Blue M-L, Morimoto C, Schlossman S. Cross-linking of T3 (CD3) with T4 (CD4) enhances the proliferation of resting T lymphocytes. J Immunol. 1987;139:679-682.
42. Eichmann K, Johnson J-I, Falk I, Emrich F. Effective activation of resting mouse T lymphocytes

- by cross-linking submitogenic concentration of the T cell antigen receptor with either Lyt-2 or L3T4. *Eur J Immunol.* 1987;17:643-650.
43. Gallagher P, Fazekas de St. Groth B, Miller J. CD4 and CD8 molecules can physically associate with the same T-cell receptor. *Proc Natl Acad Sci USA.* 1989;86:10044-10048.
 44. Rudd C, Burgess K, Barber E, Schlossman S. Monoclonal antibodies to the CD4 and CD8 antigens precipitate variable amounts of CD4/CD8-associated p56-lck activity. In: Knapp W, Duken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens.* Oxford: Oxford University Press; 1989:326-327.
 45. Terry L, Disanto J, Small T, Flomenberg N. Differential expression of the CD8 and Lyt-3 antigens on a subset of human T-cell receptor g/d-bearing lymphocytes. In: Knapp W, Duken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens.* Oxford: Oxford University Press; 1989; 345-346.
 46. Moebius U. Cluster report: CD8. In: Knapp W, Duken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigen.* Oxford: Oxford University Press; 1989:342-343.
 47. Perussia B, Starr S, Abraham S, Fanning V, Trinchieri G. Human natural killer cells analyzed by B73.1 a monoclonal antibody blocking Fc receptor functions. I. Characterization of the lymphocyte subset reactive with B73.1. *J Immunol.* 1983;130(5):2133-2141.
 48. Perussia B, Acuto O, Terhorst C, et al. Human natural killer cells analyzed by B73.1, a monoclonal antibody blocking Fc receptor functions, II. Studies of B73.1 antibody-antigen interaction on the lymphocyte membrane. *J Immunol.* 1983;130(5):2142-2148.
 49. Perussia B, Trinchieri G, Jackson A, et al. The Fc receptor for IgG on human natural killer cells: Phenotypic, functional, and comparative studies with monoclonal antibodies, *J Immunol.* 1984;133(1):180-189.
 50. Schmidt R. Non-lineage/natural killer cluster report: New and previously defined clusters. In: Knapp W, Duken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens.* Oxford: Oxford University Press; 1989:517-542.
 51. Schubert J, Lanier L, Schmidt R. Cluster report: CD56. In: Knapp W, Duken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens.* Oxford: Oxford University Press; 1989:699-702.
 52. Lanier L, Testi R, Bindl J, Phillips J. Identity of Leu-19 (CD56) leukocyte differentiation antigen and neural cell adhesion molecule. *J Exp Med.* 1989;169(June):2233-2238.
 53. Lanier L, Chang C, Azuma M, Ruitenberg J, Hemperly J, Phillips J. Molecular and functional analysis of human natural killer cell-associated neural cell adhesion molecule (N-CAM/CD56). *J Immunol.* 1991;146:4421-4426.
 54. Lanier L, Le A, Civin C, Loken M, Phillips J. The relationship of CD16 (Leu-11) and Leu-19 (NKH-1) antigen expression on human peripheral blood NK cells and cytotoxic T lymphocytes. *J Immunol.* 1986;136(12):4480-4486.
 55. Phillips J, Lanier L. Dissection of the lymphokine-activated killer phenomenon: Relative contribution of peripheral blood natural killer cells and T lymphocytes to cytotoxicity. *J Exp Med.* 1986;164:814-825.
 56. Giorgi J. Lymphocyte subset measurements: Significance in clinical medicine. In: Rose N, Friedman H, Fahey J, eds. *Manual of Clinical Laboratory Immunology.* 3rd ed. Washington DC: American Society of Microbiology. 1986:236-246.
 57. Angadi C. Lack of Leu-3a epitope on T-helper (CD4) lymphocytes. *J Clin Lab Anal.* 1990;4:193-195.
 58. Prince H, Hirji K, Waldbeser L, Plaeger-Marshall S, Kleinman S, Lanier L. Influence of racial background on the distribution of T cell subsets and Leu-11-positive lymphocytes in healthy blood donors. *Diag Immunol.* 1985;3:33-37.

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