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1. Fitzgerald-Bocarsly P, Herberman R, Hercend T, et al. A definition of natural killer cells. In: Ades E, Lopez C, eds. Natural Killer Cells and Host Defense. Basel: Karger; 1989:1.
2. Brain P, Gordon J, Willets W. Rosette formation by peripheral lymphocytes. Clin Exp Immunol. 1970;6:681-688.
3. Jondal M, Holm G, Wigzell H. Surface markers on human T and B lymphocytes. I. A large population of lymphocytes forming nonimmune rosettes with sheep red blood cells. J Exp Med. 1972;136:207-215.
4. Wolde-Mariam W, Peter J. Recent diagnostic advances in cellular immunology. Diagnost Med. 1984;7:25-32.
5. Froland S, Natvig J, Bernald P. Surface-bound immunoglobulin as a marker of B lymphocytes in man. Nature. 1971;234:251-252.
6. Bellanti J, Peters S. Diagnostic applications of immunology. In: Bellanti J, ed. Immunology III. Philadelphia: WB Saunders Co; 1985:558-569.
7. Winchester R, Fu S, Hoffman H, Kunkel H. IgG on lymphocyte surfaces: technical problems and the significance of a third cell population. J Immunol. 1975;114:1210-1212.
8. Wofsy L, Henry C, Kimura J, North J. Modification and use of antibodies to label cell surface antigens. In: Mishell B, Shiigi S, eds. Selected Methods in Cellular Immunology. New York, NY: WH Freeman and Co; 1980:287-304.
9. Henry C, Mishell B, Shiigi S. Hemolytic plaque assays. In: Mishell B, Shiigi S, eds. Selected Methods in Cellular Immunology. New York, NY: WH Freeman and Co; 1980:69-123.
10. Grabstein K, Chen Y. Cell-mediated cytolytic responses. In: Mishell B, Shiigi S, eds. Selected Methods in Cellular Immunology. New York, NY: WH Freeman and Co; 1980:124-137.
11. Oi V, Herzenberg L. Immunoglobulin-producing hybrid cell lines. In: Mishell B, Shiigi S, eds. Selected Methods in Cellular Immunology. New York, NY: WH Freeman and Co; 1980:351-372.
12. Jackson A. Basic phenotyping of lymphocytes: selection and testing of reagents and interpretation of data. Clin Immunol Newslett. 1990;10:43-55.
13. Jackson A, Warner N. Preparation, staining, and analysis by flow cytometry of peripheral blood leukocytes. In: Rose N, Friedman H, Fahey J, eds. Manual of Clinical Laboratory Immunology. 3rd ed. Washington DC: American Society for Microbiology; 1986:226-235.
14. Nicholson J, Jones B. Immunophenotyping by flow cytometry: its use in HIV infection. Labmedica. 1989;6:21-26.
15. Kidd P, Vogt R. Report of the workshop on the evaluation of T-cell subsets during HIV infection and AIDS. Clin Immunol Immunopathol. 1989;52:3-9.
16. Landay A, Muirhead K. Procedural guidelines for performing immunophenotyping by flow cytometry. Clin Immunol Immunopathol. 1989;52:48-60.
17. Schmidt RE. Monoclonal antibodies for diagnosis of immunodeficiencies. Blut. 1989;59:200-206.
18. Nicholson J. Use of flow cytometry in the evaluation and diagnosis of primary and secondary immunodeficiency diseases. Arch Pathol Lab Med. 1989;113:598-605.
19. Foucar K, Gueken JA. Clinical application of immunologic techniques to the diagnosis of lymphoproliferative and immunodeficiency disorders. Lab Medicine. 1982;13:403-413.
20. Cohen S, Weetman A. Activated interstitial and intraepithelial thyroid lymphocytes in autoimmune thyroid disease. Acta Endocr. 1988;119:161-166.
21. Smolen J, Chused T, Leiserson W, Reeves J, Alling D, Steinberg A. Heterogeneity of immunoregulatory T-cell subsets in systemic lupus erythematosus. Am J Med. 1982;72:783-790.
22. Giorgi J, Hultin L. Lymphocyte subset alterations and immunophenotyping by flow cytometry in HIV disease. Clin Immunol Newslett. 1990;10:55-62.
23. Antel J, Bania M, Noronha A, Neely S. Defective suppressor cell function mediated by T8+ cell lines from patients with progressive multiple sclerosis. J Immunol. 1986;137:3436-3439.
24. Gratama J, Naipal A, Olijans P, et al. T lymphocyte repopulation and differentiation after bone marrow transplantation: early shifts in the ratio between T4+ and T8+ T lymphocytes correlate with the occurrence of acute graft-versus-host disease. Blood. 1984;63:1416-1423.
25. Bishop G, Hall B, Duggin G, Horvath J, Sheil A, Tiller D. Immunopathology of renal allograft rejection analyzed with monoclonal antibodies to mononuclear cell markers. Kidney Internat. 1986;29:708-717.
26. Lanier L, Le A, Phillips J, Warner N, Babcock G. Subpopulations of human natural killer cells defined by expression of the Leu-7 (HNK-1) and Leu-11 (NK-15) antigens. J Immunol. 1983;131:1789-1796.
27. Landay A, Ohlsson-Wilhelm B, Giorgi J. Application of flow cytometry to the study of HIV infection. AIDS. 1990;4:479-497.
28. Cobbold S, Hale G, Waldmann H. Non-lineage, LFA-1 family, and leucocyte common antigens: new and previously defined clusters. In: McMichael A, ed. Leucocyte Typing III: White Cell Differentiation Antigens. New York, NY: Oxford University Press; 1986:788-803.
29. Beverley P. Production and use of monoclonal antibodies in transplantation immunology. In: Touraine J, Trager J, Betuel H, eds. Transplantation and Clinical Immunology XI. Amsterdam: Excerpta Medica; 1980:87-94.
30. Rowe D, Beverley P. Characterisation of breast cancer infiltrates using monoclonal antibodies to human leucocyte antigens. Br J Cancer. 1984;49:149-159.
31. Bernard A, Boumsell L, Hill C. Joint report of the First International Workshop on Human Leucocyte Differentiation Antigens by the investigators of the participating laboratories: T2 protocol. In: Bernard A, Boumsell L, Dausset J, Milstein C, Schlossman S, eds. Leucocyte Typing. Berlin: Springer-Verlag; 1984:25-60.
32. Dimitriu-Bona A, Burmester G, Waters S, Winchester R. Human mononuclear phagocyte differentiation antigens. I. Patterns of antigenic expression on the surface of human monocytes and macrophages defined by monoclonal antibodies. J Immunol. 1983;130(1):145-152.
33. Herrmann F, Komischke B, Odenwald E, Ludwig W. Use of monoclonal antibodies as a diagnostic tool in human leukemia. I. Acute myeloid leukemia and acute phase of chronic myeloid leukemia. Blut. 1983;47:157-163.
34. Schwizner R. Cluster report: CD45/CD45R. In: Knapp W, Dörken B, Gilks W, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. New York, NY: Oxford University Press; 1989:628-634.
35. Goyert S, Tesio L, Ashman L, et al. Report on the CD14 Cluster Workshop. In: Knapp W, Dörken B, Gilks W, et al, eds. Leucocyte Typing IV: White Cell Differentiation Antigens. New York, NY: Oxford University Press; 1989:789-794.
36. Bernstein I, Self S. Joint report of the myeloid section of the Second International Workshop on Human Leucocyte Differentiation Antigens. In: Reinherz E, Haynes B, Nadler L, Bernstein I, eds. Leucocyte Typing II: Human Myeloid and Hematopoietic Cells. New York, NY: Springer-Verlag;

- 1984;3:1-25.
37. Jayaram Y, Hogg N. Surface expression of CD14 molecules on human neutrophils. In: Knapp W, Dörken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:796-797.
 38. Haynes B. Summary of T cell studies performed during the Second International Workshop and Conference on Human Leukocyte Differentiation Antigens. In: Reinherz E, Haynes B, Nadler L, Bernstein I, eds. *Leucocyte Typing II: Human T Lymphocytes*. New York, NY: Springer-Verlag; 1986:1:3-30.
 39. Ledbetter J, Evans R, Lipinski M, Cunningham-Rundles C, Good R, Herzenberg L. Evolutionary conservation of surface molecules that distinguish T lymphocyte helper/inducer and T cytotoxic/suppressor subpopulations in mouse and man. *J Exp Med*. 1981;153:310-323.
 40. Kan E, Wang C, Wang L, Evans R. Non-covalently bonded subunits of 22 and 28 kd are rapidly internalized by T cells reacted with Anti-Leu-4 antibody. *J Immunol*. 1983;131:536-539.
 41. Knowles R. Immunohistochemical analysis of the T cell-specific antigens. In: Reinherz E, Haynes B, Nadler L, Bernstein I, eds. *Leucocyte Typing II: Human T Lymphocytes*. New York, NY: Springer-Verlag; 1986:259-288.
 42. Nadler L. B cell/leukemia panel workshop: summary and comments. In: Reinherz E, Haynes B, Nadler L, Bernstein I, eds. *Leucocyte Typing II: Human B Lymphocytes*. New York, NY: Springer-Verlag; 1986:2:3-43.
 43. van Dongen J, Krissansen G, Wolvers-Tettero I, et al. Cytoplasmic expression of the CD3 antigen as a diagnostic marker for immature T-cell malignancies. *Blood*. 1988;71:603-612.
 44. Brenner M, Groh V, Porcelli F, Hochstenbach H, et al. Structure and distribution of the human γ δ T-cell receptor. In: Knapp W, Dörken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:1049-1053.
 45. Clevers H, Alarcón B, Wileman T, Terhorst C. The T cell receptor/CD3 complex: a dynamic protein ensemble. *Annual Rev Immunol*. 1988;6:629-662.
 46. Dörken B, Moller P, Pezzutto A, Schwartz-Albiez R, Molderhauer G. B-cell antigens: CD19. In: Knapp W, Dörken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1986:34-36.
 47. Evans R, Wall D, Platsoucas C, et al. Thymus-dependent membrane antigens in man: inhibition of cell-mediated lympholysis by monoclonal antibodies to the TH2 antigen. *Proc Natl Acad Sci USA*. 1981;78(1):544-548.
 48. Maddon P, Dalgleish A, McDougal J, Clapham P, Weiss R, Axel R. The T4 gene encodes the AIDS virus receptor and is expressed in the immune system and the brain. *Cell*. 1986;47:333-348.
 49. Dalgleish A, Beverly P, Clapham P, Crawford D, Greaves M, Weiss R. The CD4 (T4) antigen is an essential component of the receptor for the AIDS virus. *Nature*. 1984;312:763-767.
 50. Kurrel R. Cluster report: CD3. In: Knapp W, Dörken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:290-293.
 51. Wood F, Warner N, Warnke R. Anti-Leu-3/T4 antibodies react with cells of monocyte/macrophage and Langerhans lineage. *J Immunol*. 1983;131:212-216.
 52. Anderson P, Blue M-L, Morimoto C, Schlossman S. Cross-linking of T3 (CD3) with T4 (CD4) enhances the proliferation of resting T lymphocytes. *J Immunol*. 1987;139:678-682.
 53. Eichmann K, Johnson J-I, Falk I, Emmrich F. Effective activation of resting mouse T lymphocytes by cross-linking submitogenic concentrations of the T cell antigen receptor with either Lyt-2 or L3T4. *Eur J Immunol*. 1987;17:643-650.
 54. Gallagher P, Fazekas de St. Groth B, Miller J. CD4 and CD8 molecules can physically associate with the same T-cell receptor. *Proc Natl Acad Sci USA*. 1989;86:10044-10048.
 55. Rudd C, Burgess K, Barber E, Schlossman S. Monoclonal antibodies to the CD4 and CD8 antigens precipitate variable amounts of CD4/CD8-associated p56-lck activity. In: Knapp W, Dörken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989: 326-327.
 56. Terry L, Disanto J, Small T, Flomenberg N. Differential expression of the CD8 and Lyt-3 antigens on a subset of human T-cell receptor γ/δ -bearing lymphocytes. In: Knapp W, Dörken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:345-346.
 57. Moebius U. Cluster report: CD8. In: Knapp W, Dörken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigen*. New York, NY: Oxford University Press; 1989:342-343.
 58. Perussia B, Starr S, Abraham S, Fanning V, Trinchieri G. Human natural killer cells analyzed by B73.1, a monoclonal antibody blocking Fc receptor functions. I. Characterization of the lymphocyte subset reactive with B73.1. *J Immunol*. 1983;130(5):2133-2141.
 59. Perussia B, Acuto O, Terhorst C, et al. Human natural killer cells analyzed by B73.1, a monoclonal antibody blocking Fc receptor functions. II. Studies of B73.1 antibody-antigen interaction on the lymphocyte membrane. *J Immunol*. 1983;130:2142-2148.
 60. Perussia B, Trinchieri G, Jackson A, et al. The Fc receptor for IgG on human natural killer cells: phenotypic, functional, and comparative studies with monoclonal antibodies. *J Immunol*. 1984;133:180-189.
 61. Schmidt R. Non-lineage/natural killer section report: new and previously defined clusters. In: Knapp W, Dörken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:517-542.
 62. Schubert J, Lanier L, Schmidt R. Cluster report: CD56. In: Knapp W, Dörken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:699-702.
 63. Lanier L, Testi R, Bindl J, Phillips J. Identity of Leu-19 (CD56) leukocyte differentiation antigen and neural cell adhesion molecule. *J Exp Med*. 1989;169:2233-2238.
 64. Lanier L, Chang C, Azuma M, Ruitenberg J, Hemperly J, Phillips J. Molecular and functional analysis of human natural killer cell-associated neural cell adhesion molecule (N-CAM/CD56). *J Immunol*. 1991;146:4421-4426.
 65. Lanier L, Le A, Civin C, Loken M, Phillips J. The relationship of CD16 (Leu-11) and Leu-19 (NKH-1) antigen expression on human peripheral blood NK cells and cytotoxic T lymphocytes. *J Immunol*. 1986;136:4480-4486.
 66. Knapp W, Rieber P, Dörken B, Schmidt R, Stein H, van dem Borne A. Towards a better definition of human leukocyte surface molecules. *Immunol Today*. 1989;10:253-258.
 67. Phillips J, Lanier L. Dissection of the lymphokine-activated killer phenomenon: relative contribution of peripheral blood natural killer cells and T lymphocytes to cytotoxicity. *J Exp Med*. 1986;164:814-825.
 68. National Committee for Clinical Laboratory Standards. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture (H3-A3). Villanova, PA: The National Committee for Clinical Laboratory Standards; 1991.
 69. National Committee for Clinical Laboratory Standards. Clinical Applications of Flow Cytometry: Quality Assurance and Immunophenotyping of Peripheral Blood Lymphocytes; Tentative Guideline. (H42-T). Villanova, PA: National Committee for Clinical Laboratory Standards; 1992.
 70. Giorgi J. Lymphocyte subset measurements: significance in clinical medicine. In: Rose N, Friedman H, Fahey J, eds. *Manual of Clinical Laboratory Immunology*. 3rd ed. Washington DC: American Society of Microbiology. 1986:236-246.
 71. Angadi C. Lack of Leu-3a epitope on T-helper (CD4) lymphocytes. *J Clin Lab Anal*. 1990;4:193-195.
 72. Reichert T, DeBruyère M, Deneys V, et al. Lymphocyte subset reference ranges in adult Caucasians. *Clin Immunol Immunopathol*. 1991;60:190-208.
 73. Prince H, Hirji K, Waldbeser L, Plaeger-Marshall S, Kleinman S, Lanier L. Influence of racial background on the distribution of T cell subsets and Leu-11-positive lymphocytes in healthy blood donors. *Diag Immunol*. 1985;3:33-37.
 74. Hicks C. *Fundamental Concepts in the Design of Experiments*. 3rd Edition. New York, NY: Holt, Rinehart, and Winston; 1982:222-223.

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