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1. Wolde-Mariam W, Peter J. Recent diagnostic advances in cellular immunology. *Diagnost Med.* 1984;7:25-32.
2. Henry C, Mishell B, Shiigi S. Hemolytic plaque assays. In: Mishell B, Shiigi S, ed. *Selected Methods in Cellular Immunology*. New York: WH Freeman and Co; 1980:69-123.
3. Grabstein K, Chen Y. Cell-mediated cytolytic responses. In: Mishell B, Shiigi S, ed. *Selected Methods in Cellular Immunology*. New York: WH Freeman and Co; 1980:124-137.
4. Wofsy L, Henry C, Kimura J, North J. Modification and use of antibodies to label cell surface antigens. In: Mishell B, Shiigi S, ed. *Selected Methods in Cellular Immunology*. New York: WH Freeman and Co; 1980:287-304.
5. Jackson A, Warner N. Preparation, staining, and analysis by flow cytometry of peripheral blood leukocytes. In: Rose N, Friedman H, Fahey J, ed. *Manual of Clinical Laboratory Immunology*. Washington, DC: American Society for Microbiology; 1986:226-235.
6. Jackson A. Basic phenotyping of lymphocytes: selection and testing of reagents and interpretation of data. *Clin Immunol Newsletter.* 1990;10:43-55.
7. Nicholson J, Jones B. Immunophenotyping by flow cytometry: its use in HIV infection. *Labmedica.* 1989;6:21-26.
8. Landay A, Muirhead K. Procedural guidelines for performing immunophenotyping by flow cytometry. *Clin Immunol Immunopathol.* 1989;52:48-60.
9. Schmidt R. Monoclonal antibodies for diagnosis of immunodeficiencies. *Blut.* 1989;59:200-206.
10. Ohno T, Kanoh T, Suzuki T, et al. Comparative analysis of lymphocyte phenotypes between carriers of human immunodeficiency virus (HIV) and adult patients with primary immunodeficiency using two-color immunofluorescence flow cytometry. *J Exp Med.* 1988;154:157.
11. Antel J, Bania M, Noronha A, Neely S. Defective suppressor cell function mediated by T8+ cell lines from patients with progressive multiple sclerosis. *J Immunol.* 1986;137:3436-3439.
12. Smolen J, Chused T, Leiserson W, Reeves J, Alling D, Steinberg A. Heterogeneity of immunoregulatory T-cell subsets in systemic lupus erythematosus: Correlation with clinical features. *Am J Med.* 1982;72:783-790.
13. Legendre C, Schiffrin A, Weitzner G, Colle E, Guttman R. Two color flow cytometry analysis of activated T-lymphocyte subsets in type 1 diabetes mellitus. *Diabetes.* 1988;37:792-795.
14. Tomkinson B, Wagner D, Nelson D, Sullivan J. Activated lymphocytes during acute Epstein-Barr virus infection. *J Immunol.* 1987;139:3802-3807.
15. Bishop G, Hall B, Duggin G, Horvath J, Sheil A, Tiller D. Immunopathology of renal allograft rejection analyzed with monoclonal antibodies to mononuclear cell markers. *Kidney Internat.* 1986;29:708-717.
16. Lanier L, Le A, Phillips J, Warner N, Babcock G. Subpopulations of human natural killer cells defined by expression of the Leu-7 (HNK-1) and Leu-11 (NK-15) antigens. *J Immunol.* 1983;131:1789-1796.
17. Giorgi J, Hultin L. Lymphocyte subset alterations and immunophenotyping by flow cytometry in HIV disease. *Clin Immunol Newsletter.* 1990;10:55-61.
18. Gratama J, Naipal A, Oljans P, et al. T lymphocyte repopulation and differentiation after bone marrow transplantation: Early shifts in the ratio between T4+ and T8+ T lymphocytes correlate with the occurrence of acute graft-versus-host disease. *Blood.* 1984;63:1416-1423.
19. Sides D, Casavant C, McHugh T, et al. Flow cytometric analysis of lymphocyte phenotypes in AIDS using monoclonal antibodies and simultaneous dual immunofluorescence. *Clin Immunol Immunopath.* 1986;38:161-177.
20. Lanier L, Le A, Civin C, Loken M, Phillips J. The relationship of CD16 (Leu-11) and Leu-19 (NKH-1) antigen expression on human peripheral blood NK cells and cytotoxic T lymphocytes. *J Immunol.* 1986;136:4480-4486.
21. Fitzgerald-Bocarsly P, Herberman R, Hercend T, et al. A definition of natural killer cells. In: Ates E, Lopez C, ed. *Natural Killer Cells and Host Defense*. Basel: Karger; 1989:1.
22. Aiuti F, Sirianni M, Mezzaroma I, et al. HIV-1 infection: epidemiological features and immunological alterations during the natural history of the disease. *Clin Immunol Immunopath.* 1989;50:S157-S165.
23. Jandl J. Blood and blood forming tissues. In: Jandl J, ed. *Blood Textbook of Hematology*. Boston, Mass: Little Brown and Company; 1987:1-48.
24. Beverley P. Production and use of monoclonal antibodies in transplantation immunology. In: Touraine J, Trager J, Betuel H, et al., ed. *Transplantation and Clinical Immunology XI*. Amsterdam: Excerpta Medica; 1980:87-94.
25. Cobbold S, Hale G, Waldmann H. Non-lineage, LFA-1 family, and leucocyte common antigens: new and previously defined clusters. In: McMichael A, ed. *Leucocyte Typing III*. New York, NY: Oxford University Press; 1986:788-803.
26. Rowe D, Beverley P. Characterisation of breast cancer infiltrates using monoclonal antibodies to human leucocyte antigens. *Br J Cancer.* 1984;49:149-159.
27. Bernard A, Boumsell L, Hill C. Joint report of the First International Workshop on Human Leucocyte Differentiation Antigens by the investigators of the participating laboratories: M2 protocol. In: Bernard A, Boumsell L, Dausset J, Milstein C, Schlossman S, ed. *Leucocyte Typing*. Berlin: Springer-Verlag; 1984:82-108.
28. Dimitriu-Bona A, Burmester G, Waters S, Winchester R. Human mononuclear phagocyte differentiation antigens. I. Patterns of antigenic expression on the surface of human monocytes and macrophages defined by monoclonal antibodies. *J Immunol.* 1983;130:145-152.
29. Herrmann F, Komischke B, Odenwald E, Ludwig W. Use of monoclonal antibodies as a diagnostic tool in human leukemia. I. Acute myeloid leukemia and acute phase of chronic myeloid leukemia. *Blut.* 1983;47:157-163.
30. Schwitzer R. Cluster report: CD45/CD45R. In: Knapp W, ed. *Leucocyte Typing IV*. New York, NY: Oxford University Press; 1989:628-634.
31. Goyert S, Tesio L, Ashman L, et al. Report on the CD14 Cluster Workshop. In: Knapp W, ed. *Leucocyte Typing IV*. New York, NY: Oxford University Press; 1989:789-794.
32. Bernstein I, Self S. Joint report of the myeloid section of the Second International Workshop on Human Leucocyte Differentiation Antigens. In: Reinherz E, Haynes B, Nadler L, Bernstein

- I, ed. *Leukocyte Typing II: Human Myeloid and Hematopoietic Cells*. New York: Springer-Verlag; 1984:1-25.
33. Jayaram Y, Hogg N. Surface expression of CD14 molecules on human neutrophils. In: Knapp W, Dorken B, Gilks W, et al., ed. *Leucocyte Typing IV*. New York, NY: Oxford University Press; 1989:796-797.
 34. Haynes B. Summary of T cell studies performed during the Second International Workshop and Conference on Human Leukocyte Differentiation Antigens. In: Reinherz E, Haynes B, Nadler L, Bernstein I, ed. *Leukocyte Typing II: Human T Lymphocytes*. New York: Springer-Verlag; 1986:3-30.
 35. Ledbetter J, Evans R, Lipinski M, Cunningham-Rundles C, Good R, Herzenberg L. Evolutionary conservation of surface molecules that distinguish T lymphocyte helper/inducer and T cytotoxic/suppressor subpopulations in mouse and man. *J Exp Med*. 1981;153:310-323.
 36. Kan E, Wang C, Wang L, Evans R. Non-covalently bonded subunits of 22 and 28 kd are rapidly internalized by T cells reacted with Anti-Leu-4 antibody. *J Immunol*. 1983;131:536-539.
 37. Knowles R. Immunochemical analysis of the T cell-specific antigens. In: Reinherz E, Haynes B, Nadler L, Bernstein I, ed. *Leukocyte Typing II: Human T Lymphocytes*. New York: Springer-Verlag; 1986:259-288.
 38. Nadler L. B Cell/Leukemia Panel Workshop: Summary and comments. In: Reinherz E, Haynes B, Nadler L, Bernstein I, ed. *Leukocyte Typing II: Human B Lymphocytes*. New York: Springer-Verlag; 1986:3-43.
 39. van Dongen J, Krissansen G, Wolvers-Tettero I, et al. Cytoplasmic expression of the CD3 antigen as a diagnostic marker for immature T-cell malignancies. *Blood*. 1988;71:603-612.
 40. Lanier L, Federspiel N, Ruitenberg J, et al. The T cell antigen receptor complex expressed on normal peripheral blood CD4-, CD8- T lymphocytes. *J Exp Med*. 1987;165:1076-1094.
 41. Dorken B, Moller P, Pezzutto A, Schwartz-Albiez R, Molderhauer G. B-cell antigens: CD19. In: Knapp W, ed. *Leucocyte Typing IV*. New York, NY: Oxford University Press; 1986:34-36.
 42. Bernard A, Boumsell L, Hill C. Joint report of the First International Workshop on Human Leukocyte Differentiation Antigens by the investigators of the participating laboratories: T2 protocol. In: Bernard A, Boumsell L, Dausett J, Milstein C, Schlossman S, ed. *Leucocyte Typing*. Berlin: Springer-Verlag; 1984:25-60.
 43. Evans R, Wall D, Platsoucas C, et al. Thymus-dependent membrane antigens in man: Inhibition of cell-mediated lympholysis by monoclonal antibodies to the TH2 antigen. *Proc Natl Acad Sci USA*. 1981;78:544-548.
 44. Maddon P, Dalgleish A, McDougall J, Clapham P, Weiss R, Axel R. The T4 gene encodes the AIDS virus receptor and is expressed in the immune system and the brain. *Cell*. 1986;47:333-348.
 45. Dalgleish A, Beverly P, Clapham P, Crawford D, Greaves M, Weiss R. The CD4 (T4) antigen is an essential component of the receptor for the AIDS virus. *Nature*. 1984;312:763-767.
 46. Kurlle R. Cluster report: CD3. In: Knapp W, ed. *Leucocyte Typing IV*. New York, NY: Oxford University Press; 1989:290-293.
 47. Wood G, Warner N, Warnke R. Anti-Leu-3/T4 antibodies react with cells of monocyte/macrophage and Langerhans lineage. *J Immunol*. 1983;131:212-216.
 48. Anderson P, Blue M-L, Morimoto C, Schlossman S. Cross-linking of T3 (CD3) with T4 (CD4) enhances the proliferation of resting T lymphocytes. *J Immunol*. 1987;139:678-682.
 49. Eichmann K, Johnson J-I, Falk I, Emmerich F. Effective activation of resting mouse T lymphocytes by cross-linking submitogenic concentrations of the T cell antigen receptor with either Lyt-2 or L3T4. *Eur J Immunol*. 1987;17:643-650.
 50. Gallagher P, Fazekas de St. Groth B, Miller J. CD4 and CD8 molecules can physically associate with the same T-cell receptor. *Proc Natl Acad Sci USA*. 1989;86:10044-10048.
 51. Rudd C, Burgess K, Barber E, Schlossman S. Monoclonal antibodies to the CD4 and CD8 antigens precipitate variable amounts of CD4/CD8-associated p56-lck activity. In: Knapp W, Dorken B, Gilks W, et al., ed. *Leucocyte Typing IV*. New York, NY: Oxford University Press; 1989:326-327.
 52. Terry L, Disanto J, Small T, Flomenberg N. Differential expression of the CD8 and Lyt-3 antigens on a subset of human T-cell receptor γ/δ -bearing lymphocytes. In: Knapp W, Dorken B, Gilks W, et al., ed. *Leucocyte Typing IV*. New York, NY: Oxford University Press; 1989:345-346.
 53. Moebius U. Cluster report: CD8. In: Knapp W, Dorken B, Gilks W, et al., ed. *Leucocyte Typing IV*. New York, NY: Oxford University Press; 1989:342-343.
 54. Lampson L, Levy R. Two populations of Ia-like molecules on a human B cell line. *J Immunol*. 1980;125:293-299.
 55. Brodsky F. A matrix approach to human Class II histocompatibility antigens: reactions of four monoclonal antibodies with the products of nine haplotypes. *Immunogenetics*. 1984;19:179-194.
 56. Robbins P, Evans E, Ding A, Warner N, Brodsky F. Monoclonal antibodies that distinguish between Class II antigens (HLA-DP, DQ, and DR) in 14 haplotypes. *Human Immunol*. 1987;18:301-313.
 57. Terstappen L, Hollander Z, Meiners H, Loken M. Quantitative comparison of myeloid antigens on five lineages of mature peripheral blood cells. *J Leuk Biol*. 1990;48:138-148.
 58. Warnke R, Levy R. Detection of T and B antigens with hybridoma monoclonal antibodies: a biotin-avidin-horseradish peroxidase method. *J Histochem Cytochem*. 1980;28:771-776.
 59. Engleman E, Warnke R, Fox R, Dille J, Benike C, Levy R. Studies of a human T lymphocyte antigen recognized by a monoclonal antibody. *Proc Natl Acad Sci USA*. 1981;78:1791-1795.
 60. Edwards J, Durant B, Jones D, Evans P, Smith J. Differential expression of HLA Class II antigens in fetal human spleen: relationship of HLA-DP, DQ, and DR to immunoglobulin expression. *J Immunol*. 1985;137:479-490.
 61. Perussia B, Starr S, Abraham S, Fanning V, Trinchieri G. Human natural killer cells analyzed by B73.1, a monoclonal antibody blocking Fc receptor functions. I. Characterization of the lymphocyte subset reactive with B73.1. *J Immunol*. 1983;130:2133-2141.
 62. Perussia B, Acuto O, Terhorst C, et al. Human natural killer cells analyzed by B73.1, a monoclonal antibody blocking Fc receptor functions. II. Studies of B73.1 antibody-antigen interaction on the lymphocyte membrane. *J Immunol*. 1983;130:2142-2148.
 63. Perussia B, Trinchieri G, Jackson A, et al. The Fc receptor for IgG on human natural killer cells: phenotypic, functional, and comparative studies with monoclonal antibodies. *J Immunol*. 1984;133:180-189.
 64. Schmidt R. Non-lineage/natural killer section report: new and previously defined clusters. In: Knapp W, Dorken B, Gilks W, et al., ed. *Leucocyte Typing IV*. New York, NY: Oxford University Press; 1989:517-542.
 65. Schubert J, Lanier L, Schmidt R. Cluster report: CD56. In: Knapp W, Dorken B, Gilks W, et al., ed. *Leucocyte Typing IV*. New York, NY: Oxford University Press; 1989:699-702.
 66. Knapp W, Rieber P, Dörken B, Schmidt R, Stein H, van der Borne A. Towards a better definition of human leucocyte surface molecules. *Immunol Today*. 1989;10:253-258.
 67. Lanier L, Testi R, Bindl J, Phillips J. Identity of Leu-19 (CD56) leukocyte differentiation antigen and neural cell adhesion molecule. *J Exp Med*. 1989;169:2233-2238.
 68. Slockbower J, Belgeri K, Bruck E, et al. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture (H3-A2). Villanova: The National Committee for Clinical Laboratory Standards, 1984.
 69. Reichert T, DéBruyere M, Deney V, et al. Lymphocyte subset reference ranges in adult Caucasians. *Clin Immunol Immunopath*. 1991;60:190-208.
 70. Prince H, Hirji K, Waldbeser L, Plaeger-Marshall S, et al. Influence of racial background on the distribution of T cell subsets and Leu-11-positive lymphocytes in healthy blood donors. *Diag Immunol*. 1985;3:33-37.

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