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- van Dongen JJM, Krissansen GW, Wolvers-Tettero ILM, et al. Cytoplasmic expression of the CD3 antigen as a diagnostic marker for immature T-cell malignancies. *Blood*. 1988;71:603-612.
- Brenner MB, Groh V, Porcelli SA, et al. Structure and distribution of the human  $\gamma/\delta$  T-cell receptor. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:1049-1053.
- Clevers H, Alarcón B, Wileman T, Terhorst C. The T cell receptor/CD3 complex: a dynamic protein ensemble. *Annu Rev Immunol*. 1988;6:629-662.
- Reichert T, DeBruyère M, Deneys V, et al. Lymphocyte subset reference ranges in adult Caucasians. *Clin Immunol Immunopath*. 1991;60:190-208.
- Lanier LL, Allison JP, Phillips JH. Correlation of cell surface antigen expression on human thymocytes by multi-color flow cytometric analysis: implications for differentiation. *J Immunol*. 1986;137:2501-2507.
- Garson JA, Beverley PCL, Coakham HB, Harper EJ. Monoclonal antibodies against human T lymphocytes label Purkinje neurones of many species. *Nature*. 1982;298:375-377.
- Kaneoka H, Perez-Rojas G, Sasasaki T, Benike CJ, Engleman EG. Human T lymphocyte proliferation induced by a pan-T monoclonal antibody (Anti-Leu-4): heterogeneity of response is a function of monocytes. *J Immunol*. 1983;131:158.
- Ledbetter JA, Evans RL, Lipinski M, Cunningham-Rundles C, Good RA, Herzenberg LA. Evolutionary conservation of surface molecules that distinguish T-lymphocyte helper/inducer and T cytotoxic/suppressor subpopulations in mouse and man. *J Exp Med*. 1981;153:310-323.
- Kan EAR, Wang CY, Wang LC, Evans RL. Noncovalently bonded subunits of 22 and 28 kd are rapidly internalized by T cells reacted with Anti-Leu-4 antibody. *J Immunol*. 1983;131:536-539.
- Knowles RW. Immunochemical analysis of the T-cell-specific antigens. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leucocyte Typing II: Human T Lymphocytes*. New York, NY: Springer-Verlag; 1986;1:259-288.
- Haynes BF. Summary of T-cell studies performed during the Second International Workshop and Conference on Human Leucocyte Differentiation Antigens. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leucocyte Typing II: Human T Lymphocytes*. New York, NY: Springer-Verlag; 1986;1:3-30.
- Allison JP, Lanier LL. Structure, function, and serology of the T-cell antigen receptor complex. *Ann Rev Immunol*. 1987;5:503.
- Freidrich W, O'Reilly RJ, Kozenger B, Gebhard DF, Good RA, Evans RL. T lymphocyte recipients of bone marrow transplants: analysis of regulatory T cell imbalances in GVHD. *Blood*. 1982;59:696.
- Kurrl R. Cluster report: CD3. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:290-293.
- Maecker HT, Maino VC. Analyzing T-cell responses to cytomegalovirus by cytokine flow cytometry. *Hum Immunol*. 2004;65:493-499.
- Maino VC, Picker LJ. Identification of functional subsets by flow cytometry: intracellular detection of cytokine expression. *Cytometry*. 1998;34:207-215.
- Campbell MJ, Scott J, Maecker HT, Park JW, Esserman LJ. Immune dysfunction and micrometastases in women with breast cancer. *Breast Cancer Res Treat*. 2005;91:163-171.
- Tu W, Chen S, Sharp M, et al. Persistent and selective deficiency of CD4+ T cell immunity to cytomegalovirus in immunocompetent young children. *J Immunol*. 2004;172:3260-3267.
- Wolde-Mariam W, Peter JB. Recent diagnostic advances in cellular immunology. *Diag Med*. 1984;7:25-32.
- Schmidt RE. Monoclonal antibodies for diagnosis of immunodeficiencies. *Blut*. 1989;59:200-206.
- Giorgi JV, Hultin LE. Lymphocyte subset alterations and immunophenotyping by flow cytometry

- in HIV disease. *Clin Immunol Newslett.* 1990;10:55-61.
22. Bernard A, Bousmell L, Hill C. Joint report of the first international workshop on human leucocyte differentiation antigens by the investigators of the participating laboratories. In: Bernard A, Bousmell L, Dausset J, Milstein C, Schlossman SF, eds. *Leucocyte Typing*. New York, NY: Springer-Verlag; 1984:9-108.
  23. Evans RL, Wall DW, Platsoucas CD, et al. Thymus-dependent membrane antigens in man: inhibition of cell-mediated lympholysis by monoclonal antibodies to the TH2 antigen. *Proc Natl Acad Sci USA.* 1981;78:544-548.
  24. Engleman EG, Benike CJ, Glickman E, Evans RL. Antibodies to membrane structures that distinguish suppressor/cytotoxic and helper T lymphocyte subpopulations block the mixed leukocyte reaction in man. *J Exp Med.* 1981;153:193-198.
  25. Kotzin BK, Benike CJ, Engleman EG. Induction of immunoglobulin secreting cells in the allogeneic mixed leukocyte reaction: regulation by helper and suppressor lymphocyte subsets in man. *J Immunol.* 1981;127:931-935.
  26. Wood GS, Warner NL, Warnke RA. Anti-Leu-3/T4 antibodies react with cells of monocyte/macrophage and Langerhans lineage. *J Immunol.* 1983;131:212-216.
  27. Dalgleish AG, Beverley PCL, Clapham PR, Crawford DH, Greaves MF, Weiss RA. The CD4 (T4) antigen is an essential component of the receptor for the AIDS retrovirus. *Nature.* 1984;312:763-767.
  28. Sattentau QJ, Dalgleish AG, Weiss RA, Beverley PCL. Epitopes of the CD4 antigen and HIV infection. *Science.* 1986;234:1120-1123.
  29. Lewis DE, Puck JM, Babcock GF, Rich R. Disproportionate expansion of a minor T cell subset in patients with lymphadenopathy syndrome and acquired immunodeficiency syndrome. *J Infect Dis.* 1985;151:555-559.
  30. Ohno T, Kanoh T, Suzuki T, et al. Comparative analysis of lymphocyte phenotypes between carriers of human immunodeficiency virus (HIV) and adult patients with primary immunodeficiency using two-color immunofluorescence flow cytometry. *J Exp Med.* 1988;154:157-172.
  31. Sites DP, Casavant CH, McHugh TM, et al. Flow cytometric analysis of lymphocyte phenotypes in AIDS using monoclonal antibodies and simultaneous dual immunofluorescence. *Clin Immunol Immunopathol.* 1986;38:161-177.
  32. Gatenby PA, Kansas GS, Xian CY, Evans RL, Engleman EG. Dissection of immunoregulatory subpopulations of T lymphocytes within the helper and suppressor sublineages in man. *J Immunol.* 1982;129:1997-2000.
  33. Picker LJ, Weiss LM, Medeiros LJ, Wood GS, Warnke RA. Immunophenotypic criteria for the diagnosis of non-Hodgkin's lymphoma. *Am J Pathol.* 1987;128:181-201.
  34. Maecker HT, Maino VC, Picker LJ. Immunofluorescence analysis of T-cell responses in health and disease. *J Clin Immunol.* 2000;20:391399.
  35. Knowles RW. Immunohistochemical analysis of the T-cell-specific antigens. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leucocyte Typing II: Human T Lymphocytes*. New York, NY: Springer-Verlag; 1986;1:259-288.
  36. Engleman EG, Warnke R, Fox RI, Dille J, Benike CJ, Levy R. Studies of a human T lymphocyte antigen recognized by a monoclonal antibody. *Proc Natl Acad Sci USA.* 1981;78:1791-1795.
  37. Ledbetter JA, Frankel AE, Herzenberg LA. Human Leu T-cell differentiation antigens: quantitative expression on normal lymphoid cells and cell lines. In: Hämmerling G, Hämmerling U, Kearney J, eds. *Monoclonal Antibodies and T Cell Hybridomas: Perspectives and Technical Notes*. New York, NY: Elsevier/North Holland; 1981:16-22.
  38. Warnke RA, Levy R. Detection of T and B antigens with hybridoma monoclonal antibodies: a biotinavidin-horseradish peroxidase method. *J Histochem Cytochem.* 1980;28:771-776.
  39. Warnke R, Miller R, Grogan T, Pederson M, Dille J, Levy R. Immunologic phenotype in 30 patients with diffuse large-cell lymphoma. *N Eng J Med.* 1980;303:293-300.
  40. Zipf RF, Fox R, Dille J, Levy R. Definition of the high risk ALL patient by immunologic phenotyping with monoclonal antibodies. *Cancer Res.* 1981;41:4786.
  41. Gadol N, Ault KA. Phenotypic and functional characterization of human Leu-1 (CD5) B cells. *Immunol Rev.* 1986;93:23.
  42. Royston I, Majda JA, Baird SM, Meserve BL, Griffiths JC. Human T-cell antigens defined by monoclonal antibodies: the 65,000-dalton antigen of T cells (T65) is also found on chronic lymphocytic leukemia cells bearing surface immunoglobulin. *J Immunol.* 1980;125:725.
  43. Antin JH, Emerson SG, Martin P, Gadol N, Ault KA. Leu-1+ (CD5) B cells. A major lymphoid subpopulation in human fetal spleen: Phenotypic and functional studies. *J Immunol.* 1986;136:505.
  44. Hardy RR, Hayakawa K, Shimizu M, Yamasaki K, Kishimoto T. Rheumatoid factor secretion from human Leu-1+ B cells. *Science.* 1987;236:81-83.
  45. Lanier LL, Le AM, Civin CI, Loken MR, Phillips JH. The relationship of CD16 (Leu-11) and Leu-19 (NKH-1) antigen expression on human peripheral blood NK cells and cytotoxic T lymphocytes. *J Immunol.* 1986;136:4480-4486.
  46. Lanier LL, Kipps T, Phillips JH. Functional properties of a unique subset of cytotoxic CD3+ T lymphocytes that express Fc receptors for IgG (CD16/Leu-11 antigen). *J Exp Med.* 1985;162:2089.
  47. Lanier LL, Le AM, Phillips JH, Warner NL, Babcock GF. Subpopulations of human natural killer cells defined by expression of the Leu-7 (HNK-1) and Leu-11 (NK-15) antigens. *J Immunol.* 1983;131:1789-1796.
  48. Perussia B, Trinchieri G, Jackson A, et al. The Fc receptor for IgG on human natural killer cells: phenotypic, functional, and comparative studies with monoclonal antibodies. *J Immunol.* 1984;133:180-189.
  49. Perussia B, Starr S, Abraham S, Fanning V, Trinchieri G. Human natural killer cells analyzed by B73.1, a monoclonal antibody blocking Fc receptor functions. I. Characterization of the lymphocyte subset reactive with B73.1. *J Immunol.* 1983;130:2133-2141.
  50. Gerosa F, Baldani-Guerra B, Nisii C, Marchesini V, Carra G, Trinchieri G. Reciprocal activating interaction between natural killer cells and dendritic cells. *J Exp Med.* 2002;195:327-333.
  51. Phillips JH, Babcock GF. A monoclonal antibody reactive against purified human natural killer cells and granulocytes. *Immunol Letters.* 1983;6:143.
  52. Thompson J, Goeken N, Brown S, Rhodes J. Phenotypic definition of human monocyte/macrophage subpopulations by monoclonal antibodies. II. Distinction of non-T, non-B natural killer cells (NK) from antigen-presenting stimulating cells in the mixed lymphocyte response (MLR). Presented at the American Association for Clinical Histocompatibility Testing;
  53. Perussia B, Acuto O, Terhorst C, et al. Human natural killer cells analyzed by B73.1, a monoclonal antibody blocking Fc receptor functions. II. Studies of B73.1 antibody-antigen interaction on the lymphocyte membrane. *J Immunol.* 1983;130:2142-2148.
  54. Schmidt RE. Non-lineage/natural killer section report: new and previously defined clusters. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:517-542.
  55. Cooper MD, Fehniger TA, Caligiuri MA. The biology of human natural killer-cell subsets. *Trends in Immunol.* 2001;22:633-640.
  56. Campbell JJ, Qin S, Unutmaz D, et al. Unique subpopulations of CD56+ NK and NK-T peripheral blood lymphocytes identified by chemokine receptor expression repertoire. *J Immunol.* 2001;166:6477-6482.
  57. Galandriani R, Tassi I, Mattia G. SH2-containing inositol phosphatase (SHIP-1) transiently translocates to raft domains and modulates CD16-mediated cytotoxicity in human NK cells. *Blood.* 2002;100:4581-4589.
  58. Fröland S, Natvig J, Bernald P. Surface-bound immunoglobulin as a marker of B lymphocytes in man. *Nature.* 1971;234:251-252.
  59. Warner N. Membrane immunoglobulins and antigen receptors on B and T lymphocytes. *Advances Immunol.* 1974;19:67.
  60. Kohler G, Milstein C. Continuous cultures of fused cells secreting antibody of predefined specificity. *Nature.* 1975;256:495.
  61. Dörken B, Möller P, Pezzutto A, Schwartz-Albiez R, Moldenhauer G. B-cell antigens: CD19. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. Oxford: Oxford University Press; 1989:34-36.
  62. Loken M, Shah V, Dattilio K, Civin C. Flow cytometric analysis of human bone marrow. II. Normal B-lymphocyte development. *Blood.* 1987;70:1316-1324.
  63. Ryan D, Kossover S, Mitchell S, Frantz C, Hennessy L, Cohen H. Subpopulations of common acute lymphoblastic leukemia antigen-positive lymphoid cells in normal bone marrow identified by hematopoietic differentiation antigens. *Blood.* 1986;68:417-425.
  64. Foucar K, Goeken J. Clinical application of immunologic techniques to the diagnosis of lymphoproliferative and immunodeficiency disorders. *Lab Medicine.* 1982;13:403-413.
  65. Meeker T, Miller R, Link M, Bindl J, Warnke R, Levy R. A unique human B lymphocyte antigen defined by a monoclonal antibody. *Hybridoma.* 1984;3:305-320.
  66. Clark P, Normansell D, Innes D, Hess C. Lymphocyte subsets in normal bone marrow. *Blood.* 1986;67:1600.
  67. Nadler L. B cell/leukemia panel workshop: Summary and comments. In: Reinherz E, Haynes B, Nadler L, ID B, eds. *Leucocyte Typing II: Human B Lymphocytes*. New York: Springer-Verlag; 1986:3-43.
  68. Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens on health-care settings. *MMWR.* 1988; 37(24).
  69. National Committee for Clinical Laboratory Standards. Protection of Laboratory Workers from Infectious Disease Transmitted by Blood, Body Fluids, and Tissue: Tentative Guideline (M29-T2). Villanova, PA: National Committee for Clinical Laboratory Standards; 1991.
  70. Centers for Disease Control. Guidelines for the performance of CD4+ T-cell determinations in persons with human immunodeficiency virus infection. *MMWR.* 1992; 41(No. RR-8).
  71. National Committee for Clinical Laboratory Standards. Clinical Applications of Flow Cytometry: Quality Assurance and Immunophenotyping of Peripheral Blood Lymphocytes (H42-T). Villanova, PA: National Committee for Clinical Laboratory Standards; 1992.
  72. Giorgi J. Lymphocyte subset measurements: Significance in clinical medicine. In: Rose N, Friedman H, Fahey J, eds. *Manual of Clinical Laboratory Standards*. Washington DC: American Society of Microbiology; 1986:236-246.
  73. Mishell B, Shiigi S, Henry C, et al. Preparation of mouse cell suspensions. In: Mishell B, Shiigi S, eds. *Selected Methods in Cellular Immunology*. New York: WH Freeman and Co; 1980:16-17.
  74. Jackson A, Warner N. Preparation, staining, and analysis by flow cytometry of peripheral blood leukocytes. In: Rose N, Friedman H, Fahey J, eds. *Manual of Clinical Laboratory Immunology*. 3rd ed. Washington, DC: American Society for Microbiology; 1986:226-235.
  75. Tursz T, Brouet J-C, Flandrin G, Danon F, Clauvel J-P, Seligmann M. Clinical and pathologic features of Waldenström's macroglobulinemia in seven patients with serum monoclonal IgG or IgA. *Am J Med.* 1977;63:499-502.
  76. Brown G, Greaves M. Enumeration of absolute numbers of T and B cells in human blood. *Scand J Immunol.* 1974;3:161.
  77. Dwyer J, Finklestein F, Mangi R, Fisher K, Hendler E. Assessment of the adequacy of immunosuppressive therapy using microscopy techniques to study immunologic competence. *Transplant Proc.* 1975;7:785.
  78. Dörken B, Möller P, Pezzutto A, Schwartz-Albiez R, Moldenhauer G. B-cell antigens: CD20. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. 1 ed. New York, NY: Oxford University Press; 1989:46-48.
  79. Hultin LE, Hausner MA, Hultin PM, Giorgi JV. CD20 (Pan-B cell) antigen is expressed at a low level on a subpopulation of human T lymphocytes. *Cytometry.* 1993;14:196-204.
  80. Ling NR, MacLennan ICM, Mason DY. B-cell and plasma cell antigens: new and previously defined clusters. In: McMichael AJ, ed. *Leucocyte Typing III: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1987:303-334.
  81. Clark EA, Shu G, Ledbetter JA. Role of the B p35 cell surface polypeptide in human B cell activation. *Proc Natl Acad Sci USA.* 1985;82:1766.
  82. Abts H, Emmerich M, Miltenyi S, Radbruch A, Tesch H. CD20 positive human B lymphocytes separated with the magnetic cell sorter (MACS) can be induced to proliferation and antibody secretion in vitro. *J Immunol Methods.* 1989;125:19-28.
  83. Marti GE, Faguet G, Bertin P, et al. CD20 and CD5 expression in B-chronic lymphocytic leukemia. *Ann NY Acad Sci.* 1992;651:480-483.
  84. Deans JP, Li H, Polyak MJ. CD20-mediated apoptosis: signalling through lipid rafts. *Immunology.* 2002;107:176-182.
  85. Lanier LL, Testi R, Bindl J, Phillips JH. Identity of the Leu-19 (CD56) leukocyte differentiation antigen and neural cell adhesion molecule. *J Exp Med.* 1989;169:2233-2238.
  86. Lanier LL, Chang C, Azuma M, Ruitenberg JJ, Hemperly JJ, Phillips JH. Molecular and functional analysis of human natural killer cell-associated neural cell adhesion molecule (N-CAM/CD56). *J Immunol.* 1991;146:4421-4426.
  87. Schubert J, Lanier LL, Schmidt RE. Cluster report: CD56. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York: Oxford University Press; 1989:699-702.
  88. Hercend T, Griffin J, Benussan A, et al. Generation of monoclonal antibodies to a human natural

killer clone. Characterization of two natural killer-associated antigens NKH1 and NKH2 expressed on subsets of large granular lymphocytes. *J Clin Invest.* 1985;75:932.

89. Edelman GM. Cell adhesion molecules. *Science.* 1983;219:450-457.

90. Phillips JH, Lanier LL. Dissection of the lymphokine-activated killer phenomenon: Relative contribution of peripheral blood natural killer cells and T lymphocytes to cytotoxicity. *J Exp Med.* 1986;164:814-825.

91. Abo T, Balch CMA. A differentiation antigen of human NK and K cells identified by a monoclonal antibody (HNK-1). *J Immunol.* 1981;127:1024-1029.

92. Manara DC, de Panfilis B, Ferrari C. Ultrastructural characterization of human large granular lymphocyte subsets defined by the expression of HNK-1 (Leu-7), Leu-11, or both HNK-1 and Leu-11 antigens. *J Histochem Cytochem.* 1985;33:1129-1133.

93. Mechtersheimer G, Staudter M, Möller P. Expression of the natural killer cell-associated antigens CD56 and CD57 in human neural and striated muscle cells and in their tumors. *Cancer Res.* 1991;51:1300-1307.

94. Kruse J, Mailhammer R, Wennecke H, et al. Neural cell-adhesion molecules and MAG share a common carbohydrate moiety recognized by monoclonal antibodies L2 and HNK-1. *Nature.* 1984;311:153-155.

95. McGarry RC, Helfand SL, Quarles RH, Roder JC. Recognition of myelin-associated glycoprotein by the monoclonal antibody HNK-1. *Nature.* 1983;306:376.

96. Bunn PA, Gazdar AF, Carney DN, Minna J. Small cell lung carcinoma and natural killer cells share an antigenic determinant, Leu-7. *Clin Res.* 1984;32:416A.

97. Schubert J, Lanier LL, Schmidt RE. Cluster report: CD57. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens.* Oxford: Oxford University Press; 1989:711-714.

98. Lewis DE, Puck JM, Babcock GM, Rich RR. Disproportionate expansion of a minor T cell subset in patients with lymphadenopathy syndrome and acquired immunodeficiency syndrome. *J Infect Disease.* 1985;151:555.

99. Prince HE, Kreiss JK, Kasper CK, et al. Distinctive lymphocyte subpopulation abnormalities in patients with congenital coagulation disorders who exhibit lymph node enlargement. *Blood.* 1985;66:64.

100. Stites D, Casavant C, McHugh T, et al. Flow cytometric analysis of lymphocyte phenotypes in AIDS using monoclonal antibodies and simultaneous dual immunofluorescence. *Clin Immunol Immunopathol.* 1986;38:161-177.

101. Wood GS, Warner N, Warnke R. Anti-Leu-3/T4 antibodies react with cells of monocyte/macrophage and Langerhans lineage. *J Immunol.* 1983;131(1):212-216.

102. Giorgi J, Hultin L. Lymphocyte subset alterations and immunophenotyping by flow cytometry in HIV disease. *Clin Immunol Newslett.* 1990;10(4):55-62.

103. Bernard A, Boumsell L, Hill C. Joint report of the First International Workshop on Human Leucocyte Differentiation Antigens by the investigators of the participating laboratories: T2 protocol. In: Bernard A, Boumsell L, Dausset J, Milstein C, Schlossman S, eds. *Leucocyte Typing.* Berlin: Springer-Verlag; 1984:25-60.

104. Picker LF, Weiss LM, Medeiros LJ, Wood GS, Warnke RA. Immunophenotypic criteria for the diagnosis of non-Hodgkin's lymphoma. *Am J Pathol.* 1987;128:181-201.

105. Schwinger R. Cluster report: CD45/CD45R. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens.* New York, NY: Oxford University Press; 1989:628-634.

106. Jackson A. Basic phenotyping of lymphocytes: selection and testing of reagents and interpretation of data. *Clin Immunol Newslett.* 1990;10:43-55.

107. Terstappen LWMM, Levin J. Bone marrow cell differential counts obtained by multidimensional flow cytometry. *Blood Cells.* 1992;18:311-330.

108. Loken MR, Brosnan J, Bach B, Ault K. Establishing optimal lymphocyte gates for immunophenotyping by flow cytometry. *Cytometry.* 1990;11:453-459.

109. Hermiston ML, Xu Z, Weiss A. CD45: a critical regulator of signaling thresholds in immune cells. *Annu Rev Immunol.* 2003;21:107-137.

110. Pellat-Deceunynck C, Bataille R. Normal and malignant human plasma cells: proliferation, differentiation, and expansions in relation to CD45 expression. *Blood Cells Mol Dis.* 2004;32:293-301.

111. Link M, Warnke R, Finlay J, et al. A single monoclonal antibody identifies T-cell lineage of childhood lymphoid malignancies. *Blood.* 1983;62:722-728.

112. Palker TJ, Searce RM, Hensley LL, Ho W, Haynes BF. Comparison of the CD7 (3A1) group of T cell workshop antibodies. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leucocyte Typing II: Human T Lymphocytes.* New York, NY: Springer-Verlag; 1986:303-313.

113. Weiss LM, Crabtree GS, Rouse RV, Warnke RA. Morphologic and immunologic characterization of 50 peripheral T-cell lymphomas. *Am J Pathol.* 1985;118:316-324.

114. Weiss LM, Wood GS, Warnke RA. Immunophenotypic differences between dermatopathic lymphadenopathy and lymph node involvement in mycosis fungoides. *Am J Pathol.* 1985;120:179-185.

115. Foon KA, Todd RF. Immunologic classification of leukemia and lymphoma. *Blood.* 1986;68:1-31.

116. Wood GS, Abel EA, Hoppe RT, Warnke RA. Leu-8 and Leu-9 antigen phenotypes: immunologic criteria for the distinction of mycosis fungoides from cutaneous inflammation. *J Am Acad Dermatol.* 1986;14:1006-1013.

117. Lazarovits AI, White MJ, Karsh J. CD7- T cells in rheumatoid arthritis. *Arthritis Rheum.* 1992;35:615-624.

118. Lazarovits AI, Karsh J. Decreased expression of CD7 occurs in rheumatoid arthritis. *Clin Exp Immunol.* 1988;72:470-475.

119. Greaves MF. Monoclonal antibodies as probes for leukemic heterogeneity and hematopoietic differentiation. In: Knapp W, ed. *Leukemia Markers.* New York, NY: Academic Press; 1981:19.

120. Zola H. CD10 Workshop Panel report. In: Schlossman SF, Boumsell L, Gilks W, et al, eds. *Leucocyte Typing V: White Cell Differentiation Antigens.* New York, NY: Oxford University Press; 1995:505-507.

121. Letarte M, Vera S, Tran R, et al. Common acute lymphocytic leukemia antigen is identical to neutral endopeptidase. *J Exp Med.* 1988;168:1247-1253.

122. Le Bien TW, McCormack RT. The common acute lymphoblastic leukemia antigen (CD10): emancipation from a functional enigma. *Blood.* 1989;73:625-635.

123. Consolini R, Legitimo A, Rondelli R, et al. Clinical relevance of CD10 expression in childhood ALL. *Haematologica.* 1998;83:967-973.

124. Nadler LM, Korsmeyer SJ, Anderson KC, et al. B cell origin of non-T cell acute lymphoblastic

leukemia: a model for discrete stages of neoplastic and normal pre-B cell differentiation. *J Clin Invest.* 1984;74:332-340.

125. Hann IM, Richards SM, Eden OB, Hill FGH. Analysis of the immunophenotype of children treated on the Medical Research Council United Kingdom Acute Lymphoblastic Leukaemia Trial XI (MRC UKALLXI). *Leukemia.* 1998;12:1249-1255.

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