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1. Schmidt RE. Monoclonal antibodies for diagnosis of immunodeficiencies. *Blut*. 1989;59:200-206.
2. Nicholson JKA. Use of flow cytometry in the evaluation and diagnosis of primary and secondary immunodeficiency diseases. *Arch Pathol Lab Med*. 1989;113:598-605.
3. Foucar K, Goeken JA. Clinical application of immunologic techniques to the diagnosis of lymphoproliferative and immunodeficiency disorders. *Lab Med*. 1982;13:403-413.
4. Cohen SB, Weetman AP. Activated interstitial and intraepithelial thyroid lymphocytes in autoimmune thyroid disease. *Acta Endocrinol*. 1988;119:161-166.
5. Smolen JS, Chused TM, Leiserson WM, Reeves JP, Alling D, Steinberg AD. Heterogeneity of immunoregulatory T cell subsets in systemic lupus erythematosus: correlation with clinical features. *Am J Med*. 1982;72:783-790.
6. Giorgi JV, Hultin LE. Lymphocyte subset alterations and immunophenotyping by flow cytometry in HIV disease. *Clin Immunol Newslett*. 1990;10:55-61.
7. Landay A, Ohlsson-Wilhelm B, Giorgi JV. Application of flow cytometry to the study of HIV infection. *AIDS*. 1990;4:479-497.
8. Antel J, Bania M, Noronha A, Neely S. Defective suppressor cell function mediated by T8+ cell lines from patients with progressive multiple sclerosis. *J Immunol*. 1986;137:3436-3439.
9. Gratama JW, Naipal A, Olijans P, et al. T-lymphocyte repopulation and differentiation after bone marrow transplantation: early shifts in the ratio between T4+ and T8+ T lymphocytes correlate with the occurrence of acute graft-versus-host disease. *Blood*. 1984;63:1416-1423.
10. Bishop GA, Hall BM, Duggin GG, Horvath JS, Sheil AGR, Tiller DJ. Immunopathology of renal allograft rejection analyzed with monoclonal antibodies to mononuclear cell markers. *Kidney Internat*. 1986;29:708-717.
11. Wolde-Mariam W, Peter JB. Recent diagnostic advances in cellular immunology. *Diag Med*. 1984;7:25-32.
12. Centers for Disease Control. 1997 Revised guidelines for performing CD4+ T-cell determinations in persons with human immunodeficiency virus (HIV). *MMWR*. 1997;46(No. RR-2):1-29.
13. Nicholson JKA, Jones BM, Hubbard M. CD4 T lymphocyte determinations on whole blood specimens using a single-tube three-color assay. *Cytometry*. 1993;14:685-689.
14. Nicholson J, Kidd P, Mandy F, Livnat D, Kagan J. Threecolor supplement to the NIAID DAIDS guideline for flow cytometric immunophenotyping. *Cytometry*. 1996;26:227-230.
15. Nicholson JKA, Hubbard M, Jones BM. Use of CD45 fluorescence and side-scatter characteristics for gating lymphocytes when using the whole blood lysis procedure and flow cytometry. *Cytometry*. 1996;26: 16-21.
16. Evans RL, Wall DW, Platsoucas CD, et al. Thymusdependent membrane antigens in man: inhibition of cellmediated lympholysis by monoclonal antibodies to the TH2 antigen. *Proc Natl Acad Sci USA*. 1981;78:544-548.
17. Bernard A, Bounsell L, Hill C. Joint report of the first international workshop on human leucocyte differentiation antigens by the investigators of the participating laboratories: T2 protocol. In: Bernard A, Bounsell L, Dausset J, Milstein C, Schlossman SF, eds. *Leucocyte Typing*. New York, NY: Springer-Verlag; 1984:25-60. 18. Wood GS, Warner NL, Warnke RA. Anti-Leu-3/T4 antibodies react with cells of monocyte/macrophage and Langerhans lineage. *J Immunol*. 1983;131:212-216.
18. Haynes BF. Summary of T-cell studies performed during the Second International Workshop and Conference on Human Leucocyte Differentiation Antigens. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leucocyte Typing II: Human T Lymphocytes*. New York, NY: Springer-Verlag; 1986:3-30.
19. Kan EAR, Wang CY, Wang LC, Evans RL. Noncovalently bonded subunits of 22 and 28 kd are rapidly internalized by T cells reacted with Anti-Leu-4 antibody. *J Immunol*. 1983;131:536-539.
20. Knowles RW. Immunohistochemical analysis of the T-cell-specific antigens. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leucocyte Typing II: Human T Lymphocytes*. New York, NY: Springer-Verlag; 1986;1:259-288.
21. Organizing Committee of the Fourth International Workshop on Human Leucocyte Differentiation Antigens. Appendix A: CD guide. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press Inc; 1989:1074-1093.
22. Dalgleish AG, Beverley PCL, Clapham PR, Crawford DH, Greaves MF, Weiss RA. The CD4 (T4) antigen is an essential component of the receptor for the AIDS retrovirus. *Nature*. 1984;312:763-767.
23. Maddon PJ, Dalgleish AG, McDougal JS, Clapham PR, Weiss RA, Axel R. The T4 gene encodes the AIDS virus receptor and is expressed in the immune system and the brain. *Cell*. 1986;47:333-348.
24. Rudd CE, Burgess KE, Barber EK, Schlossman SF. Monoclonal antibodies to the CD4 and CD8 antigens precipitate variable amounts of CD4/CD8-associated p56lck activity. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:326-327.
25. Schwinger R. Cluster report: CD45/CD45R. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:628-634.
26. Moebius U. Cluster report: CD8. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:342-343.
27. Gallagher PF, Fazekas de St. Groth B, Miller JFAP. CD4 and CD8 molecules can physically associate with the same T-cell receptor. *Proc Natl Acad Sci, USA*. 1989;86:10044-10048.
28. Anderson P, Blue M-L, Morimoto C, Schlossman SF. Cross-linking of T3 (CD3) with T4 (CD4) enhances the proliferation of resting T lymphocytes. *J Immunol*. 1987;139:678-682.
29. Eichmann K, Jönsson J-I, Falk I, Emrlich F. Effective activation of resting mouse T lymphocytes by cross-linking submitogenic concentrations of the T cell antigen receptor with either Lyr-2 or L3T4. *Eur J Immunol*. 1987;17:643-650.
30. Wood GS, Warner NL, Warnke RA. Anti-Leu-3/T4 antibodies react with cells of monocyte/macrophage and Langerhans lineage. *J Immunol*. 1983;131:212-216.
31. Moebius U. Cluster report: CD8. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:342-343.
32. Schwinger R. Cluster report: CD45/CD45R. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:628-634.
33. Cobbold SP, Hale G, Waldmann H. Non-lineage, LFA-1 family, and leucocyte common antigens: new and previously defined clusters. In: McMichael AJ, ed. *Leucocyte Typing III: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1987:788-803.
34. van Dongen JJM, Krissansen GW, Wolvers-Tettero ILM, et al. Cytoplasmic expression of the CD3 antigen as a diagnostic marker for immature T-cell malignancies. *Blood*. 1988;71:603-612.
35. Brenner MB, McClean J, Dialynas DP, et al. Identification of a putative second T cell receptor. *Nature*. 1986;322:145-149.
36. Clevers H, Alarcón B, Wileman T, Terhorst C. The T cell receptor/CD3 complex: a dynamic protein ensemble. *Annu Rev Immunol*. 1988;6:629-662.
37. Fitzgerald-Bocarsly P, Herberman R, Hercend T, et al. A definition of natural killer cells. In: Ades E, Lopez C, eds. *Natural Killer Cells and Host Defense*. Basel: Karger; 1989:1.
38. Lanier LL, Le AM, Civin CI, Loken MR, Phillips JH. The relationship of CD16 (Leu-11) and Leu-19 (NKH-1) antigen expression on human peripheral blood NK cells and cytotoxic T lymphocytes. *J Immunol*. 1986;136:4480-4486.
39. Perussia B, Acuto O, Terhorst C, et al. Human natural killer cells analyzed by B73.1, a monoclonal antibody blocking Fc receptor functions, II: studies of B73.1 antibody-antigen

- interaction on the lymphocyte membrane. *J Immunol.* 1983;130:2142-2148.
40. Perussia B, Starr S, Abraham S, Fanning V, Trinchieri G. Human natural killer cells analyzed by B73.1, a monoclonal antibody blocking Fc receptor functions, I: characterization of the lymphocyte subset reactive with B73.1. *J Immunol.* 1983;130:2133-2141.
 41. Schmidt RE. Non-lineage/natural killer section report: new and previously defined clusters. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens.* New York, NY: Oxford University Press; 1989:517-542.
 42. Ritz J, Trinchieri G, Lanier LL. NK-cell antigens: section report. In: Schlossman SF, Bousmell L, Gilks W, et al, eds. *Leucocyte Typing V: White Cell Differentiation Antigens.* New York, NY: Oxford University Press; 1995;2:1367-1372.
 43. Nadler LM. B Cell/Leukemia Panel Workshop: summary and comments. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leucocyte Typing II: Human B Lymphocytes.* New York, NY: Springer-Verlag; 1986;2:3-43.
 44. Dörken B, Möller P, Pezzutto A, Schwartz-Albiez R, Moldenhauer G. B-cell antigens: CD19. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens.* New York, NY: Oxford University Press; 1989:34-36.
 45. Berti E, Parravicini C, Cattoretti G, Delia D, de Braud F, Cusini M. Immunohistochemical reactivity of anti-B cell monoclonal antibodies in thymus, lymph node, and normal skin. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leucocyte Typing II: Human B Lymphocytes.* New York, NY: Springer-Verlag; 1986;2:313-318.
 46. Centers for Disease Control. Update: universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens in health-care settings. *MMWR.* 1988;37:377-388.
 47. Protection of Laboratory Workers from Infectious Disease Transmitted by Blood, Body Fluids, and Tissue: Tentative Guideline. Villanova, PA: National Committee for Clinical Laboratory Standards; 1991. NCCLS document M29-T2.
 48. Nicholson JK, Browning SW, Orloff SL, McDougal JS. Inactivation of HIV-infected H9 cells in whole blood preparations by lysing/fixing reagents used in flow cytometry. *J Immunol Methods.* 1993;160:215-218.
 49. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture: Approved Standard. Villanova, PA: National Committee for Clinical Laboratory Standards; 1991. NCCLS document H3-A3.
 50. Clinical Applications of Flow Cytometry: Quality Assurance and Immunophenotyping of Peripheral Blood Lymphocytes; Tentative Guideline. Villanova, PA: National Committee for Clinical Laboratory Standards; 1992. NCCLS document H42-T.
 51. Giorgi JV. Lymphocyte subset measurements: significance in clinical medicine. In: Rose NR, Friedman H, Fahey JL, eds. *Manual of Clinical Laboratory Immunology.* 3rd ed. Washington, DC: American Society for Microbiology; 1986:236-246.
 52. Jackson AL, Warner NL. Preparation, staining, and analysis by flow cytometry of peripheral blood leukocytes. In: Rose NR, Friedman H, Fahey JL, eds. *Manual of Clinical Laboratory Immunology.* 3rd ed. Washington, DC: American Society for Microbiology; 1986:226-235.
 53. Prince HE, Hirji K, Waldbeser LS, Plaeger-Marshall S, Kleinman S, Lanier LL. Influence of racial background on the distribution of T-cell subsets and Leu 11-positive lymphocytes in healthy blood donors. *Diagn Immunol.* 1985;3(1):33-37.
 54. Angadi CV. Lack of Leu-3a epitope on T-helper (CD4) lymphocytes. *J Clin Lab Anal.* 1990;4:193-195.
 55. How to Define and Determine Reference Intervals in the Clinical Laboratory: Approved Guideline. Villanova, PA: National Committee for Clinical Laboratory Standards; 1992. NCCLS document C28-A.
 56. Lanier LL, Le AM, Phillips JH, Warner NL, Babcock GF. Subpopulations of human natural killer cells defined by expression of the Leu-7 (HNK-1) and Leu-11 (NK-15) antigens. *J Immunol.* 1983;131:1789-1796.

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