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1. Terstappen LWMM, Huang S, Picker LJ. Flow cytometric assessment of human T cell differentiation in thymus and bone marrow. Blood.1992;79:666-677.
2. Porwir-MacDonald, A. Multiparameter flow cytometry detection of minimal residual disease in T-cell acute lymphoblastic leukemia. Presented at the International Symposium on Minimal Residual Disease: From Methodological Problems to Clinical Goals; October 31-November 1,1997; Salamanca, Spain.
3. Knowles DM. Immunophenotypic and antigen receptor gene rearrangement analysis in T cell neoplasia. Am J Pathol. 1989;134:761-785.
4. Ginaldi L, Matutes E, Farahat N, De Martinis M, Morilla R, Catovsky D. Differential expression of CD3 and CD7 in T-cell malignancies: a quantitative study by flow cytometry. Br J Haematol. 1996;93:921-927.
5. Hastrup N, Ralfkiaer E, Pallesen G. Aberrant phenotypes in peripheral T cell lymphomas. J Clin Pathol. 1989;42:398-402.
6. Ferrara F, Cimino R, Antinolfi I, et al. Clinical relevance of immunological dissection in TALL: a report on 20 cases with stem cell (CD7+, CD4-, CD8-, CD1-) phenotype. Am J Hematol. 1992;40:98-102.
7. Khalidi HS, Medeiros LJ, Chang KL, Brynes RK, Slovak ML, Arber DA. The immunophenotype of adult acute myeloid leukemia: high frequency of lymphoid antigen expression and comparison of immunophenotype, French-American-British classification, and karyotypic abnormalities. Am J Clin Pathol. 1998;109:211-220.
8. Lauder TM, Bray RA, Stempora L, Chenggis ML, Farhi DC. Lymphoid-associated antigen expression by acute myeloid leukemia. Am J Clin Pathol. 1996;106:185-191.
9. Macedo A, Orfao A, Vidriales MB, et al. Characterization of aberrant phenotypes in acute myeloblastic leukemia. Ann Hematol. 1995;70:189-194.
10. Bieber CP, Howard FD, Pennock J, Wong J, Shorthouse R, Stinson EB. Preparation, characterization, and primate testing of monoclonal antithymocyte globulin. Transplantation. 1981;31:283-289.
11. Link M, Warnke R, Finlay J, et al. A single monoclonal antibody identifies T-cell lineage of

- childhood lymphoid malignancies. *Blood*. 1983;62:722-728.
12. Palker TJ, Scearce RM, Hensley LL, Ho W, Haynes BF. Comparison of the CD7 (3A1) group of T cell workshop antibodies. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human T Lymphocytes*. New York, NY: Springer-Verlag; 1986:303-313.
 13. van Dongen JJM, Krissansen GW, Wolvers-Tettero ILM, et al. Cytoplasmic expression of the CD3 antigen as a diagnostic marker for immature T-cell malignancies. *Blood*. 1988;71:603-612.
 14. Clevers H, Alarcón B, Wileman T, Terhorst C. The T cell receptor/CD3 complex: a dynamic protein ensemble. *Annual Rev Immunol*. 1988;6:629-662.
 15. Howard FD, Ledbetter JA, Wong J, Bieber CP, Stinson EB, Herzenberg LA. A human T lymphocyte differentiation marker defined by monoclonal antibodies that block E-rosette formation. *J Immunol*. 1981;126:2117-2122.
 16. Lanier LL, Phillips JH. A map of the cell surface antigens expressed on resting and activated human natural killer cells. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human Myeloid and Hematopoietic Cells*. New York, NY: Springer-Verlag; 1986:3:157-170.
 17. Weiss LM, Crabtree GS, Rouse RV, Warnke RA. Morphologic and immunologic characterization of 50 peripheral T-cell lymphomas. *Am J Pathol*. 1985;118:316-324.
 18. Reiter C. Cluster report: CD7. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leukocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:341-342.
 19. Ledbetter JA, Evans RL, Lipinski M, Cunningham-Rundles C, Good RA, Herzenberg LA. Evolutionary conservation of surface molecules that distinguish T lymphocyte helper/inducer and T cytotoxic/suppressor subpopulations in mouse and man. *J Exp Med*. 1981;153:310-323.
 20. Haynes BF. Summary of T-cell studies performed during the Second International Workshop and Conference on Human Leukocyte Differentiation Antigens. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human T Lymphocytes*. New York, NY: Springer-Verlag; 1986:3-30.
 21. Kan E, Wang C, Wang L, Evans R. Non-covalently bonded subunits of 22 and 28 kd are rapidly internalized by T cells reacted with Anti-Leu-4 antibody. *J Immunol*. 1983;131:536-539.
 22. Knowles RW. Immunohistochemical analysis of the T cell-specific antigens. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human T Lymphocytes*. New York, NY: Springer-Verlag; 1986:259-288.
 23. Hoffkes HG, Schmidtke G, Uppenkamp M, Schmucker U. Multiparametric immunophenotyping of B cells in peripheral blood of healthy adults by flow cytometry. *Clin Diagn Lab Immunol*. 1996;3:30-36.
 24. Loken MR, Shah VO, Dattilio KL, Civin CI. Flow cytometric analysis of human bone marrow. II. Normal B-lymphocyte development. *Blood*. 1987;70:1316-1324.
 25. Freedman AS, Nadler LM. Immunologic markers in non-Hodgkin's lymphoma. *Hematol Oncol Clin North Am*. 1991;5(5):871-889.
 26. Jennings CD, Foon KA. Recent advances in flow cytometry: application to the diagnosis of hematologic malignancy. *Blood*. 1997;90:2863-2892.
 27. Zola H. CD10 Workshop Panel report. In: Schlossman SF, Bousmell L, Gilks W, et al, eds. *Leukocyte Typing V: White Cell Differentiation Antigens*. New York: Oxford University Press; 1995:505-507.
 28. Gadol N, Ault K. Phenotypic and functional characterization of human Leu-1 (CD5) B cells. *Immunol Rev*. 1986;93:23.
 29. Schmidt CJ, Domenico L, Ward P, Barcos MP, Stewart CC. Aberrant antigen expression detected by multiparameter three color flow cytometry in intermediate and high grade B-cell lymphomas [In Process Citation]. *Leuk Lymphoma*. 1999;34(5-6):539-544.
 30. Cabezudo E, Carrara P, Morilla R, Matutes E. Quantitative analysis of CD79b, CD5 and CD19 in mature B-cell lymphoproliferative disorders. *Haematologica*. 1999;84(5):413-418.
 31. Almasri NM, Iturraspe JA, Braylan RC. CD10 expression in follicular lymphoma and large cell lymphoma is different from that of reactive lymph node follicles. *Arch Pathol Lab Med*. 1998;122:539-544.
 32. Pui CH, Rivera GK, Hancock ML, et al. Clinical significance of CD10 expression in childhood acute lymphoblastic leukemia. *Leukemia*. 1993;7(1):35-40.
 33. De Rossi G, Grossi C, Foa R, et al. Immunophenotype of acute lymphoblastic leukemia cells: the experience of the Italian Cooperative Group (Gimema). *Leuk Lymphoma*. 1993;9(3):221-228.
 34. Geisler CH, Larsen JK, Hansen NE, et al. Prognostic importance of flow cytometric immunophenotyping of 540 consecutive patients with B-cell chronic lymphocytic leukemia. *Blood*. 1991;78:1795-1802.
 35. DiGiuseppe JA, Borowitz MJ. Clinical utility of flow cytometry in the chronic lymphoid leukemias. *Semin Oncol*. 1998;25(1):6-10.
 36. Baldini L, Cro L, Cortezzi A, et al. Immunophenotypes in "classical" B-cell chronic lymphocytic leukemia. Correlation with normal cellular counterpart and clinical findings. *Cancer*. 1990;66(8):1738-1742.
 37. Spier CM, Grogan TM, Fielder K, Richter L, Rangel C. Immunophenotypes in "well-differentiated" lymphoproliferative disorders, with emphasis on small lymphocytic lymphoma. *Hum Pathol*. 1986;17(11):1126-1136.
 38. Greaves M. Monoclonal antibodies as probes for leukemic heterogeneity and haematopoietic differentiation. In: Knapp W, ed. *Leukemia Markers*. New York: Academic Press; 1981:19.
 39. Letarte M, Vera S, Tran R, et al. Common acute lymphocytic leukemia antigen is identical to neutral endopeptidase. *J Exp Med*. 1988;168:1247-1253.
 40. Moldenhauer G, Dörken B, Schwartz R, Pezzutto A, Knops J, Hämmerling GJ. Analysis of ten B lymphocyte-specific workshop monoclonal antibodies. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human B Lymphocytes*. New York: Springer-Verlag; 1986:2:61-67.
 41. Nadler LM. B Cell/Leukemia Panel Workshop: summary and comments. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human B Lymphocytes*. New York: Springer-Verlag; 1986:2:3-43.
 42. Engleman EG, Warnke R, Fox RI, Dilley J, Benike CJ, Levy R. Studies of a human T lymphocyte antigen recognized by a monoclonal antibody. *Proc Natl Acad Sci USA*. 1981;78:1791-1795.
 43. Ledbetter JA, Frankel AE, Herzenberg LA. Human Leu T-cell differentiation antigens: Quantitative expression on normal lymphoid cells and cell lines. In: Hämmerling G, Hämmerling U, Kearney J, eds. *Monoclonal Antibodies and T Cell Hybridomas: Perspectives and Technical Notes*. New York: Elsevier/North Holland; 1981:16-22.
 44. Warnke RA, Levy R. Detection of T and B antigens with hybridoma monoclonal antibodies: a biotin-avidin-horseradish peroxidase method. *J Histochem Cytochem*. 1980;28:771-776.
 45. Warnke R, Miller R, Grogan T, Pederson M, Dilley J, Levy R. Immunologic phenotype in 30 patients with diffuse large-cell lymphoma. *N Engl J Med*. 1980;303:293-300.
 46. Zipf RF, Fox R, Dilley J, Levy R. Definition of the high risk ALL patient by immunologic phenotyping with monoclonal antibodies. *Cancer Res*. 1981;41:4786.
 47. Royston I, Majda JA, Baird SM, Meserve BL, Griffiths JC. Human T-cell antigens defined by monoclonal antibodies: the 65,000-dalton antigen of T cells (T65) is also found on chronic lymphocytic leukemia cells bearing surface immunoglobulin. *J Immunol*. 1980;125:725.
 48. Kilo MN, Dorfman DM. The utility of flow cytometric immunophenotypic analysis in the distinction of small lymphocytic lymphoma/chronic lymphocytic leukemia from mantle cell lymphoma. *Am J Clin Pathol*. 1996;105:451-457.
 49. LeBien T, McCormack R. The common acute lymphoblastic leukemia antigen (CD10)-emancipation from a functional enigma. *Blood*. 1989;73:625-635.
 50. Reichert T, DeBruyère M, Deneys V, et al. Lymphocyte subset reference ranges in adult Caucasians. *Clin Immunol Immunopath*. 1991;60:190-208.
 51. Tedder T, Zhou L-J, Engel P. The CD19/CD21 signal transduction complex of B lymphocytes. *Immunol Today*. 1994;15(9):437-442.
 52. Weinberg DS, Pinkus GS, Ault KA. Cytofluorometric detection of B cell clonal excess: a new approach to the diagnosis of B-cell lymphoma. *Blood*. 1984;63:1080-1087.
 53. Tétu B, Manning JT, Ordóñez NG. Comparison of monoclonal and polyclonal antibodies directed against immunoglobulin light and heavy chains in non-Hodgkin's lymphoma. *Am J Clin Pathol*. 1985;85:25-31.
 54. Molot RJ, Meeker TC, Wittwer CT, et al. Antigen expression and polymerase chain reaction amplification of mantle cell lymphomas. *Blood*. 1994;83:1626-1631.
 55. Almasri NM, Duque RE, Iturraspe J, Everett E, Braylan RC. Reduced expression of CD20 antigen as a characteristic marker for chronic lymphocytic leukemia. *Am J Hematol*. 1992;40:259-263.
 56. Zander DS, Iturraspe JA, Everett ET, Massey JK, Braylan RC. Flow cytometry. In vitro assessment of its potential application for diagnosis and classification of lymphoid processes in cytologic preparations from fine-needle aspirates. *Am J Clin Pathol*. 1994;101:577-586.
 57. Matutes E, Morilla R, Owusu-Ankomah K, Houlihan A, Meeus P, Catovsky D. The immunophenotype of hairy cell leukemia (HCL). Proposal for a scoring system to distinguish HCL from B-cell disorders with hairy or villous lymphocytes. *Leuk Lymphoma*. 1994;14 Suppl 1:57-61.
 58. Tworek JA, Singleton TP, Schnitzer B, Hsi ED, Ross CW. Flow cytometric and immunohistochemical analysis of small lymphocytic lymphoma, mantle cell lymphoma, and plasmacytoid small lymphocytic lymphoma. *Am J Clin Pathol*. 1998;110:582-589.
 59. Huh YO, Pugh WC, Kantarjian HM, et al. Detection of subgroups of chronic B-cell leukemias by FMC7 monoclonal antibody. *Am J Clin Pathol*. 1994;101:283-289.
 60. Hubl W, Iturraspe J, Braylan RC. FMC7 antigen expression on normal and malignant B-cells can be predicted by expression of CD20. *Cytometry*. 1998;34:71-74.
 61. Marti GE, Zenger V, Caproaso NE, et al. Antigenic expression of B-cell chronic lymphocytic leukemic lymphocytes. *Anal Quant Cytol Histol*. 1989;11:315-323.
 62. Hassan IB, Hagberg H, Sundstrom C. Immunophenotype of hairy-cell leukemia. *Eur J Haematol*. 1990;45:172-176.
 63. Ault KA. Flow cytometric evaluation of normal and neoplastic B cells. In: Rose NR, Friedman H, Fahey JL, eds. *Manual of Clinical Laboratory Immunology*. Washington, DC: American Society for Microbiology; 1986:247-253.
 64. Harris NL, Data RE. The distribution of neoplastic and normal B-lymphoid cells in nodular lymphomas: use of an immunoperoxidase technique on frozen sections. *Hum Pathol*. 1982;13:610-617.
 65. Meis JM, Osborne BM, Butler JJ. A comparative marker study of large cell lymphoma, Hodgkin's disease, and true histiocytic lymphoma in paraffin-embedded tissue. *Am J Clin Pathol*. 1986;86:591-599.
 66. Picker LJ, Weiss LM, Medeiros LJ, Wood GS, Warnke RA. Immunophenotypic criteria for the diagnosis of non-Hodgkin's lymphoma. *Am J Pathol*. 1987;128:181-201.
 67. Foon KA, Todd RF. Immunologic classification of leukemia and lymphoma. *Blood*. 1986;68:1-31.
 68. Tubbs RR, Sheibani K, Weiss RA, Sebek BA, Deodhar SD. Tissue immunocytometric evaluation of monoclonality of B-cell lymphomas: comparison with cell suspension studies. *Am J Clin Pathol*. 1981;76:24-28.
 69. Smith BR, Weinberg DS, Robert NJ, et al. Circulating monoclonal B lymphocytes in non-Hodgkin's lymphoma. *N Engl J Med*. 1984;311:1476-1481.

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